

STEMdiff™ Hematopoietic Kit

For the Generation of Hematopoietic Progenitor Cells from hPSCs

The **STEMdiff™ Hematopoietic Kit** consists of defined, serum-free basal medium and supplements designed for the generation of hematopoietic progenitor cells (HPCs) from human embryonic stem (hES) cells and induced pluripotent stem (hiPS) cells. Optimized for a standardized, 12-day differentiation protocol, this kit supports robust differentiation of human pluripotent stem cells (hPSCs) into HPCs that can be identified by the expression of CD34 and CD45 (Figures 1-3) and by the ability to form hematopoietic colonies of multiple lineages in colony-forming unit (CFU) assays with MethoCult™ medium (Figure 4). This kit is formulated for use in feeder-free conditions, optimized for the differentiation of hPSCs maintained in TeSR™ media and compatible with multiple hES and hiPS cell lines. After differentiation, the resulting HPCs may be used for additional downstream assays.



Advantages:

DEFINED. Serum-free and feeder-free formulation.

EASY-TO-USE. Simple monolayer protocol produces HPCs in suspension for easy harvest.

RAPID. Generation of HPCs in 12 days.

HIGH YIELD. One kit typically generates 4 - 18 million CD34⁺CD45⁺ HPCs.

FLEXIBLE. Robust generation of HPCs across multiple hES and hiPS cell lines.

PRODUCT	SIZE	CATALOG #
STEMdiff™ Hematopoietic Kit*	1 Kit	05310
STEMdiff™ Hematopoietic Basal Medium	120 mL	05311
STEMdiff™ Hematopoietic Supplement A (200X)	225 µL	05312
STEMdiff™ Hematopoietic Supplement B (200X)	375 µL	05313

* Kit includes basal medium and supplements A and B.




Figure 1. Hematopoietic Differentiation Protocol

One day before the differentiation protocol, hPSC colonies are harvested and seeded as small aggregates (100 - 200 µm in diameter) at 10 - 20 aggregates/cm² in mTeSR™1 or TeSR™-E8™. After one day, TeSR™ medium is replaced with Medium A (STEMdiff™ Hematopoietic Basal Medium containing Supplement A) to begin inducing the cells towards a mesoderm-like state (day 0). On day 2, a half medium change is performed with fresh Medium A. On day 3, the medium is changed to Medium B (STEMdiff™ Hematopoietic Basal Medium containing Supplement B) with half medium changes on days 5, 7 and 10, to promote further differentiation into hematopoietic cells. Typically, by day 12, large numbers of HPCs can be harvested from the culture supernatant.

STEMdiff™ Hematopoietic Kit

For the Generation of Hematopoietic Progenitor Cells from hPSCs




Figure 2. Morphology of hPSC-Derived HPCs

Representative images of (A) hES (H1) cells and (B) hiPS (WLS-1C) cells on day 12 of differentiation to HPCs using the STEMdiff™ Hematopoietic Kit. Differentiated cells exhibit typical HPC morphology as round cells that float freely in suspension.




Figure 3. Efficient and Robust Generation of CD34⁺CD45⁺/CD43⁺ HPCs

hES and hiPS cells were cultured for 12 days in single wells of 12-well plates using the STEMdiff™ Hematopoietic Kit. At the end of the culture period, cells in suspension were harvested and analyzed by flow cytometry for expression of hematopoietic cell surface markers: CD34, CD45 and CD43. (A,B) Example flow cytometry plots for hematopoietic cell surface-marker analysis of cultures of hES (H1 and H9) and hiPS (STiPS-M001) cells. (C,D) Percentages and total numbers of CD34⁺CD45⁺ cells in cultures of hES (H1 and H9) or hiPS (WLS-1C, STiPS-F016, STiPS-M001 and STiPS-B004) cells are shown. Data shown as mean ± SEM; n ≥ 3.





Figure 4. hPSC-Derived HPCs Produce Colonies of Multiple Lineages

Cells in suspension were harvested from cultures on day 12 of the hematopoietic differentiation protocol and assessed in colony-forming unit (CFU) assays using MethoCult™ H4435 Enriched (Catalog #04435) methylcellulose-based medium. Shown are (A) CFU frequencies observed in cultures of HPCs-derived from 6 different hPSC lines (data shown as mean ± SEM; n ≥ 3.) CFU frequencies were variable between different cell lines, with on average approximately 120 CFU per 10,000 hPSC-derived HPCs plated. The progenitor cell types observed included granulocyte/macrophage (CFU-M, CFU-G and CFU-GM), erythroid (BFU-E and CFU-E) and occasional mixed (CFU-GEMM) colonies. (B) A selection of representative colony images are shown at 40X magnification. MethoCult™ SF H4636 (Catalog #04636) has also been tested for use in CFU assays of hPSC-derived HPCs in serum-free conditions.

For a complete list of related products, including specialized cell culture and storage media, matrices, antibodies, cytokines and small molecules, visit www.stemcell.com/hPSCworkflow or contact us at techsupport@stemcell.com.