

# MOUSE T CELL RESEARCH

Cell Isolation and Culture



# EasySep™

## The Easy Choice for Simple Cell Isolation

EasySep™ is a powerful immunomagnetic cell isolation platform that combines the specificity of monoclonal antibodies with the simplicity of a column-free magnetic system, offering easy and fast isolation of highly purified cell populations that are immediately ready for a wide range of downstream applications. EasySep™ protocols are optimized to minimize hands-on time, making cell isolation simpler and more efficient. This versatile technology supports positive selection, negative selection, or cell depletion and can be used to isolate cells in as little as 8 minutes from a variety of species (e.g. human, mouse, non-human primate, and rat) and sample sources, including splenocytes, lymph nodes, peripheral blood mononuclear cells (PBMCs), whole blood, leukopaks, cord blood, and bone marrow.



### Why Use EasySep™ to Isolate Mouse Immune Cells?





- Isolate cells in as little as 15 minutes, with a simple pour.
- Achieve up to 99% cell purities with high recoveries.
- Obtain viable, functional cells without the need for columns and washes.
- Isolate cells from virtually any sample source.
- Choose a method backed by 20+ years of published research.



### Request a Sample

Try EasySep™ Column-Free Cell Isolation  
[stemcell.com/Try-EasySep](https://stemcell.com/Try-EasySep)

## EasySep™ Magnets

	EasySep™ Magnet <sup>1</sup>	"The Big Easy" EasySep™ Magnet <sup>2</sup>	Easy 50 EasySep™ Magnet	Easy 250 EasySep™ Magnet	EasyEights™ EasySep™ Magnet	EasyPlate™ EasySep™ Magnet
						
Catalog #	18000	18001	18002	100-0821	18103	18102
Number of Samples	1	1	1	1	8 on each side = 16 total	96
Start Sample Cell Number Range <sup>3</sup>	0.1 - 2.5 x 10 <sup>8</sup> cells per 5 mL tube	0.2 - 10 x 10 <sup>8</sup> cells per 14 mL tube	0.5 - 20 x 10 <sup>8</sup> cells per 50 mL tube	2.00 - 12.5 x 10 <sup>9</sup> cells	0.125 - 2.0 x 10 <sup>8</sup> cells per 5 mL tube 0.25 - 8.0 x 10 <sup>8</sup> cells per 14 mL tube	0.025 - 0.2 x 10 <sup>8</sup> cells per well
Collection Method	Pour off	Pour off	Pipette off	Pipette off	Pipette off	Pipette off



1. Multiple EasySep™ magnets can be used together with the EasySep™ Multistand (Catalog #18010) for processing up to 4 samples simultaneously or with up to 6 EasySep™ EasyStands™ (Catalog #18130) for processing up to 6 samples.  
2. Multiple "The Big Easy" EasySep™ magnets can be used together with the EasySep™ Multistand (Catalog #18010) for processing up to 4 samples simultaneously.  
3. Minimum and maximum cell number range and volumes depend on the cell isolation kit, sample source, cell type being isolated, and flasks used during isolation.

# EasySep™ for Mouse Cells

## Immunomagnetic Cell Isolation in As Little As 15 Minutes


EasySep™ kits are available for the isolation of mouse cells in as little as 15 minutes from a variety of sample sources, including spleen, bone marrow, lymph nodes, and whole blood. Our EasySep™ negative selection kits use biotinylated antibodies to target unwanted cells, while positive selection kits use an antibody complex to select cells of interest.

### Typical EasySep™ Mouse Cell Isolation Protocol (Negative Selection)



**Figure 1.** EasySep™ Mouse CD4<sup>+</sup> T Cell Isolation Kit (Catalog #19852)

Starting with mouse splenocytes, the CD4<sup>+</sup> T cell content (CD3<sup>+</sup>CD4<sup>+</sup>) of the isolated fraction is  $95.4 \pm 3\%$  (mean  $\pm$  SD using the purple EasySep™ magnet). In the above example, the purities of the start and isolated CD3<sup>+</sup>CD4<sup>+</sup> fractions are 16.4% and 95.1%, respectively.



**Figure 2.** EasySep™ Mouse F4/80 Positive Selection Kit (Catalog #100-0659)

Starting with a naïve mouse lung single-cell suspension, the F4/80<sup>+</sup> cell content of the isolated fraction is typically  $94.3 \pm 2.8\%$  (mean  $\pm$  SD using the purple EasySep™ magnet). In the above example, the purities of the start and final isolated fractions are 26.3% and 95.3%, respectively.

## Mouse Cell Isolation Products

### Mouse T Cell and Subset Isolation by Negative Selection

Cell Type	Source	Product	Purity <sup>1</sup>	For Processing	Compatible Staining Antibodies	Catalog #
T Cells	Spleen	EasySep™ Mouse T Cell Isolation Kit	96.6 ± 2.0%	1 x 10 <sup>9</sup> cells	CD3e (Catalog #60015) CD90 (Catalog #60024)	19851 19851RF
Pan-Naïve T Cells	Spleen	EasySep™ Mouse Pan-Naïve T Cell Isolation Kit	90.0 - 97.0%	1 x 10 <sup>9</sup> cells		19848 19848RF
CD4 <sup>+</sup> T Cells	Spleen	EasySep™ Mouse CD4 <sup>+</sup> T Cell Isolation Kit	95.4 ± 3%	1 x 10 <sup>9</sup> cells	CD3e (Catalog #60015) CD4 (Catalog #60017)	19852 19852RF
Naïve CD4 <sup>+</sup> T Cells	Spleen	EasySep™ Mouse Naïve CD4 <sup>+</sup> T Cell Isolation Kit	90.0 - 95.0%	1 x 10 <sup>9</sup> cells	CD4 (Catalog #60017) CD44 (Catalog #60068) CD62L (Catalog #60109)	19765 19765RF
Memory CD4 <sup>+</sup> T Cells	Spleen	EasySep™ Mouse Memory CD4 <sup>+</sup> T Cell Isolation Kit	78.0 - 96.0%	1 x 10 <sup>9</sup> cells		19767 19767RF
CD8 <sup>+</sup> T Cells	Spleen	EasySep™ Mouse CD8 <sup>+</sup> T Cell Isolation Kit	94.4 ± 0.7%	1 x 10 <sup>9</sup> cells	CD3e (Catalog #60015) CD8a (Catalog #60023)	19853 19853RF
Naïve CD8 <sup>+</sup> T Cells	Spleen	EasySep™ Mouse Naïve CD8 <sup>+</sup> T Cell Isolation Kit	92.0 - 98.0%	1 x 10 <sup>9</sup> cells	CD8a (Catalog #60023) CD44 (Catalog #60068) CD62L (Catalog #60109)	19858 19858RF

### Mouse T Cell and Subset Isolation by Positive Selection

Cell Type	Source	Product	Purity <sup>1</sup>	For Processing	Compatible Staining Antibodies	Catalog #
CD90.2+ (Thy 1.2) Cells	Spleen	EasySep™ Mouse CD90.2 Positive Selection Kit II	97.0 - 99.0%	2 x 10 <sup>9</sup> cells	CD90.2	18951 18951RF
CD90.1+ (Thy 1.1) Cells	Spleen, Lymph Node, Whole Blood	EasySep™ Mouse CD90.1 Positive Selection Kit	93.5 ± 3.9%	2 x 10 <sup>9</sup> cells	CD90 (Catalog #60024)	18958 18958RF
CD4 <sup>+</sup> Cells	Spleen	EasySep™ Mouse CD4 Positive Selection Kit II	94.8 ± 3.5 %	2 x 10 <sup>9</sup> cells	CD3e (Catalog #60015) CD4 (Catalog #60029)	18952 18952RF
CD4+CD62L+ Cells	Spleen	EasySep™ Mouse CD4+CD62L+ T Cell Isolation Kit	92.0 - 97.0%	1 x 10 <sup>9</sup> cells	CD4 (Catalog #60017) CD44 (Catalog #60068)	18765 18765RF
CD4+CD304+ Cells	Spleen	EasySep™ Release Mouse CD4+CD304+ Regulatory T Cell Isolation Kit	90.4 ± 3.0%	1 x 10 <sup>9</sup> cells	CD4 (Catalog #60029) CD4 (Catalog #60017) FOXP3, clone FJK-16s	100-1570 100-1564
	Lymph Nodes		92.1 ± 0.7%			
	Spleen		94.1 ± 2.9%			
CD4+CD25+ Cells	Spleen	EasySep™ Mouse CD4+CD25+ Regulatory T Cell Isolation Kit II	70.0 - 93.0%	1 x 10 <sup>9</sup> cells	CD4 (Catalog #60029) CD4 (Catalog #60017)	18783
CD25+ Cells	Spleen	EasySep™ Mouse CD25 Regulatory T Cell Positive Selection Kit	80.0 - 93.0%	1 x 10 <sup>9</sup> cells	CD4 (Catalog #60017)	18782 18782RF
CD8a+ Cells	Spleen	EasySep™ Mouse CD8a Positive Selection Kit II	96.3 ± 1.4%	2 x 10 <sup>9</sup> cells	CD3e (Catalog #60015) CD8a (Catalog #60023)	18953 18953RF

RoboSep™-S Reagent Kits (RF) contain an EasySep™ Selection Kit with RoboSep™ Buffer and 1 - 2 boxes of RoboSep™ Tip Racks.

1. Purities shown as either a range or mean ± SD.

# Culture Cells with ImmunoCult™

## Activate, Expand, Maintain, and Differentiate Immune Cells

Ensure optimal activation, expansion, and differentiation of immune cell subsets by using STEMCELL Technologies' ImmunoCult™ cell culture media and supplements. ImmunoCult™ products allow you to culture various cell types, including monocytes, T cells, NK cells, B cells, dendritic cells, and macrophages, under defined stimulatory conditions for consistent and reliable results.

ImmunoCult™ products are part of a wider and complete immunology workflow of STEMCELL products.

### ImmunoCult™ Mouse T Cell Products

Product	Catalog #	Size
ImmunoCult™ Mouse T Cell Activator Kit	100-1572	3 x 1 mL
ImmunoCult™ -XF T Cell Expansion Medium	10981	500 mL
ImmunoCult™ Mouse Th1 Differentiation Supplement	10953	1 mL
ImmunoCult™ Mouse Th2 Differentiation Supplement	10955	1 mL
ImmunoCult™ Mouse Treg Differentiation Supplement	10957	1 mL

### Why Use ImmunoCult™?

- Activate, expand, or differentiate immune cells in culture conditions optimized to promote high yield and frequency.
- Minimize variation by using serum-free culture conditions.
- Consistently achieve high yields of immune cells with the desired phenotype and function.
- Mix and match media, activators, and supplements to suit your specific research needs.



### Wallchart

Frequencies & Percentages of Mouse Immune Cell Types  
[stemcell.com/MouseCellFreq](https://www.stemcell.com/MouseCellFreq)





Figure 3. Products for Your Immune Cell Culture Workflow




## ImmunoCult™ for Mouse T Cell Activation and Expansion

Activate mouse T cells with confidence using a single kit that allows for flexible design and activation by CD3/CD28 or CD3/CD28/CD2. The ImmunoCult™ Mouse T Cell Activator Kit provides a gentle activation stimulus that encourages high viability of activated T cells, which can then be further expanded in common mouse T cell growth media. For efficient isolation of mouse T cells, the ImmunoCult™ Mouse T Cell Activator Kit can be used with EasySep™ as a part of a complete T cell workflow.



**Figure 4.** Activation of EasySep™-Isolated Mouse T Cells Stimulated with ImmunoCult™ Mouse T Cell Activator Kit

Mouse T cells were isolated using EasySep™ Mouse T Cell Isolation Kit, stimulated with ImmunoCult™ Mouse T Cell Activator Kit, and cultured in IMDM + FBS formulation. Following 3 days of culture, the mean  $\pm$  SD frequency of CD25+ cells was  $91.9 \pm 5.1\%$  ( $n = 11$ ) or  $99.9 \pm 0.1\%$  ( $n = 5$ ) when stimulated with ImmunoCult™ Mouse CD3/CD28 T Cell Activator or ImmunoCult™ Mouse CD3/CD28/CD2 T Cell Activator, respectively. Stimulated mouse T cells maintained expression levels of CD25 throughout the 7-day culture period.




**Figure 5.** Robust Expansion of EasySep™-Isolated Mouse T Cells Can Be Achieved Following Stimulation with ImmunoCult™ Mouse T Cell Activator Kit

EasySep™-isolated mouse T cells were expanded with ImmunoCult™ Mouse T Cell Activator Kit in IMDM + FBS formulation over 7 days. The number of viable cells was assessed every 2 - 3 days, when fresh medium supplemented with IL-2 was added. No additional ImmunoCult™ Mouse T Cell Activator was added during the 7-day culture period. After 7 days in culture with ImmunoCult™ Mouse CD3/CD28 T Cell Activator or ImmunoCult™ Mouse CD3/CD28/CD2 T Cell Activator, stimulation resulted in a fold expansion of  $23 \pm 3.4$  or  $29.3 \pm 4.8$  (mean  $\pm$  SEM,  $n = 6$ ), respectively.

### Why Use ImmunoCult™ Activation and Expansion Reagents?

- Achieve optimal activation without the use of magnetic beads, feeder cells, or antigens.
- Maintain functional T cells with high viability by using a gentle activation stimulus.
- Rely on highly stable, filter-sterilized, and soluble reagents that are compatible with common T cell growth medium.



**Figure 6.** High Cell Proliferation Is Observed in EasySep™-Isolated T cells After Stimulation with ImmunoCult™ Mouse T Cell Activator

EasySep™-isolated mouse T cells were labeled with CFDA-SE (Catalog #75003), stimulated with ImmunoCult™ Mouse T Cell Activator Kit, and cultured in IMDM + FBS formulation. On Day 3, cells were harvested, stained with anti-mouse CD4 and CD8a antibodies, then measured by flow cytometry. Shown are CFDA-SE-labeled mouse T cells, gated on viable (A) CD4+ or (B) CD8a+ cells, cultured with no activator (top panel), with ImmunoCult™ Mouse CD3/CD28 T Cell Activator (middle panel), or with ImmunoCult™ Mouse CD3/CD28/CD2 T Cell Activator (bottom panel). Due to cell proliferation, the intensity of CFDA-SE signal is reduced by 50% for each cell division.




### Protocol Video

How to Prepare a Single-Cell Suspension from Mouse Spleen

[stemcell.com/Mouse-Sample-Processing](https://www.stemcell.com/Mouse-Sample-Processing)


## ImmunoCult™ for Mouse T Cell Differentiation

Differentiate mouse T cells with ImmunoCult™ supplements. These differentiation supplements contain recombinant mouse cytokines and blocking monoclonal antibodies formulated to promote the differentiation of spleen-derived mouse naïve CD4<sup>+</sup> T cells into Th1, Th2, or regulatory T (Treg) cells. ImmunoCult™ differentiation supplements are intended for use with RPMI 1640 Medium containing fetal bovine serum and other additives, as well as anti-mouse CD3 and CD28 monoclonal antibodies as activating agents.



**Figure 7. ImmunoCult™ Mouse Th1 Differentiation Supplement Produces CD4<sup>+</sup>IFNγ<sup>+</sup>T-Bet<sup>+</sup> Cells Under Th1 Polarizing Conditions**

Naïve CD4<sup>+</sup> T cells were isolated from mouse splenocytes using the EasySep™ Mouse Naïve CD4<sup>+</sup> T Cell Isolation Kit (Catalog #19765), activated with plate-bound anti-CD3 and soluble anti-CD28, and cultured in medium alone (non-polarized cultures) or medium supplemented with ImmunoCult™ Mouse Th1 Differentiation Supplement (polarized cultures) for 5 days. Cells were subsequently stimulated with PMA/Ionomycin, in conjunction with Brefeldin A treatment, and stained with anti-CD4, anti-IFNγ, anti-Tbet, and a viability dye, and analyzed by flow cytometry. Shown is the expression of IFNγ and Tbet, back-gated on viable CD4<sup>+</sup> cells from (A) non-polarized or (B) polarized cultures. The mean proportion of CD4<sup>+</sup>IFNγ<sup>+</sup>Tbet<sup>+</sup> cells is significantly higher in cells cultured in ImmunoCult™ Mouse Th1 Differentiation Supplement (44 ± 6%) compared to non-polarized cells (2 ± 1%;  $p < 0.001$ ;  $n = 13$ ). Data from experimental groups were compared using a paired T-test.




**Figure 8. ImmunoCult™ Mouse Th2 Differentiation Supplement Produces CD4<sup>+</sup>GATA3<sup>+</sup>IL-4<sup>+</sup> Cells Under Th2 Polarizing Conditions**

Naïve CD4<sup>+</sup> T cells were isolated from mouse splenocytes using the EasySep™ Mouse Naïve CD4<sup>+</sup> T Cell Isolation Kit (Catalog #19765), activated with plate-bound anti-CD3 and soluble anti-CD28, and cultured in medium alone (non-polarized cultures) or in medium supplemented with ImmunoCult™ Mouse Th2 Differentiation Supplement (polarized cultures) for 6 days. Cells were subsequently stimulated with PMA/Ionomycin, in conjunction with monensin treatment, and stained with anti-CD4, anti-IL-4, anti-GATA3, and a viability dye, and analyzed by flow cytometry. Shown is the expression of GATA3 and IL-4 back-gated on viable CD4<sup>+</sup> cells from (A) non-polarized or (B) polarized cultures. The mean proportion of CD4<sup>+</sup>IL-4<sup>+</sup>GATA3<sup>+</sup> cells is significantly higher in cells cultured in ImmunoCult™ Mouse Th2 Differentiation Supplement (25 ± 3%) compared to non-polarized cells (4 ± 1%;  $p < 0.001$ ;  $n = 10$ ). Data from experimental groups were compared using a paired T-test.

### Why Use ImmunoCult™ Differentiation Supplements?

- Induce Th1, Th2, or Treg cells from naïve CD4<sup>+</sup> T cells isolated from mouse spleen.
- Add 100x concentrated supplements directly to the medium.
- Use cell culture reagents that are compatible with upstream EasySep™ cell isolation.



**Figure 9. ImmunoCult™ Mouse Treg Differentiation Supplement Produces CD4<sup>+</sup>CD25<sup>+</sup>FOXP3<sup>+</sup> Cells Under Treg Polarizing Conditions**

Naïve CD4<sup>+</sup> T cells were isolated from mouse splenocytes using the EasySep™ Mouse Naïve CD4<sup>+</sup> T Cell Isolation Kit (Catalog #19765), activated with plate-bound anti-CD3 and soluble anti-CD28, and cultured in medium alone (non-polarized cultures) or in medium supplemented with ImmunoCult™ Mouse Treg Differentiation Supplement (polarized cultures) for 6 days. Cells were subsequently stained with anti-CD4, anti-CD25, anti-FOXP3, and a viability dye, and analyzed by flow cytometry. Shown is the expression of CD25 and FOXP3 back-gated on viable CD4<sup>+</sup> cells from (A) non-polarized or (B) polarized cultures. The mean proportion of CD4<sup>+</sup>FOXP3<sup>+</sup> cells is significantly higher in cells cultured in ImmunoCult™ Mouse Treg Differentiation Supplement (91 ± 2%) compared to non-polarized cells (2 ± 0.4%;  $p < 0.001$ ;  $n = 14$ ). Data from experimental groups were compared using a paired T-test.

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, ArciTeCt, EasyEights, EasyPlate, EasySep, EasyStand, ImmunoCult, RapidSpheres, and RoboSep are trademarks of STEMCELL Technologies Canada Inc. CryoStor and ThawSTAR are registered trademarks of BioLife Solutions, Inc. Neon is a registered trademark of Thermo Fisher Scientific. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

UNLESS OTHERWISE STATED, PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES. FOR PRODUCT-SPECIFIC COMPLIANCE AND INTENDED USE INFORMATION, REFER TO THE PRODUCT INFORMATION SHEET. GENERAL INFORMATION ON QUALITY AT STEMCELL MAY BE FOUND AT [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).

# MOUSE T CELL RESEARCH

Cell Isolation and Culture



TOLL FREE PHONE 1 800 667 0322

PHONE +1 604 877 0713

INFO@STEMCELL.COM

TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

DOCUMENT #27256 VERSION 1.2.0 MAR 2025