

CELL CULTURE MATRICES FOR HUMAN PLURIPOTENT STEM CELLS

Vitronectin XF™ and CellAdhere™ Laminin-521

Cell culture matrices support the growth and differentiation of human pluripotent stem cells (hPSCs), including human embryonic stem (ES) and induced pluripotent stem (iPS) cells, by mimicking the in vivo extracellular matrix. When used with TeSR™ maintenance media, our defined cell culture surfaces, [Vitronectin XF™](#) and [CellAdhere™ Laminin-521](#), provide a robust culture system for cell maintenance under feeder-free conditions.

Vitronectin XF™

For Growth and Differentiation of hPSCs Under Serum-Free, Feeder-Free Conditions

Vitronectin XF™, developed and manufactured by Nucleus Biologics, is a defined, xeno-free cell culture matrix that supports the growth and differentiation of hPSCs. Use with mTeSR™1 ([Catalog #85850](#)), mTeSR™ Plus ([Catalog #100-0276](#)), TeSR™-E8™ ([Catalog #05990](#)), or TeSR™-AOF ([Catalog #100-0401](#)) medium to provide a defined culture system for the maintenance of ES and iPS cells and greater control over the culture environment, resulting in more consistent cell populations and reproducible results in downstream applications. Human ES and iPS cells cultured on Vitronectin XF™ retain pluripotency and normal colony morphology, without the need for an adaptation step (Figure 1). Pair with Gentle Cell Dissociation Reagent (GCDR; [Catalog #07174](#)) or ReLeSR™ ([Catalog #05872](#)) when passaging to maintain high-quality cultures.

Product	Size	Catalog #
Vitronectin XF™	2 mL	07180
	10 x 2 mL	100-0763
CellAdhere™ Dilution Buffer	100 mL	07183



Why Use Vitronectin XF™?

- Decrease sources of variability in your experiment with a recombinant human protein matrix
- Handle at room temperature without matrix gelling
- Use with any TeSR™ family medium to maintain hPSCs
- Create a completely xeno-free system when used with TeSR™-E8™ or TeSR™-AOF




Figure 1. Morphology of Human ES and iPS Cells Cultured on Vitronectin XF™ Cell Culture Matrix in TeSR™-E8™

Undifferentiated human ES (H9) and iPS (WLS-1C) cell cultures exhibit normal morphology when cultured on Vitronectin XF™. Colonies are round, tightly packed and multilayered, with a high nucleus-to-cytoplasm ratio. Cells were transferred directly from Matrigel® hESC-Qualified Matrix without an adaptation step.

Note: Colonies grown in TeSR™-E8™ have a more condensed and round morphology when grown on Vitronectin XF™ matrix, compared to colonies grown on Matrigel® hESC-Qualified Matrix, which are more diffuse and irregularly shaped.

CellAdhere™ Laminin-521

For Long-Term Maintenance in Feeder-Free Conditions

CellAdhere™ Laminin-521 is a defined and xeno-free cell culture matrix that supports the growth and differentiation of ES and iPSC cells under feeder-free conditions. Laminins are membrane proteins that bind to cell surface receptors, which results in the activation of signaling cascades and leads to cells that are more physiologically relevant. Laminin-521 is expressed and secreted by hPSCs in the inner cell mass of the embryo and can therefore be used to create a biologically relevant hPSC culture environment in vitro.

For consistent cell populations and reproducible results in downstream applications, use CellAdhere™ Laminin-521 with TeSR™ maintenance media to provide a defined culture substrate for cell maintenance. Compared to other matrices, CellAdhere™ Laminin-521 increases single-cell attachment and survival and does not require the addition of apoptotic inhibitors for long-term culture. CellAdhere™ Laminin-521 is recommended for use with eTeSR™ (Catalog #100-1215) maintenance medium for single-cell passaging. Pair with Gentle Cell Dissociation Reagent (GCDR; Catalog #07174) or ReLeSR™ (Catalog #05872) for routine passaging of PSC aggregates, or ACCUTASE™ (Catalog #07920) for single-cell passaging workflows.

Note: Single-cell passaging of human ES and iPSC cells can result in selective pressure and lead to genetic aberrations. If passaging as single cells, check the karyotype frequently.

Why Use CellAdhere™ Laminin-521?

- Decrease sources of variability in your experiment with a recombinant human protein matrix
- Passage cells without the need for apoptotic inhibitors
- Use with any TeSR™ family medium to maintain hPSCs
- Increase cell attachment and survival when using with eTeSR™ for single-cell passaging
- Mimic the stem cell niche with this physiologically relevant matrix

Product	Size	Catalog #
CellAdhere™ Laminin-521	100 µg	77003
	500 µg	200-0117

Combine Vitronectin XF™ or CellAdhere™ Laminin-521 and TeSR™ media with human iPSC lines, such as iPSCdirect™ (Catalog #100-1028) and SCTi003-A (Catalog #200-0511), for a fully defined culture system.

Table 1. Comparison of the Different Matrices

Features	Corning® Matrigel®	Vitronectin-XF™	CellAdhere™ Laminin-521
Key Advantage	hESC-qualified alternative to feeders	Defined matrix that supports the growth and differentiation of hPSCs	Defined matrix found to increase single-cell attachment and survival of hPSCs without apoptotic inhibitors
Grade	Feeder-free	Xeno-free	Xeno-free
Source	Solubilized Engelbreth-Holm-Swarm (mouse) sarcoma extract	Recombinant human protein	Recombinant human protein
Coating Temperature	2 - 8°C	15 - 25°C (room temperature)	2 - 8°C
Plasticware	Tissue culture-treated	Non-tissue culture-treated	Tissue culture-treated

For a complete list of related products, including specialized cell culture and storage media, matrices, antibodies, cytokines, and small molecules, visit www.stemcell.com/hPSCworkflow or contact us at techsupport@stemcell.com.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, CellAdhere, iPSCdirect, and ReLeSR are trademarks of STEMCELL Technologies Canada Inc. CellAdhere Laminin-521 is manufactured by BioLamina. ACCUTASE is a trademark of Innovative Cell Technologies Inc. TeSR, E8, mTeSR, and eTeSR are trademarks of WARF. Vitronectin-XF is a trademark of, and developed and manufactured by Nucleus Biologics. Corning and Matrigel are registered trademarks of Corning Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

UNLESS OTHERWISE STATED, PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES. FOR PRODUCT-SPECIFIC COMPLIANCE AND INTENDED USE INFORMATION, REFER TO THE PRODUCT INFORMATION SHEET. GENERAL INFORMATION ON QUALITY AT STEMCELL MAY BE FOUND AT WWW.STEMCELL.COM/COMPLIANCE.