A microscopic view of a cell culture dish showing a large, dense cluster of small, round cells in the center, with many smaller, individual cells scattered throughout the field of view. The cells have a distinct nucleus and a thin cytoplasmic layer.

HEMATOPOIETIC STEM AND PROGENITOR CELLS

Products for Your Research



TABLE OF CONTENTS

Introduction

- 3 [Your Ideas. Our Tools.](#)
- 4 [Hematopoietic Cell Therapy Development Workflow](#)

Source and Isolate

- 5 [Primary Human Hematopoietic Cells](#)
- 8 [Cell Isolation Products & Platforms](#)
- 8 [EasySep™](#)
- 8 [RoboSep™](#)
- 9 [SepMate™](#)
- 10 [RosetteSep™](#)
- 14 [EasySep™ RBC Depletion Agent](#)
- 15 [HetaSep™](#)
- 16 [Cryopreservation Media](#)

Expand & Differentiate

- 17 [StemSpan™ Expansion Media & Supplements](#)
- 18 [StemSpan™-AOF](#)
- 20 [CD34 Expansion Supplements](#)
- 23 [Small Molecules](#)
- 25 [STEMdiff™ Kits](#)
- 26 [STEMdiff™ Megakaryocyte Kit](#)
- 26 [STEMdiff™ Erythroid Kit](#)
- 27 [Recombinant Cytokines](#)

Characterize

- 28 [MyeloCult™](#)
- 30 [MethoCult™](#)
- 31 [The Colony Forming Unit \(CFU\) Assay](#)
- 34 [STEMvision™](#)
- 37 [SmartDish™ and STEMgrid™-6](#)
- 38 [MegaCult™ Collagen-Based Media](#)
- 39 [Antibodies & ELISA Kits](#)
- 40 [ALDH^{br} Assay Kit](#)

Support Products & Services

- 41 [Tissue Culture Reagents and Supplies](#)
- 42 [Proficiency Testing and Quality Control Kits](#)
- 43 [Training Courses and Instructional Materials](#)
- 44 [Contract Assay Services](#)
- 45 [Services for Cell Therapy Program](#)

References

- 46 [References](#)

Your Ideas. Our Tools.

Products for Every Step of Your Hematopoietic Stem and Progenitor Cell Research

STEMCELL Technologies is the world leader in developing tools for hematology ranging from fundamental research to cell therapy development. Our portfolio includes a comprehensive range of products for cell sourcing and isolation, expansion and differentiation, and analysis of hematopoietic stem and progenitor cells (HSPCs). To help ensure standardization throughout your HSPC research, use STEMCELL products from the beginning to the end of your workflow.




STEMCELL Products for Every Step of Your HSPC Research

| Cell Sourcing & Isolation | Expansion & Differentiation | Analysis |
|---|---|---|
| Human Primary Hematopoietic Cells Cell Isolation Products & Platforms EasySep™ RoboSep™ SepMate™ RosetteSep™ HetaSep™ Cryopreservation Media | StemSpan™ Media & Kits StemSpan™ Supplements STEMdiff™ Kits for Differentiation of hPSCs ¹ to Hematopoietic Progenitor Cells Recombinant Cytokines Small Molecules | MyeloCult™ Media MethoCult™ Media STEMvision™ Instrument MegaCult™ Media Antibodies and ELISA Kits Proficiency Testing Programs Contract Assay Services |

1. hPSCs: human pluripotent stem cells

Visit www.stemcell.com/HSPCworkflow for a full listing of products for your HSPC research.

Hematopoietic Cell Therapy Development Workflow



1. Certain products are only available in select territories. Please contact your regional Sales representative or Product and Scientific Support at techsupport@stemcell.com for further information.

2. MNCs: Mononuclear cells

3. Explore the [Human CD34+ Cell Isolation Product Selection infographic](#) for a complete list of recommendations.

4. For large-scale isolation of human CD34+ cells from whole blood and normal or mobilized leukopaks, see the large-format (1x10¹⁰ cells) kit (Catalog # [100-1569](#)).

5. Additional validation for specific application may be required.

Human Primary Hematopoietic Cells

It All Starts with the Right Cells

Cell Sourcing & Isolation

Expansion & Differentiation

Analysis

Starting with the right primary cells builds a strong foundation for your experiments and is the first step toward success in your research. Choose from a wide range of fresh and cryopreserved human primary cells for your downstream applications.*

STEMCELL's cryopreserved hematopoietic cells are isolated from human cord blood, bone marrow, peripheral blood, and mobilized peripheral blood leukopaks. Whole bone marrow, whole peripheral blood, and leukapheresis (leukopak) preparations are also available for users requiring fresh, unprocessed tissue. All products are verified for purity and viability—ensuring reproducibility across multiple experiments.

Visit www.stemcell.com/primarycells for a complete list of mononuclear cells, isolated subsets, plasma, serum, and unprocessed tissues.

Donor Criteria

All human primary cell products are ethically sourced using consent forms and protocols approved by either an Institutional Review Board, the U.S. Food and Drug Administration, the U.S. Department of Health and Human Services, and/or an equivalent regulatory authority. Donations performed in the United States comply with applicable federal, state, and local laws, regulations, and guidance. Cells sourced in the United Kingdom are collected using protocols and informed consent forms (ICFs) approved by the National Health Service (and Health Research Authority) Research Ethics Committee (REC) or equivalent agency. Donors are prescreened for general health and viral status, including HIV-1, HIV-2, hepatitis B, and hepatitis C**. Additional screening or analysis is available upon request. Our Quality Assurance, Quality Control, and Regulatory Affairs departments are ready to assist you with any necessary documentation to meet specific institutional requirements.



Resource

Explore Frequently Asked Questions on Primary Cells
www.stemcell.com/primarycellsfaqs

Why Use Human Primary Cells from STEMCELL Technologies?

EFFICIENT. Reduce time spent collecting and culturing primary cells.

CONVENIENT. Reserve large numbers of cryopreserved cells and start experiments on your schedule with cells you've already tested.

CUSTOMIZABLE. Request custom products for non-standard cell types or collections with specific requirements.

PHYSIOLOGICALLY RELEVANT. Choose cells that are more physiologically representative of cells in vivo.

ETHICALLY SOURCED. Access donor samples collected using regulatory authority-approved consent forms and protocols.



Figure 1. Fresh Whole Bone Marrow

Whole Bone Marrow (Catalog #70502) is collected using heparin as the anticoagulant and supplied in a 100 mL bottle.

*Certain fresh and cryopreserved products are only available in select territories. Please contact your regional Sales representative or Product and Scientific Support at techsupport@stemcell.com for further information.

** Leukopak, Whole Blood, and Bone Marrow (LP, WB, and BM) Donor Screening: Donors are screened for HIV-1, HIV-2, hepatitis B, and hepatitis C. Donors in the U.K. are also screened for HTLV III and syphilis. For fresh mobilized leukopaks donors are screened for HIV-1, HIV-2, hepatitis B, hepatitis C, HTLV-II-II, syphilis, and WNV. Cryopreserved and fresh LP, WB, and BM: If the donor has been screened within 90 days prior to donation and the results are negative, the product will be shipped with the negative test result and date of the most recent viral testing on the Certificate of Analysis (CoA). If the donor has not been screened within 90 days prior to collection, a test sample will be taken at the time of collection (with the exception of cancer patient donors) and the product will be shipped before the screening results are available. In the event that a test result is positive, the customer will be contacted as soon as possible (usually within 2 - 4 business days from the time of shipment). Fresh Cord Blood (CB) Donor Screening: Testing for HIV-1, HIV-2, hepatitis B, and hepatitis C is performed on a maternal blood sample and/or on a sample of the donated CB. Cryopreserved CB: Products with negative test results from donor screening are shipped with the CoA.

Product Listing

Fresh Human Bone Marrow Products*

| Description | Anticoagulant | Quantity | Catalog # |
|-------------------|---------------|----------|-----------|
| Whole Bone Marrow | Heparin | > 25 mL | 70502.2 |
| | | > 50 mL | 70502.1 |
| | | > 100 mL | 70502 |

Leukopaks*

| Description | Anticoagulant | Quantity | Catalog # |
|--|-------------------|--------------|-----------|
| Fresh Peripheral Blood Leukopak ¹ | ACDA ² | Tenth Size | 200-0092 |
| | | Quarter Size | 70500.2 |
| | | Half Size | 70500.1 |
| | | Full Size | 70500 |
| Frozen Peripheral Blood Leukopak | ACDA ² | Tenth Size | 200-0470 |
| | | Quarter Size | 200-0132 |
| | | Half Size | 200-0131 |
| | | Full Size | 200-0130 |

Cryopreserved Human Umbilical Cord Blood Products*

| Description | Quantity | Catalog # |
|---|-------------------|-----------|
| Mononuclear Cells | 15 million cells | 70007.1 |
| | 50 million cells | 70007.2 |
| | 150 million cells | 70007 |
| CD34+ Cells (Mixed Donor) | 0.2 million cells | 70008.1 |
| | 0.5 million cells | 70008.3 |
| | 1 million cells | 70008 |
| CD34+ Cells (Single Donor) ³ | 0.2 million cells | 70008.2 |
| | 0.5 million cells | 70008.4 |
| | 1 million cells | 70008.5 |

Cryopreserved Human Bone Marrow Products*

| Description | Quantity | Catalog # |
|--------------------------|-------------------|-----------|
| Mononuclear Cells | 5 million cells | 70001.1 |
| | 15 million cells | 70001.2 |
| | 25 million cells | 70001 |
| | 50 million cells | 70001.3 |
| | 100 million cells | 70001.4 |
| | | |
| CD34+ Cells | 0.1 million cells | 70002.1 |
| | 0.3 million cells | 70002.2 |
| | 0.5 million cells | 70002.3 |
| | 1 million cells | 70002 |
| | 2 million cells | 70002.4 |
| | 5 million cells | 70002.5 |
| CD36+ Cells ⁴ | 1 million cells | 70003 |

Cryopreserved Human Peripheral Blood Cells*

| Description | Quantity | Catalog # |
|-------------------|-------------------|-----------|
| Mononuclear Cells | 15 million cells | 70025.1 |
| | 25 million cells | 70025.2 |
| | 50 million cells | 70025.3 |
| | 100 million cells | 70025 |

Cryopreserved Mobilized Human Peripheral Blood Products*

| Description | Quantity | Catalog # |
|--|-------------------|-----------|
| G-CSF Mobilized Mononuclear Cells | 5 million cells | 70049.4 |
| | 15 million cells | 70049.2 |
| | 25 million cells | 70049.3 |
| | 50 million cells | 70049.1 |
| | 100 million cells | 70049 |
| G-CSF Mobilized CD34+ Cells | 0.2 million cells | 70060.2 |
| | 1 million cells | 70060.1 |
| | 5 million cells | 70060 |
| | 10 million cells | 70060.3 |
| | 20 million cells | 70060.4 |
| G-CSF and Plerixafor Mobilized Mononuclear Cells | 5 million cells | 70072.4 |
| | 15 million cells | 70072.2 |
| | 25 million cells | 70072.3 |
| | 50 million cells | 70072.1 |
| | 100 million cells | 70072 |
| G-CSF and Plerixafor Mobilized CD34+ Cells | 0.2 million cells | 70073.2 |
| | 1 million cells | 70073.1 |
| | 5 million cells | 70073 |
| | 10 million cells | 70073.3 |
| | 20 million cells | 70073.4 |
| Plerixafor Mobilized Mononuclear Cells | 5 million cells | 70074.4 |
| | 15 million cells | 70074.2 |
| | 25 million cells | 70074.3 |
| | 50 million cells | 70074.1 |
| | 100 million cells | 70074 |
| Plerixafor Mobilized CD34+ Cells | 0.2 million cells | 70075.2 |
| | 1 million cells | 70075.1 |
| | 5 million cells | 70075 |
| | 10 million cells | 70075.3 |
| | 20 million cells | 70075.4 |

For a complete listing of fresh and cryopreserved products, visit www.stemcell.com/human-primary-cell-products.

* Certain fresh and cryopreserved products are only available in select territories. Please contact your Sales representative or Product and Scientific Support at techsupport@stemcell.com for further information.

1. A full size leukopak typically contains $1.1 \pm 0.3 \times 10^{10}$ cells and has a volume of approximately 120 mL.
2. ACDA: acid citrate dextrose solution A.
3. Additional sizes such as 0.6 million cells, 0.7 million cells, and 0.8 million cells are also available for umbilical cord blood-sourced CD34+ cells (Single Donor).
4. Cultured cell product.

















Mobilized Leukopaks^{1,2}


| Description | Anticoagulant | Quantity | Apheresis | Catalog # |
|---|-------------------|----------|--------------------------------|--------------------------------|
| Fresh Human Mobilized Peripheral Blood Leukopak, G-CSF | ACDA ³ | 1 bag | First Collection (Day 5) | 200-0602 |
| | | 2 bags | Second Collection (Day 6) | 200-0603 |
| | | | Both Collections (Day 5 and 6) | 100-1101 (200-0602 & 200-0603) |
| Fresh Human Mobilized Peripheral Blood Leukopak, Plerixafor | ACDA ³ | 1 bag | First Collection | 200-0604 |
| Fresh Human Mobilized Peripheral Blood Leukopak, G-CSF and Plerixafor | ACDA ³ | 1 bag | First Collection (Day 5) | 200-0607 |
| | | 2 bags | Second Collection (Day 6) | 200-0608 |
| | | | Both Collections (Day 5 and 6) | 100-1103 (200-0607 & 200-0608) |

¹ Certain cryopreserved products are only available in select territories. Please contact your Sales representative or Product and Scientific Support at techsupport@stemcell.com for further information.


² The cytokine G-CSF or filgrastim (Neupogen®) is the most commonly used mobilizing agent and is used as a gold standard in the clinic. The bicyclam plerixafor (Mozobil®) is a rapid-acting mobilization agent. The combination of both mobilization agents works synergistically, increasing the CD34+ mobilization compared to single-agent mobilization.

³ ACDA: acid citrate dextrose solution A.


| Mobilizing Agent | Mobilizing Regimen | Mobilization and Collection Schedule | | | | | | Quantity | Apheresis | Catalog # |
|---|--------------------|---|---|---|---|---|--|----------|--------------------------------|--------------------------------|
| | | Day 1 | Day 2 | Day 3 | Day 4 | Day 5 | Day 6 | | | |
| Granulocyte Colony-Stimulating Factor (G-CSF) | 5-Day |  |  |  |  |  | | 1 Bag | First Collection (Day 5) | 200-0602 |
| | | | | | |  |  | | Second Collection (Day 6) | 200-0603 |
| | | | | | | First Collection | Second Collection | 2 Bags | Both Collections (Day 5 and 6) | 100-1101 (200-0602 & 200-0603) |
| Plerixafor | 1-Day |  |  | | | | | 1 Bag | First Collection | 200-0604 |
| G-CSF + Plerixafor (Combo) | 5-Day |  |  |  |  |  | | 1 Bag | First Collection (Day 5) | 200-0607 |
| | | | | | |  |  | | Second Collection (Day 6) | 200-0608 |
| | | | | | | First Collection | Second Collection | 2 Bags | Both Collections (Day 5 and 6) | 100-1103 (200-0607 & 200-0608) |



G-CSF (NEUPOGEN®):
a maximum of 10 µg/kg/day



Plerixafor (MOZOBIL®):
a maximum of 0.24 mg/kg/day



One leukopak bag collected on the
specified day of the collection schedule

Figure 2. Mobilized Leukopak Selection Guide

Normal donors are mobilized with specified doses of G-CSF, plerixafor, or both mobilizing agents, based on the regimen shown above, prior to collection. The time interval between injection and apheresis can be found listed in the Certificate of Analysis.

Cell Isolation Products & Platforms

Fast and Easy Hematopoietic Stem and Progenitor Cell Isolation

We have combined our years of technical expertise in hematopoietic stem cell research with our powerful cell separation systems to develop a wide range of optimized hematopoietic stem and progenitor cell (HSPC) isolation products. Our innovative cell separation platforms provide fast, easy, and effective methods for isolating HSPCs with high purity and recovery.

EasySep™

The Easy Choice for Simple Cell Isolation

Isolating HSPCs can be a challenge due to their low frequency in tissues. EasySep™ column-free immunomagnetic cell isolation technology is the ideal method for the gentle, high-purity isolation of HSPCs and can be used with a variety of sources, including cord blood, whole blood, and bone marrow.

Using EasySep™, HSPCs can be easily and quickly isolated based on lineage-specific markers, such as expression of CD34 (for human cells), or c-KIT and SCA1 (for mouse cells). EasySep™ can be used either manually or with RoboSep™, the fully automated cell separator.

Why Use EasySep™ to Isolate Cells?


EASY & EFFICIENT. Isolate cells in as little as 8 minutes with a simple pour.

HIGH PURITY & RECOVERY. Achieve up to 99% cell purities with high recoveries.

COLUMN-FREE. Obtain viable, functional cells without the need for columns and washes.

VERSATILE. Isolate cells from virtually any sample source, including whole blood and leukopaks.

PROVEN. Choose a method widely used in published research for over 20 years to support diverse downstream applications.



**Times are typical for next-generation negative selection kits. Time for each kit will vary depending on the exact isolation protocol and magnet used. No particle incubation step is required for next-generation negative selection protocols.*

Figure 3. Typical EasySep™ Human Cell Isolation Protocol

RoboSep™

Fully Automated Cell Isolation System

Automate your cell processing and isolations using RoboSep™ cell isolation instruments. The RoboSep™-S and RoboSep™-16 instruments are true walk-away automated cell separation systems that retain the speed and simplicity of EasySep™ cell separation technology, with as little as 5 minutes of hands-on time. The column-free system minimizes sample handling and eliminates the risk of cross-contamination between samples. You can isolate cells from up to 16 samples at once for flow cytometry, functional studies, or other downstream applications.



SepMate™

Hassle-Free PBMC Isolation

SepMate™ is a specialized tube for fast and easy peripheral blood mononuclear cell (PBMC) isolation in just 15 minutes. The SepMate™ tube contains a unique insert that prevents the density gradient medium (e.g. Lymphoprep™) and blood sample from mixing. The density gradient medium is pipetted through a central hole in the insert, and the sample is rapidly pipetted or poured on top of the insert. This eliminates the need to carefully layer the sample directly onto the density gradient medium, an otherwise time-consuming and highly laborious step. Only 10 minutes of centrifugation are required, and this step can be carried out with the brake on, further reducing the total time necessary for separation. After centrifugation, plasma and PBMCs are simply poured into a new tube.

SepMate™ can be used on its own to isolate PBMCs in 15 minutes or combined with RosetteSep™ to enrich specific cell subsets, including CD34+ cells, directly from whole blood in 25 minutes. SepMate™ is available in 15 mL and 50 mL sizes for isolating individual samples of 0.5 - 17 mL in volume.




Figure 4. Typical SepMate™ Protocol

Why Use SepMate™?

- EASY.** Avoid slow and laborious sample layering over the density gradient medium.
- FAST.** Centrifuge for just 10 minutes with the brake on and simply pour off PBMCs into a new tube.
- CONSISTENT.** Eliminate errors and minimize variability between users.
- VERSATILE.** Combine with RosetteSep™ to isolate purified cell subsets from whole blood in 25 minutes.



Product Listing

| Product Name | Catalog # | Blood Volume Processed | Unit Size |
|---------------------------------|-----------|-------------------------|------------|
| SepMate™-15 (IVD ¹) | 85415 | 0.5 - 5 mL | 100 tubes |
| | 85420 | | 500 tubes |
| SepMate™-15 (RUO ²) | 86415 | 4 - 17 mL | 100 tubes |
| | 86420 | | 500 tubes |
| SepMate™-50 (IVD ¹) | 85450 | | 100 tubes |
| | 85460 | | 500 tubes |
| SepMate™-50 (RUO ²) | 86450 | | 100 tubes |
| | 86460 | | 500 tubes |
| Product Name | Catalog # | Density | Unit Size |
| Lymphoprep™ ⁴ | 18060 | 1.077 g/mL ³ | 250 mL |
| | 18061 | | 500 mL |
| | 07811 | | 4 x 250 mL |
| | 07861 | | 6 x 500 mL |


- SepMate™ is available as an in vitro diagnostic device for the isolation of mononuclear cells from human whole blood or bone marrow by density gradient centrifugation in Canada, the United States, Europe, and Australia. This product is also available in China where it is considered a non-medical device by the China Food and Drug Administration (CFDA) and should therefore be used as general laboratory equipment.
- SepMate™ RUO is available in other regions for research use only.
- Lymphoprep™ has the same density as Ficoll-Paque® and can be substituted for Ficoll-Paque® without any need to change your existing protocols.
- Lymphoprep™ is for Research Use Only (RUO).

RosetteSep™

Unique Immunodensity Cell Isolation

RosetteSep™ is a fast and easy immunodensity procedure for the isolation of untouched cells directly from whole blood. By crosslinking unwanted cells to red blood cells (RBCs) present in the sample, RosetteSep™ eliminates the need for a separate magnetic separation step because cells are purified during standard density gradient centrifugation. This approach significantly reduces sample handling time and maximizes convenience.

Typical RosetteSep™ Protocol



*Use SepMate™ to reduce centrifugation time to 10 minutes with brake on.




Figure 5. Micrograph of a Blood Sample After Addition of the RosetteSep™ Cocktail and Prior to Density Gradient Centrifugation




Figure 6. RosetteSep™ Human NK Cell Enrichment Cocktail (Catalog #15025)

Starting with whole peripheral blood, the NK cell content (gated on CD45+ cells) of the isolated fraction typically ranges from 80 - 98%. In the example above, the purities of the start and isolated fractions are 7.2% and 92.4%, respectively.

RosetteSep™ and SepMate™

Simplified and Standardized Cell Isolation

RosetteSep™ is easily combined with SepMate™ to rapidly and reproducibly isolate PBMC subsets from whole blood. By using the unique SepMate™ tube, sample throughput is increased and errors associated with improper sample layering are eliminated. This allows even users with minimal training to consistently perform cell isolation by density gradient centrifugation in a busy laboratory environment.

Human CD34+ Hematopoietic Progenitor Cell Isolation by Positive Selection

Applications

- Isolate human CD34+ cells from a variety of tissues.
- Obtain highly purified CD34+ cell populations for use in downstream assays.

Product Listing

| Source | Product Name | Purity ¹ | Capacity | Catalog # |
|-------------------------------|---|---------------------|---|------------------------------|
| Mobilized PBMCs, CBMCs, BMMCs | EasySep™ Human CD34 Positive Selection Kit II | 93.5 ± 1.1% | For labeling up to 5 x 10 ⁹ cells | 17856 ² , 17856RF |
| Mobilized Leukopak | EasySep™ Human CD34 Positive Selection Kit II | 94.6 ± 2.4% | For labeling up to 1 x 10 ¹⁰ cells | 100-1569 |
| Whole Blood, Buffy Coat | EasySep™ Human Whole Blood CD34 Positive Selection Kit | 26 - 41% | For labeling 75 mL whole blood (37 mL buffy coat) | 18086 |
| | EasySep™ Human Whole Blood CD34 Positive Selection Kit II | | | 17879, 17879RF |
| | RoboSep™ Human Whole Blood CD34 Positive Selection Kit | | | 18086RF |
| Whole Blood | Complete Kit for Human Whole Blood CD34+ Cells | 79 - 95% | For labeling 120 mL whole blood | 15086 |
| | RoboSep™ Complete Kit for Human Whole Blood CD34+ Cells | | | 15086RF |
| Fresh Cord Blood | EasySep™ Human Cord Blood CD34 Positive Selection Kit II | 91 ± 9% | For labeling 1000 mL cord blood | 17896 ^{3,4} |
| | RoboSep™ Human Cord Blood CD34 Positive Selection Kit II | | | 17896RF ^{3,4} |
| | EasySep™ Human Cord Blood CD34 Positive Selection Kit III | 87 ± 12% | | 17897 ⁴ |
| | RoboSep™ Human Cord Blood CD34 Positive Selection Kit III | | | 17897RF ⁴ |

PBMC: peripheral blood mononuclear cell; CBMC: cord blood mononuclear cell; BMMC: bone marrow mononuclear cell; hESC: human embryonic stem cell; hiPSC: human induced pluripotent stem cell.

For staining after positive selection with all Human CD34 Positive Selection kits, please use Anti-Human CD34, Clone 581 (Catalog #60013).

1. Purities shown as either a range or mean ± SD. Purity data for 18056, 18086, and 15086 are reported as a percentage of viable CD45+ cells.
2. These kits are for use with fresh or previously frozen PBMC, BMMC, and previously frozen cord blood mononuclear cells. For isolation of CD34+ cells from fresh cord blood, please use 18096RF, 17896, 17896RF, 17897, or 17897RF.
3. These kits are new versions of 18096 and 18096RF with improved performance.
4. For more information about choosing a CD34+ cell isolation kit for use with cord blood samples, please see our Technical Bulletin (Document #27003).

For a full listing of cell isolation tools and reagents, visit www.stemcell.com/HSPCworkflow under the "Cell Sourcing & Isolation" tab.

Human Hematopoietic Progenitor Cell Isolation by Negative Selection

Applications

- Enrich human hematopoietic progenitors (i.e. CD34+ cells) by using monoclonal antibodies targeted to specific cell surface antigens to deplete unwanted mature cell types (i.e. lineage depletion).
- Obtain unlabeled progenitor-enriched cell populations, i.e. lineage negative (Lin-) cells.

Product Listing

| Source | Product Name | Purity ¹ | Capacity | Catalog # |
|------------------------|--|---|--|----------------------|
| Mobilized PBMCs, BMMCs | EasySep™ Human Progenitor Cell Enrichment Kit | 42 ± 5 fold CD34+ cell enrichment (bone marrow) | For labeling up to 1 x 10 ⁹ cells | 19056 |
| | RoboSep™ Human Progenitor Cell Enrichment Kit | | | 19056RF |
| PBMCs, CBMCs | EasySep™ Human Progenitor Enrichment Kit with Platelet Depletion | 50 - 75% | For labeling up to 1 x 10 ⁹ cells | 19356 ² |
| | RoboSep™ Human Progenitor Enrichment Kit with Platelet Depletion | | | 19356RF ² |
| Bone Marrow | RosetteSep™ Human Bone Marrow Progenitor Cell Pre-Enrichment Cocktail | 25 ± 10 fold CD34+ cell enrichment | For labeling 40 mL bone marrow | 15027 |
| | | | For labeling 200 mL bone marrow | 15067 |
| Cord Blood | RosetteSep™ Human Hematopoietic Progenitor Cell Enrichment Cocktail II | 77.5 ± 16.0% | 1 x 10 ⁹ cells | 17936, 17936RF |
| | RosetteSep™ Human Hematopoietic Progenitor Cell Enrichment Cocktail | 29 ± 9% | For labeling 40 mL cord blood | 15026 |
| | | | For labeling 200 mL cord blood | 15066 |
| | Complete RosetteSep™ Human Cord Blood Progenitor Enrichment Kit | 29 ± 9% | For processing 500 mL cord blood | 15276 |
| Cord Blood | RosetteSep™ Human Cord Blood Debulking Cocktail | 5 ± 1% (CD34+ cells) | For labeling 40 mL cord blood | 15126 ³ |
| | | | For labeling 200 mL cord blood | 15166 ³ |

PBMC: peripheral blood mononuclear cell; BMMC: bone marrow mononuclear cell; CBMC: cord blood mononuclear cell.

For compatible staining antibodies, use 60013 Anti-Human CD34, 60018 Anti-Human CD45, 60026 Anti-Dextran.

1. Purities of CD34+ cells shown as either a range or mean ± SD. Purity of CD34+ cells for 19056 is reported relative to viable CD45+ cells in the start sample.

2. This product is designed for use with samples that contain large numbers of platelets.

3. This product is recommended for debulking cord blood of lineage positive cells prior to freezing.

For a full listing of cell isolation tools and reagents, visit www.stemcell.com/HSPCworkflow under the "Cell Sourcing & Isolation" tab.

Mouse Hematopoietic Progenitor Cell Isolation by Positive Selection

Applications

- Easily isolate mouse hematopoietic progenitor cells by positive selection.
- Select mouse hematopoietic progenitor cells using monoclonal antibodies to target specific cell surface antigens associated with progenitor cell phenotypes (e.g. SCA1+, cKIT+, AA4.1+).
- Achieve up to 97% purity using column-free immunomagnetic technology.

| Cell Type | Source | Product Name | Purity | Capacity | Compatible Anti-Mouse Staining Antibodies | Catalog # |
|---------------------|-------------|--|----------|--|--|-----------|
| SCA1+ Cells | Bone Marrow | EasySep™ Mouse SCA1 Positive Selection Kit | 87 - 97% | For labeling up to 2 x 10 ⁹ cells | CD3e Antibody, Clone 145-2C11 (Catalog #60015); CD11b Antibody, Clone M1/70 (Catalog #60001); CD19 Antibody, Clone 6D5 (Catalog #60006); CD45R Antibody, Clone RA3-6B2 (Catalog #60019); Gr-1 Antibody, Clone RB6-8C5 (Catalog #60028); TER119 Antibody, Clone TER-119 (Catalog #60033) | 18756 |
| | | RoboSep™ Mouse SCA1 Positive Selection Kit | | | | 18756RF |
| CD117+ (cKIT) Cells | Bone Marrow | EasySep™ Mouse CD117 (cKIT) Positive Selection Kit | 88 - 95% | For labeling up to 2 x 10 ⁹ cells | CD3e Antibody, Clone 145-2C11 (Catalog #60015); CD11b Antibody, Clone M1/70 (Catalog #60001); CD19 Antibody, Clone 6D5 (Catalog #60006); CD45R Antibody, Clone RA3-6B2 (Catalog #60019); Gr-1 Antibody, Clone RB6-8C5 (Catalog #60028); TER119 Antibody, Clone TER-119 (Catalog #60033) | 18757 |
| | | RoboSep™ Mouse CD117 (cKIT) Positive Selection Kit | | | | 18757RF |

Mouse Hematopoietic Progenitor Cell Isolation by Negative Selection

Applications

- Easily isolate untouched mouse hematopoietic progenitors (including Lin-SCA1+cKIT+ or LSK cells) by depleting unwanted mature cell types using monoclonal antibodies targeted to specific cell surface antigens expressed on these cells (i.e. lineage depletion).
- Efficiently obtain functional, unlabeled, progenitor-enriched cell populations, i.e. lineage negative (Lin-) cells.
- Achieve up to 84% purity using column-free immunomagnetic technology.

| Cell Type | Source | Product Name | Purity | Capacity | Compatible Anti-Mouse Staining Antibodies | Catalog # |
|--------------------------------|-------------|--|----------|--|---|-----------|
| Hematopoietic Progenitor Cells | Bone Marrow | EasySep™ Mouse Hematopoietic Progenitor Cell Isolation Kit | 60 - 84% | For labeling up to 1 x 10 ⁹ cells | CD117 Antibody, Clone 2B8 (Catalog #60025); Sca1 Antibody, Clone E13-161.7 (Catalog #60032); CD3e Antibody, Clone 145-2C11 (Catalog #60015); CD11b Antibody, Clone M1/70 (Catalog #60001); CD19 Antibody, Clone 6D5 (Catalog #60006); CD45R Antibody, Clone RA3-6B2 (Catalog #60019); Gr-1 Antibody, Clone RB6-8C5 (Catalog #60028); TER119 Antibody, Clone TER-119 (Catalog #60033) | 19856 |
| | | RoboSep™ Mouse Hematopoietic Progenitor Cell Isolation Kit | | | | 19856RF |

For a full listing of cell isolation tools and reagents for mouse cells, visit www.stemcell.com/hspc-mouse-workflow.

EasySep™ RBC Depletion Agent

Gentle, Lysis-Free Leukocyte Purification

In many laboratories, the standard protocols for obtaining leukocytes from human whole blood samples involve density gradient centrifugation or lysing red blood cells (RBCs) with ammonium chloride. However, these methods can be time consuming, be difficult to automate, and leave residual cell debris that may alter cellular function or interfere with downstream assays.

EasySep™ RBC Depletion Reagent immunomagnetically depletes RBCs without lysis, washes, or centrifugation steps. The resulting highly purified leukocytes are untouched and ready for downstream applications, including cell culture, RNA isolation, or enzyme activity testing. EasySep™ RBC Depletion Reagent can also be used on cord blood, bone marrow, buffy coats, and leukapheresis products to meet all of your laboratory’s needs for isolating leukocytes.

Why Use EasySep™ RBC Depletion Reagent?

GENTLE. Avoid lysis buffer, centrifugation, or additional washing steps.

RELIABLE. Deplete 99.9% of RBCs immunomagnetically without leaving additional debris that may interfere with downstream assays.

FAST. Obtain leukocytes in as little as 9 minutes.

CONVENIENT. Automate blood sample processing with RoboSep™ instruments to increase laboratory throughput.

EasySep™ RBC Depletion Reagent

| Product Name | Volume Processed | Catalog # |
|---|------------------|-----------|
| EasySep™ RBC Depletion Reagent | 100 mL | 18170 |
| EasySep™ RBC Depletion Reagent for RoboSep™ | 100 mL | 18170RF |

To learn more, visit www.stemcell.com/RBCdepletion.

HetaSep™

Erythrocyte Aggregation



HetaSep™

| Product Name | Catalog # | Size |
|--------------|-----------|--------|
| HetaSep™ | 07806 | 20 mL |
| | 07906 | 100 mL |

The presence of large numbers of RBCs in a colony-forming unit (CFU) assay prevents hematopoietic colonies from being accurately visualized either manually or using STEMvision™ (Figure 7). RBCs must be removed from fresh cord blood, bone marrow, and mobilized peripheral blood samples (whether whole or processed), before performing the CFU assay.

HetaSep™ is an erythrocyte aggregation agent used to quickly separate nucleated cells from RBCs. It is based on the principle that aggregated erythrocytes settle much faster than dispersed cells.

The HetaSep™ procedure does not affect the number of progenitor cells; 97% of CFUs are recovered in the RBC-cleared sample (Figure 8). HetaSep™-mediated RBC depletion requires only 50 µL of sample and is quick, making it easy to incorporate into an institution's workflow.

For more information, see the HetaSep™ Protocol Technical Bulletin (Document #29541) or visit www.stemcell.com/hetasep_protocol.

Why Use HetaSep™?

ACCURATE. Increases the accuracy of colony counting by removing RBC background.

CONSISTENT. Recover > 97% of colonies.

FAST. Easy to perform, no centrifuge needed. Can be performed with only 50 µL of sample.




Figure 7. STEMvision™ Images of 7-Day CFU Assays of Fresh Cord Blood Samples Plated in MethoCult™ Express

(A) Cord blood samples without RBC removal using HetaSep™. There is unacceptable background for a CFU assay. Note that fewer colonies are visible due to increased RBCs in culture. (B) Cord blood sample with prior removal of RBCs using HetaSep™. Acceptable background (minimal RBCs) for a CFU assay.




Figure 8. An Average of 97% of Colony-Forming Cells Are Recovered Following RBC Depletion with HetaSep™

Cord blood samples ($n = 10$) were split into two volumes, one of which was plated in a CFU assay without RBC depletion, and the other was plated after undergoing RBC depletion using the HetaSep™ protocol. Each sample type was plated in duplicate. CFU assays were counted manually and the percent recovery of colonies in each RBC-depleted fraction was calculated relative to results of CFU assays of non-depleted cells from the same donor.

Cryopreservation Media

cGMP Freezing and Preservation Media Formulated with USP-Grade Components

The cryopreservation and subsequent storage of hematopoietic cells is an important step in hematopoietic stem and progenitor cell research. STEMCELL Technologies' suite of cGMP, protein-free, and serum-free cryopreservation products are designed to maintain high viability and maximize cell recovery after long-term storage. For short-term storage and shipping of cells, HypoThermosol® products preserve cells at 2 - 8°C.



CryoStor® Freezing Media

- **Cell Types:** blood cells derived from peripheral blood.
- Designed to mitigate temperature-induced molecular stress responses during freezing and thawing.
- CryoStor® is pre-formulated with 2%, 5%, or 10% USP-grade DMSO.
- U.S. FDA Drug Master File



HypoThermosol® FRS Preservation Media

- **Cell Types:** all cells and tissues, including hematopoietic stem and progenitor cells.
- HypoThermosol® FRS is designed for short-term storage and/or shipment of cells at 2 - 8°C rather than at cryogenic temperatures.
- U.S. FDA Drug Master File



Video

Freeze Human PBMCs with CryoStor® CS10
www.stemcell.com/Cryopreserve-Primary-Cells-With-Cryostor



BloodStor® Freezing Media

- **Cell Types:** blood cells derived from peripheral blood, and bone marrow.
- BloodStor® 55-5 is preformulated with 55% (w/v) DMSO USP, 5% (w/v) Dextran-40 USP, and water for injection- (WFI) quality water.
- BloodStor® 100 contains 100% (w/v) DMSO USP.

Product Listing

| Product Name | Catalog # | Unit Size |
|--------------------|-----------|-------------------|
| CryoStor® CS10 | 100-1061 | 100 mL |
| | 07955 | 100 mL bag |
| | 07940 | 1000 mL bag |
| | 07931 | 5 x 16 mL vials |
| | 07959 | 5 x 10 mL vials |
| | 07952 | 16 x 10 mL vials |
| CryoStor® CS5 | 07933 | 100 mL |
| | 07953 | 100 mL bag |
| | 07949 | 5 x 10 mL vials |
| CryoStor® CS2 | 07932 | 100 mL |
| HypoThermosol® FRS | 07935 | 100 mL |
| | 07936 | 500 mL |
| | 07945 | 500 mL bag |
| | 07934 | 16 x 10 mL vials |
| BloodStor® 55-5 | 07937 | 16 x 7.2 mL vials |
| BloodStor® 100 | 07951 | 50 mL |
| | 07939 | 100 mL |
| | 07938 | 5 x 100 mL |

StemSpan™ Expansion Media & Supplements

Expansion and Differentiation of Hematopoietic Stem and Progenitor Cells

Cell Sourcing & Isolation

Expansion & Differentiation

Analysis

StemSpan™ Expansion Media include serum-free, xeno-free, and animal component-free formulations, as well as StemSpan™-AOF, the first commercially available animal origin-free medium for culturing hematopoietic stem and progenitor cells (HSPCs), manufactured under relevant cGMPs. StemSpan™ media do not contain any cytokines, allowing the flexibility to prepare a medium to meet the specific requirements of an experiment. StemSpan™ Expansion Supplements are pre-mixed cocktails of recombinant human cytokines and other additives formulated to selectively expand CD34+ HSPCs and/or stimulate their differentiation into mature cells of specific lineages when added to a StemSpan™ medium.

Applications

- Ex vivo expansion of HSPCs.^{1,2}
- Identification of novel regulators of HSPCs.³⁻⁷
- Production of large numbers of mature blood cells in vitro.^{3,8-9}
- Generation of target cells for reprogramming to make induced pluripotent stem cells.¹⁰
- Gene transfer into HSPCs.^{11,12}

Serum-Free Expansion Media

| Medium | Product Name | Catalog # (Size) | Recommended For | Components |
|-----------------------|--|----------------------------------|--|---|
| Serum-Free | StemSpan™ SFEM | 09600 (100 mL) 09650 (500 mL) | Culture of human HSPCs Culture of mouse, rat, and non-human primate HSPCs | Pre-tested BSA, insulin, transferrin, and supplements in IMDM |
| | StemSpan™ SFEM II | 09605 (100 mL) 09655 (500 mL) | Culture and expansion of human HSPCs Production of mature blood cells by expansion and lineage-specific differentiation of human HSPCs | |
| Animal Origin-Free | StemSpan™-AOF | 100-0130 (500 mL) | Culture and expansion of human HSPCs in the absence of human- or animal-derived components when higher compliance is required, e.g. for cell therapy manufacturing | Recombinant and synthetic components in IMDM |
| Xeno-Free | StemSpan™-XF | 100-0073 (500 mL) | Culture and expansion of human HSPCs in the absence of animal-derived components* | Human-derived or recombinant human proteins in IMDM |
| Animal Component-Free | StemSpan™-ACF Erythroid Expansion Medium | 09860 (500 mL) | Culture and expansion of human erythroid cells in the absence of human- or animal-derived components | Recombinant proteins and synthetic components only |

BSA: bovine serum albumin; IMDM: Iscove's Modified Dulbecco's Medium.

*Contains pre-tested human-derived and recombinant human proteins




Figure 9. StemSpan™ Media Support Greater Expansion of Human CD34+ and CD34^{bright} Cells Than Other Commercial Media

Purified CB-derived CD34+ cells were cultured for 7 days in select StemSpan™ media (StemSpan™ SFEM, StemSpan™ SFEM II, StemSpan™-XF, or StemSpan™-AOF, orange bars), and in five xeno-free media formulations from other suppliers (Commercial Alternative, gray bars) including (in random order) CTS™ StemPro™ HSC (Thermo), SCGM (Cellgenix), X-VIVO™ 15 (Lonza), Stemline™ II (Sigma), and StemPro™-34 (Thermo). All media were supplemented with StemSpan™ CD34+ Expansion Supplement and UM171*. Cell expansion of viable CD34+ and CD34^{bright} cells in culture were measured based on viable cell counts and flow cytometry results. StemSpan™ showed significantly higher expansion of CD34+ and CD34^{bright} cells ($P < 0.05$ when comparing StemSpan™ SFEM II to five media from other suppliers, calculated using a one-way ANOVA followed by Dunnett's post hoc test) and StemSpan™-AOF, the only animal origin-free formulation, showed equivalent performance to all xeno-free commercial alternatives tested. Data shown are mean \pm SEM ($n = 8$).

Note: Data for StemSpan™-AOF shown were generated with the original phenol red-containing version, StemSpan™-ACF (Catalog #09855). However, internal testing showed that the performance of the new phenol red-free, cGMP-manufactured version, StemSpan™-AOF (Catalog #100-0130), was comparable.

**Similar results are expected when using UM729 (Catalog #72332) prepared to a final concentration of 1 μ M. For more information, including data comparing UM171 and UM729, see Fares et al., 2014.*

StemSpan™-AOF

GMP-Compliant Medium for Cell Therapy Development

Hematopoietic stem and progenitor cells (HSPCs) are widely used in cell and gene therapy applications. When culturing HSPCs for cell therapy development, it is important to minimize risk and variability in your cell culture medium to ensure consistent, reproducible performance and safety. Whether you are performing fundamental research or ready to transition to the clinic, StemSpan™-AOF medium helps minimize the risk of viral contamination in your cell therapy development. Choosing animal origin-free (AOF) cell culture conditions can facilitate a smoother pathway to the clinic by helping you avoid regulatory roadblocks. StemSpan™-AOF is also manufactured under relevant cGMPs, ensuring the highest quality and consistency for reproducible results.

StemSpan™-AOF for Gene Editing Applications

The ability to genetically manipulate HSPCs has significantly advanced our understanding of the mechanisms that regulate hematopoiesis and is contributing to the development of novel cellular therapies. Using a medium that supports genome editing of hematopoietic cells can help take your cell therapy development to the next level. See how StemSpan™-AOF supports optimal culture conditions for HSPC maintenance and expansion in CRISPR-Cas9 gene editing applications:




Figure 10. StemSpan™-AOF Supports Equivalent or Greater Expansion of Human CD34⁺ and CD34^{bright} Cells Compared to Other Commercial Media

Purified cord blood-derived CD34⁺ cells were cultured for 7 days in StemSpan™-AOF (orange bar) and in four alternative commercial media (gray bars). Each medium was supplemented with StemSpan™ CD34⁺ Expansion Supplement and 175 nM UM171*. The (A) frequency and (B) cell expansion of viable CD34⁺ and CD34^{bright} cells in culture were measured based on viable cell counts and flow cytometry results. StemSpan™-AOF, the only animal origin-free formulation, showed equivalent performance to all xeno-free alternative media tested.

Why Use StemSpan™-AOF for Cell Therapy Development?

SAFE. Minimize the risk of viral contamination by using a medium that does not contain any primary or secondary raw materials derived from animals.

ROBUST. Ensure consistency in your experiments by using serum-free and animal origin-free culture conditions.

FLEXIBLE. Customize your cell culture conditions by adding StemSpan™ Expansion Supplements, individual cytokines, or additives to suit your specific cell therapy development needs.




Figure 11. Human CD34⁺ Cells Cultured in StemSpan™-AOF Show Equivalent Gene Editing Efficiency Compared with Alternative Media

Cells cultured for 2 days in either StemSpan™-AOF (orange bar) or xeno-free alternative commercial media (gray bars), each supplemented with StemSpan™ CD34⁺ Expansion Supplement and 175 nM UM171*, were electroporated with CRISPR-Cas9 RNP complexes containing crRNA:tracrRNA targeting $\beta 2$ Microglobulin (B2M). Non-electroporated (No EP) cells and cells electroporated with Cas9 without gRNA (Cas9 Only) were cultured in StemSpan™ SFEM II supplemented with StemSpan™ CD34⁺ Expansion Supplement plus 175 nM UM171*. B2M knockout efficiency (% MHC-I- viable cells) was monitored by flow cytometry using a fluorophore-conjugated anti-MHC-I antibody.

**Similar results are expected when using UM729 (Catalog #72332) prepared to a final concentration of 1 μ M. For more information, including data comparing UM171 and UM729, see Fares et al. Science, 2014.*

In Vivo Engraftment of HSPCs

One of the best assays to determine the quality of a hematopoietic cell therapy product is evaluation of its engraftment and multilineage differentiation potential after intravenous injection into immunodeficient mice (e.g. NOD scid gamma (NSG) mice). The “stemness” of HSPCs can be affected by many parameters, such as cell processing methods and culture conditions used for expansion and gene editing, which may impact the ability of the cells to successfully engraft. StemSpan™-AOF, the only animal origin-free cGMP medium on the market, supports multilineage engraftment of CD34+ cells at equivalent or higher levels when compared to uncultured cells.

See how StemSpan™-AOF supports the engraftment and expansion of cord blood-derived CD34+ cells in NSG mouse recipients:




Figure 12. StemSpan™-AOF-Expanded Cord Blood CD34+ Cells Engraft in NSG Mouse Recipients

Purified cord blood-derived CD34+ cells were cultured for 7 days in StemSpan™-AOF supplemented with StemSpan™ CD34+ Expansion Supplement and UM729 (1 μ M). After 7 days of expansion, progeny of 10,000 fresh or uncultured CD34+ cells were transplanted in sub-lethally irradiated NSG mice. (A-D) The number of human platelets and the frequency of human cells expressing the pan-leukocyte marker CD45 were measured in peripheral blood at 3 and 18 weeks post-transplantation. Data shown are mean \pm SEM (n = 3 - 5 mice). (A) At 3 weeks, engraftment of human platelets was lower in recipients of cells cultured in StemSpan™-AOF than in recipients of uncultured cells. (C) At week 18, there were no significant differences in platelet engraftment between the expanded and uncultured cells. (B,D) Human CD45+ cell frequencies in recipients of cells expanded in StemSpan™-AOF were similar to those in recipients of uncultured cells. (E-H) At week 20, long-term multilineage engraftment was measured in bone marrow of transplanted NSG mice. Data shown are mean \pm SEM (n = 3 - 4 mice). (E,F) Recipients of StemSpan™-AOF expanded cells showed similar frequencies of human CD45+ and CD34+ cells in the mouse bone marrow compared to recipients of uncultured cells. (G,H) Cells expanded with StemSpan™-AOF showed similar levels of myeloid (CD45+CD33+) and lymphoid (CD45+19+ B cells and CD45+CD3+ T cells) engraftment relative to uncultured cells.

Taking Your Research to the Clinic?

STEMCELL's Services for Cell Therapy program has a team of experts who can help support your regulatory filing by providing custom solutions such as quality documentation, additional product testing, and customized product manufacturing. To learn more about how we can support your preclinical and clinical research needs, visit us at www.stemcell.com/services/cell-therapy.html.



Resource

Visit the Cell Therapy Learning Center
www.stemcell.com/hsc-learning

CD34 Expansion Supplements

Investigators studying hematopoiesis require standardized culture media and cytokines to promote the proliferation and/or lineage-specific differentiation of HSPCs from human bone marrow (BM), cord blood (CB), and other tissues. StemSpan™ media require the addition of cytokines to promote HSPC proliferation and differentiation. The choice of cytokines and supplements depends on the objective of the experiment, i.e. on the desired numbers of specific cell types to be generated in vitro.




Figure 13. StemSpan™ SFEM II Serum-Free Expansion Medium Containing CC100 Cytokine Cocktail Supports Greater Expansion of Human CD34+ Cells than Other Media Tested

Expansion of CD34+ cells, normalized relative to the values obtained in StemSpan™ SFEM medium (dark gray bars), after culturing purified CD34+ CB (A, n = 6) or bone marrow (BM) (B, n = 3) cells for 7 days in StemSpan™ SFEM, SFEM II (blue bars) and AOF (orange bars), and six media from other commercial suppliers (light gray bars). All media were supplemented with StemSpan™ CC100 Cytokine Cocktail (Catalog #02690). Vertical lines indicate 95% confidence limits, i.e. the range within which 95% of results fall. The numbers of CB and BM cells produced in StemSpan™ SFEM II were significantly higher than in all other media, with the exception of CB cells cultured in StemSpan™-AOF (*p<0.05, paired t-test).

Note: Data for StemSpan™-AOF shown were generated with the original phenol red-containing version, StemSpan™-ACF (Catalog #09855). However, internal testing showed that the performance of the new phenol red-free, cGMP version, StemSpan™-AOF (Catalog #100-0130), was comparable.

CD34+ Cell Expansion Supplements

| Product Name | Catalog # (Size) | Product Features | Components |
|--------------------------------------|--------------------|---|--|
| StemSpan™ CC100 | 02690 (1 mL, 100X) | Contains both early- and late-acting cytokines Stimulates the production of large numbers of human hematopoietic cells, including CD34+ progenitor cells | rhFlt3L, rhSCF, rhIL-3, rhIL-6 |
| StemSpan™ CC110 | 02697 (1 mL, 100X) | Contains early-acting cytokines Stimulates similar expansion of CD34+ cells as CC100, but with higher purity | rhFlt3L, rhSCF, rhTPO |
| StemSpan™ CD34+ Expansion Supplement | 02691(10 mL, 10X) | Recommended for selective expansion of human CD34+ HSPCs Stimulates greater CD34+ cell expansion compared to CC100 and CC110 | rhFlt3L, rhSCF, rhIL-6, rhTPO Other additives |

rh: recombinant human. For a full listing of products for expansion and differentiation of HSPCs, visit www.stemcell.com/HSPCworkflow.

Expansion and Lineage-Specific Differentiation Supplements and Kits

| Cell Type | Product Name | Catalog # (Size) | Product Features | Components |
|----------------------------|--|----------------------|--|---|
| Erythroid Progenitor Cells | StemSpan™ Erythroid Expansion Supplement | 02692 (1 mL, 100X) | Stimulates the production of human erythroid cells by expansion and lineage-specific differentiation of human HPCs | rhSCF, rhIL-3, rhEPO |
| Megakaryocytes | StemSpan™ Megakaryocyte Expansion Supplement | 02696 (1 mL, 100X) | Stimulates the production of human megakaryocytes by expansion and lineage-specific differentiation of human HPCs | rhSCF, rhTPO, rhIL-6, rhIL-9 |
| Granulocytes | StemSpan™ Myeloid Expansion Supplement | 02693 (1 mL, 100X) | Stimulates the production of human granulocytes by expansion and lineage-specific differentiation of human HPCs | rhSCF, rhTPO, rhG-CSF, rhGM-CSF |
| Monocytes | StemSpan™ Myeloid Expansion Supplement II | 02694 (1 mL, 100X) | Stimulates the production of human monocytes by expansion and lineage-specific differentiation of human HPCs | rhFlt3L, rhSCF, rhTPO, rhM-CSF, rhGM-CSF and supplements |
| T Cells | StemSpan™ T Cell Generation Kit | 09940 (Kit, various) | Stimulates the production of human T cells by expansion and lineage-specific differentiation of human HSPCs in stroma-free conditions | <ul style="list-style-type: none"> • SFEM II • Lymphoid Progenitor Expansion Supplement (10X) • Lymphoid Differentiation Coating Material (100X) • T Cell Progenitor Maturation Supplement (10X) |
| NK Cells | StemSpan™ NK Cell Generation Kit | 09960 | Stimulates the production of human NK cells by expansion and differentiation of human CD34+ HSPCs in stroma-free conditions | <ul style="list-style-type: none"> • SFEM II • Lymphoid Progenitor Expansion Supplement (10X) • Lymphoid Differentiation Coating Material (100X) • NK Cell Differentiation Supplement (100X) |
| B Cells | StemSpan™ B Cell Generation Kit | 100-1250 | Stimulates the production of B cells and antibody-secreting cells (ASCs) of human CD34+ HSPCs in serum- and feeder-free conditions | <ul style="list-style-type: none"> • StemSpan™ SFEM II • StemSpan™ B Cell Differentiation Supplement 1 (20X) • StemSpan™ B Cell Differentiation Supplement 2 (20X) • StemSpan™ B Cell Differentiation Supplement 3 (20X) • StemSpan™ B Cell Differentiation Supplement 4 (20X) |
| Myeloid Leukemic Cells | StemSpan™ Leukemic Cell Culture Kit | 09720 (Kit, various) | Stimulates the production of human myeloid leukemia cells, allowing users to expand, culture, and use malignant cells for drug screening | <ul style="list-style-type: none"> • SFEM II • CD34+ Expansion Supplement • UM729 |

rh: recombinant human. For a full listing of products for expansion and differentiation of HSPCs, visit www.stemcell.com/HSPCworkflow.

Differentiation of Human Progenitor Cells in Culture




Figure 14. Production of Erythroid Cells from Human CB-Derived CD34+ Cells Cultured in StemSpan™ Media Containing StemSpan™ Erythroid Expansion Supplement

Average numbers of erythroid cells generated after culturing purified CD34+ CB cells (n = 5) for 14 days in StemSpan™ SFEM gray bar), SFEM II (blue bar), or StemSpan™-ACF Erythroid Expansion Medium (ACF-E, orange bar) media containing StemSpan™ Erythroid Expansion Supplement (Catalog #02692). Shown are the number of erythroid cells that express CD71 and/or Glycophorin A (GlyA) per input CD34+ cell.




Figure 15. Comparison of Megakaryocyte Expansion in Different StemSpan™ Media Containing StemSpan™ Megakaryocyte Expansion Supplement

Average numbers of CD41+ megakaryocytic cells normalized relative to the values obtained in StemSpan™ SFEM (gray bar) after culturing purified CD34+ cord blood cells (n=6) for 14 days in StemSpan™ SFEM, SFEM II (blue bar), and AOF (orange bar) media containing StemSpan™ Megakaryocyte Expansion Supplement. Vertical lines indicate 95% confidence limits, the range within which 95% of results typically fall.

*The numbers of CD41a+ cells were significantly higher in SFEM II (p<0.01, paired t-test, n=6) compared to SFEM and AOF medium.

Note: Data for StemSpan™-AOF shown were generated with the original phenol red-containing version StemSpan™-ACF (Catalog #09855). However internal testing showed that the performance of the new phenol red-free, cGMP-manufactured version, StemSpan™-AOF (Catalog #100-0130) was comparable.

Lymphoid Expansion and Differentiation




Figure 16. Frequency and Yield of CD4 Immature Single-Positive and CD4+CD8+ Double-Positive Cells After 42 Days of Culture

CB-derived CD34+ cells were cultured with the StemSpan™ T Cell Generation Kit (Catalog #09940) for 42 days. The average (A) frequency and (B) yield of CD4 immature single-positive (ISP), CD4+CD8+ double-positive, and CD3+TCRαβ+-expressing double-positive cells (CD4+CD8+CD3+TCRαβ+) are shown. Vertical lines indicate 95% confidence interval (n = 31). On average, 38% of the total viable population were CD4+CD8+ double-positive, of which 35% co-expressed CD3 and TCRαβ. The yield of total double-positive cells per input CD34+ cell was ~23,000.




Figure 17. Frequency and Yield of CD8 Single-Positive T Cells After 49 Days of Culture

CD4+CD8+ double-positive cells were further matured into CD8 single-positive T cells by culturing for an additional 7 days in StemSpan™ SFEM II with T Cell Progenitor Maturation Supplement (Catalog #09930), IL-15 (Catalog #78031), and ImmunoCult™ CD3/CD28/CD2 T Cell Activator (Catalog #10970). The average (A) frequency and (B) yield of CD3+TCRαβ+-expressing cells and their subsets are shown. Vertical lines indicate 95% confidence interval (n = 12). On average, 54% of the CD3+TCRαβ+ cells were CD4+CD8+ double-positive and 38% were CD8 single-positive (CD4-CD8+).

Small Molecules for Human Hematopoietic Stem and Progenitor Cell Research

UM729 and StemRegenin 1

UM729 and StemRegenin 1 (SR1) are two small molecules that have been found to enhance the self-renewal and expansion of human hematopoietic stem and progenitor cells in vitro.¹³⁻¹⁵ UM729 is a pyrimido-[4,5-b]-indole derivative that acts differently than other small molecule stimulators of hematopoiesis, such as the aryl hydrocarbon receptor (AhR) antagonist SR1.¹³⁻¹⁵ UM729 was originally discovered in a screen of compounds capable of promoting human CD34+ cell expansion and later underwent structure-activity relationship optimization to develop UM171.¹⁴⁻¹⁵ Additionally, UM729 and SR1 have been shown to cooperate in culture, resulting in an additive effect in preventing the differentiation of primary human acute myeloid leukemia (AML) cells.¹⁵




Figure 18. Chemical Structure of UM729 and StemRegenin 1

(A) UM729: Methyl 4-((3-(piperidin-1-yl)propyl)amino)-9H-pyrimido[4,5-b] indole-7-carboxylate

(B) StemRegenin 1 (SR1): 4-[2-[[2-benzo[b]thien-3-yl]-9-(1-methylethyl)-9H-purin-6-yl]amino]ethyl]-phenol

Product Listing

| Product Name | Catalog # | Size | Pathway / Target | Applications |
|-------------------------------|----------------|--------------|--|---|
| UM729 | 72332 | 250 µg | Pyrimido-indole derivative that enhances HSC self-renewal in vitro | <ul style="list-style-type: none"> Expansion of human HSCs in culture Maintenance of LSCs in combination with SR1 |
| StemRegenin 1 | 72342 72344 | 1 mg 5 mg | Aryl hydrocarbon receptor (AhR) antagonist | <ul style="list-style-type: none"> Maintenance and expansion of human HSPCs in culture Differentiation of human CD34+ HPCs into functional dendritic cells Maintenance of LSC activity in culture when combined with UM729 |
| StemRegenin 1 (Hydrochloride) | 72352 72354 | 1 mg 5 mg | | |
| StemSpan™ HSC Plus Supplement | 100-1694 | 1 mL | Small molecule cocktail that enhances HSPC expansion in vitro | <ul style="list-style-type: none"> Selective expansion of primitive HSPCs, including CD34+CD45RA-CD90+CD133+EPCR+ cells in liquid cultures initiated with CD34+ cells or purified CD34+ subsets isolated from cord blood, mobilized peripheral blood, or bone marrow |

HPCs: hematopoietic progenitor cells; HSCs: hematopoietic stem cells; HSPCs: hematopoietic stem and progenitor cells; LSCs: leukemic stem cells.

StemSpan™ HSC Plus Expansion Supplement

StemSpan™ HSC Plus Supplement contains a combination of small molecules formulated to selectively promote the expansion of hematopoietic stem and progenitor cells (HSPCs) in cultures of human CD34+ cells or purified CD34+ subsets enriched for hematopoietic stem cells (HSCs), when added to the culture medium containing appropriate cytokines.




Figure 19. Expansion of Primitive HSPCs Is Enhanced by the Addition of StemSpan™ HSC Plus Supplement to Cytokine-Containing Medium

Purified CD34+ cells derived from cord blood (CB) were cultured at a concentration of 10,000 cells per mL in StemSpan™ SFEM II medium supplemented with one of the StemSpan™ cytokine-based supplements alone (CD34 Expansion Supplement, CC100, or CC110; gray bars), or together with StemSpan™ HSC Plus supplement (orange bars). After 7 days of culture, flow cytometry was used to analyze the (A) total nucleated cell (TNC) expansion, (B) expansion of CD34+ cell subsets, and (C) frequency of CD34+ HSPC subsets: CD34+, CD34+CD45RA-CD90+, and CD34+CD45RA-CD90+EPCR+CD133+. StemSpan™ HSC Plus Supplement promoted greater expansion of CD34+ CD45RA-CD90+ and CD34+CD45RA-CD90+CD133+EPCR+ HSPC subsets compared to cytokine-based supplements alone. Each subsequent cell subset was progressively more enriched for phenotypic stem/progenitor cells. Data shown are \pm SEM ($n = 24$ for the CD34 Expansion Supplement dataset and $n = 9$ for both the CC100 and CC110 datasets).

STEMdiff™ Kits for Hematopoietic Differentiation of Pluripotent Stem Cells

Generate Hematopoietic Progenitor Cells, Immune Cells, and Blood Cells from Human Embryonic Stem and Induced Pluripotent Stem Cells

The STEMdiff™ Hematopoietic Kit consists of serum-free basal medium and supplements designed for the generation of hematopoietic progenitor cells (HPCs). Optimized for a standardized, 12-day differentiation protocol, this kit supports robust differentiation of human pluripotent stem cells (hPSCs), including human embryonic stem cells (hESCs) and induced pluripotent stem cells (hiPSCs). HPCs generated in these cultures can be identified by the expression of CD34 and CD45, and by the ability to form hematopoietic colonies of multiple lineages in colony-forming unit (CFU) assays with MethoCult™ medium.

The resulting HPCs may be used for downstream assays or quantified in a CFU assay with MethoCult™ SF H4636 (Catalog # 04636) medium, designed specifically for use with hPSC-derived HPCs, or MethoCult™ H4435 (Catalog #04435) Enriched medium. HPCs generated with the STEMdiff™ Hematopoietic Kit may be further differentiated using the STEMdiff™ Microglia Differentiation Kit or STEMdiff™ Monocyte Kit. HPCs and downstream cells in the erythroid or megakaryocyte lineages may be obtained directly using the STEMdiff™ Erythroid or Megakaryocyte Kits, and HPC and immune cell types in the lymphoid lineages may be obtained using the STEMdiff™ NK and T Cell Kits.




Figure 20. Efficient and Robust Generation of CD34+CD45+ HPCs

hESCs and hiPSCs were cultured for 12 days in single wells of 12-well plates using the STEMdiff™ Hematopoietic Kit. At the end of the culture period, cells in suspension were harvested, stained, and analyzed by flow cytometry for the expression of hematopoietic cell surface markers CD34 and CD45. (A) Percentages and (B) total numbers of CD34+CD45+ cells in cultures of hESCs or hiPSCs are shown for 6 cell lines. Data shown as mean + SEM; n ≥ 3.

Why Use the STEMdiff™ Hematopoietic Kit?

CONSISTENT. Reduce variability with a serum- and feeder-free formulation.

EASY-TO-USE. Produce HPCs in suspension for easy harvest with a simple monolayer protocol.

RAPID. Generate HPCs in 12 days.

HIGH YIELD. Generate 4 - 18 million CD34+CD45+ HPCs with one kit.

FLEXIBLE. Generate HPCs across multiple hESC and hiPSC lines.




Figure 21. Schematic for Hematopoietic Differentiation

One day before the differentiation protocol, hPSC colonies are harvested and seeded as small aggregates (100 - 200 µm in diameter) at 10 - 20 aggregates/cm² in a TeSR™ maintenance medium (mTeSR™ Plus, mTeSR™1, or TeSR™-E8™). After one day, TeSR™ medium is replaced with Medium A to begin inducing the cells toward a mesoderm-like state (Day 0). On Day 2, a half-medium change is performed with fresh Medium A. On Day 3, the medium is changed to Medium B with half-medium changes on Days 5, 7, and 10 to promote further hematopoietic differentiation. Typically, by Day 12, large numbers of HPCs can be harvested from the culture supernatant.

Learn more at www.stemcell.com/STEMdiffHeme

STEMdiff™ Megakaryocyte Kit

STEMdiff™ Megakaryocyte Kit is designed for the serum-free and feeder-free differentiation of hESCs and hiPSCs to megakaryocytes expressing CD41a and CD42b. This optimized two-dimensional and two-stage protocol is capable of generating high yields of megakaryocytes per hESC or hiPSC in 17 days. The resulting megakaryocytes show high ploidy and platelet-shedding ability and are also amenable to large-scale culture.

Why Use STEMdiff™ Megakaryocyte Kit?

CONSISTENT. Reduce variability with a serum- and feeder-free formulation.

HIGH YIELD. Expand hESCs and hiPSCs by > 400-fold during megakaryocyte differentiations, with > 70% yield for CD41a and CD42b megakaryocytes.

EASY-TO-USE. Optimized two-dimensional protocol produces megakaryocytes in 17 days.




Figure 22. hESC- and hiPSC-Derived Hematopoietic Progenitor Cells Efficiently Expand and Differentiate to CD41a+CD42b+ Megakaryocytes

Frequencies and numbers of CD41a+CD42b+ MKs per input cell for two hESC lines (H1 and H9) and two hiPSC lines (WLS-1C and STiPS-R038). The average frequency of viable CD41a+CD42b+ cells on day 17 ranged between 56% and 77%. The average yield of CD41a+CD42b+ MKs generated per input cell ranged between 223 and 425. Data are shown as mean ± SEM (n = 12 for H1, n = 29 for H9, n = 27 for WLS-1C, n = 12 for STiPS-R038).

Learn more at www.stemcell.com/megakaryocyte-diff

STEMdiff™ Erythroid Kit

STEMdiff™ Erythroid Kit is designed for the serum-free and feeder-free differentiation of hESCs and hiPSCs to erythroid cells expressing CD71 and CD235ab. This optimized two-dimensional and two-stage protocol supports the generation of high yields of erythroid cells in 24 days.

Differentiate hPSCs to erythroid cells (erythroblasts) expressing Glycophorin A and CD71. hESCs or hiPSCs are induced toward erythroid-biased hematopoietic progenitor cells and then further differentiated to erythroid progenitor cells (Day 10 - 24). Cells generated using the STEMdiff™ Erythroid Kit can be further matured into normoblasts and reticulocytes once moved to appropriate culture conditions for maturation.

Why Use the STEMdiff™ Erythroid Kit?

CONSISTENT. Differentiate to erythroblasts in two stages, without feeders or serum.

HIGH YIELD. Expand hematopoietic progenitor cells by > 200-fold during erythroid differentiation, with > 70% yield for CD71+GlyA+ erythroblasts.

VERSATILE. hESC- and hiPSC-derived erythroblasts can be matured further using StemSpan™ SFEM II medium with appropriate supplements.




Figure 23. hESC- and hiPSC-Derived Erythroid Cells Are Hemoglobinized and Display Typical Erythroid Morphology

(A) Erythroid cells generated with the STEMdiff™ Erythroid Kit express a mix of primitive (embryonic) and definitive (fetal, adult) hemoglobin. Shown are the results of qPCR analysis for globin gene expression after 24 days of culture. (B) A picture of the cell pellet shows that cells produced in culture are hemoglobinized. (C) Cells display typical basophilic erythroblast morphology after 24 days of culture using the STEMdiff™ Erythroid Kit (40X magnification; May-Grunwald Giemsa stain).

Learn more at www.stemcell.com/erythro-diff

Recombinant Cytokines

Recombinant Human Cytokines

| Cytokine | Catalog # | Unit Size (µg) |
|-----------------------------------|-----------|----------------|
| bFGF* | 78003 | 50 |
| BMP-2* | 78004 | 50 |
| BMP-4 | 78211 | 20 |
| EGF* | 78006 | 500 |
| EPO | 78007 | 50 |
| EPO (HEK293-expressed) | 100-1716 | 100 |
| | 100-1717 | 1000 |
| FGF-4 | 78103 | 25 |
| FGF-7 (KGF)* | 78046 | 10 |
| FGF-8B* | 78008.1 | 50 |
| FGF-10 (KGF-2)* | 78037 | 10 |
| FGF-18 | 78041 | 10 |
| Flt3/Flk-2 Ligand* | 78009 | 100 |
| Flt3/Flk-2 Ligand (CHO-expressed) | 100-1710 | 50 |
| | 100-1711 | 1000 |
| G-CSF* | 78012 | 100 |
| GM-CSF (CHO-expressed) | 78190 | 10 |
| GM-CSF (E. coli-expressed) | 78015 | 100 |
| IFN-γ* | 78020 | 100 |
| IGF-I* | 78022 | 500 |
| IGF-II | 78023 | 50 |
| IL-1β* | 78034 | 100 |
| IL-2 (CHO-expressed)* | 78036 | 50 |
| IL-2 (E. coli-expressed)* | 78220 | 10 |
| IL-3 (CHO-expressed)* | 78194 | 10 |
| IL-3 (E. coli-expressed)* | 78040 | 100 |
| IL-4* | 78045 | 100 |
| IL-5 | 78048.1 | 10 |
| IL-6 | 78050 | 100 |
| IL-6Ra | 78083.1 | 50 |
| IL-7* | 78053 | 100 |
| IL-10* | 78024 | 50 |
| IL-11 | 78025 | 100 |
| IL-12 | 78027 | 25 |
| IL-13 | 78029 | 100 |
| IL-15* | 78031 | 100 |
| IL-21* | 78082 | 10 |
| IL-31 | 78216 | 25 |
| M-CSF* | 78057 | 100 |
| MIP-1β (CCL4) | 78090 | 5 |
| Oncostatin M* | 78094 | 10 |
| PDGF-AB | 78096 | 10 |
| PDGF-BB* | 78097 | 10 |
| PDGF-CC | 78168 | 10 |
| PDGF-DD | 78222 | 25 |
| R-Spondin-1 | 78213 | 25 |
| SCF* | 78062 | 100 |
| SCF (P. pastoris-expressed) | 100-1712 | 50 |
| | 100-1713 | 1000 |

| Cytokine | Catalog # | Unit Size (µg) |
|-----------|-----------|----------------|
| TGF-β1 | 78067 | 5 |
| TNF-α* | 78068 | 50 |
| TPO | 78210 | 25 |
| VEGF-121 | 78127 | 10 |
| VEGF-165* | 78073 | 50 |

*Animal Component-Free (ACF) version available.

Recombinant Mouse Cytokines

| Cytokine | Catalog # | Unit Size (µg) |
|-----------------------------|-----------|----------------|
| BMP-2, ACF | 78135 | 10 |
| EGF | 78016 | 10 |
| FGF-8B, ACF | 78204 | 25 |
| FGF-21 | 78108.1 | 50 |
| Flt3/Flk-2 Ligand | 78011 | 100 |
| G-CSF | 78014 | 100 |
| GM-CSF (CHO-expressed) | 78206 | 10 |
| IFN-γ | 78021 | 100 |
| IL-1β | 78035 | 50 |
| IL-2 | 78081 | 20 |
| IL-3 | 78042 | 100 |
| IL-4 | 78047 | 100 |
| IL-5 | 78049.1 | 25 |
| IL-6 | 78052 | 100 |
| IL-7 | 78054 | 50 |
| IL-10 | 78079 | 10 |
| IL-11 | 78026 | 100 |
| IL-12 | 78028.1 | 10 |
| IL-13 | 78030.1 | 10 |
| M-CSF | 78059 | 100 |
| MIP-1α (CCL3) | 78089 | 10 |
| MIP-1β (CCL4) | 78091 | 10 |
| PDGF-BB | 78178 | 10 |
| SCF | 78064 | 100 |
| SCF (P. pastoris-expressed) | 100-1714 | 50 |
| | 100-1715 | 1000 |
| SDF-1α (CXCL12) | 78121 | 5 |
| TNF-α | 78069 | 100 |
| TPO | 78072 | 50 |
| VEGF-164 | 78102 | 20 |

These high-quality cytokines ensure reproducibility across a variety of hematopoietic cell culture applications. Choose from a wide selection of cytokines, also available in other sizes, to incorporate into your research workflow. For a complete listing of cytokines, visit www.stemcell.com/cytokines.

MyeloCult™ Long-Term Culture Media

Detection of Primitive Progenitor Cells

Cell Sourcing & Isolation

Expansion & Differentiation

Analysis

Long-Term Culture Media

The long-term culture system for hematopoietic stem and progenitor cells (HSPCs), developed first for mouse marrow in the late 1970s, and then successfully adapted for human cells, establishes the essential cell types involved in hematopoiesis, *in vitro*.¹⁶⁻¹⁷ When initiated with a relatively high density of bone marrow cells ($> 10^6$ cells/mL), long-term cultures are characterized by the formation of an adherent stromal layer of mesenchymal cells, including endothelial cells, fibroblasts, and adipocytes. Primitive hematopoietic cells associated with this stromal layer typically generate myeloid clonogenic progenitor cells and mature granulocytes for many weeks, provided that appropriate medium and supplements, incubation conditions, and feeding schedules are used.

MyeloCult™ is a unique long-term culture medium that promotes the formation of human or mouse primary stromal layers and allows the proliferation and differentiation of primitive hematopoietic progenitor cells.

Long-Term Culture-Initiating Cell Assay

The unique features of the long-term culture system have allowed the development of the long-term culture-initiating cell (LTC-IC) assay to detect and quantitate primitive hematopoietic cells, which share phenotypic and functional properties with mouse or human *in vivo* repopulating cells.¹⁸⁻²⁰ In human long-term cultures, colony-forming units (CFUs) detected after more than 5 weeks represent the progeny of LTC-ICs, as CFUs present in the input cell suspension will have undergone terminal differentiation by this time.

Quantitation of LTC-ICs in a test cell suspension requires culturing the cells on a supportive feeder layer of irradiated marrow cells or suitable human or mouse fibroblast cell lines.²¹⁻²² Limiting dilution analysis is used to determine the frequency of LTC-ICs, as well as the average number of CFUs produced per LTC-IC. Once the average number of CFUs per LTC-IC is established, the LTC-IC content of a sample can be determined by a bulk culture LTC-IC assay, provided that the same source of test cells (e.g. bone marrow, mobilized peripheral blood, or cord blood) is used and the assay conditions are identical. The LTC-IC content is then calculated by dividing the total output of CFUs by the average number of CFUs produced per LTC-IC.¹⁸




Figure 24. Bulk Culture LTC-IC Assay Procedure

Published Applications

Quantitate Frequency of LTC-ICs and Study Their Phenotypic and Functional Properties. Experiments using MyeloCult™ and related products show that LTC-ICs are a heterogeneous population of cells that can differ in phenotype, cell cycling characteristics, and expansion potential.²⁴⁻²⁹

Facilitate Gene Transfer. Culturing primitive hematopoietic progenitor cells in MyeloCult™ facilitates retroviral gene transfer and expansion of these transfected cells.^{30,31}

Expand Multi-Potential Hematopoietic Cells In Vitro. Human colony-forming units (CFUs) and LTC-ICs have been expanded using MyeloCult™ in stirred suspension cultures.^{32,33} Mouse totipotent hematopoietic cells expanded in static long-term cultures with MyeloCult™ can sustain lymphomyelopoiesis in irradiated recipients.³⁰

Evaluate Factors Regulating Myelopoiesis. The role of stroma-derived factors (positive and/or negative regulators and adhesion molecules) in the regulation of myelopoiesis can be evaluated in long-term cultures using MyeloCult™.³⁴⁻³⁷

Examine Differences Between Normal and Malignant Cells. MyeloCult™ has been used to culture LTC-ICs from patients with chronic myeloid leukemia, acute myeloid leukemia, and aplastic anemia.³⁸⁻⁴²

Study Differentiation of CD34+ Cells into Natural Killer (NK) Cells. In the presence of IL-2 and IL-7, or SCF and IL-15, a subpopulation of CD34+ cord blood cells can be induced to differentiate into NK cells when cultured in MyeloCult™.⁴³

Culture NK Cell Lines. MyeloCult™ has been used for the culture of the human NK-92 cell line.⁴⁴

MyeloCult™ Media for Initiation and Maintenance of Myeloid Long-Term Cultures

| Product Name | Catalog # | Size | Description |
|------------------|-----------|--------|--|
| MyeloCult™ H5100 | 05150 | 500 mL | Myeloid long-term cultures of human hematopoietic progenitor cells and stromal feeder layers |

MethoCult™ Media to Evaluate the Number of CFU per LTC-IC

| Product Name | Catalog # | Size | Description |
|---------------------------|----------------|---------------------|--|
| MethoCult™ H4435 Enriched | 04435 04445 | 100 mL 24 x 3 mL | Detection of human colonies derived from LTC-ICs |
| MethoCult™ GF M3434 | 03434 03444 | 100 mL 24 x 3 mL | Detection of mouse colonies derived from LTC-ICs |

Support Products

| Product Name | Catalog # | Size | Description |
|------------------|-----------|--------|--|
| Hydrocortisone | 07904 | 100 mg | Supplementation of MyeloCult™ media to a final concentration of 10 ⁻⁶ M. Suitable for supplementing both human (Catalog #05150) and mouse (Catalog #05350) MyeloCult™ media for LTC-ICs assays. |
| L-Calc™ software | 28600 | N/A | Limiting dilution analyses for determination of frequencies and other applications NOTE: Free download of L-Calc™ Software is available at www.stemcell.com/l-calc-software.html |

MethoCult™ Media

For Performing Colony-Forming Unit Assays

The colony-forming unit (CFU) assay is an in vitro functional assay for counting multipotential and lineage-committed hematopoietic progenitor cells (HPCs) in bone marrow, blood, and other hematopoietic tissues. MethoCult™ is a line of methylcellulose-based media formulated to promote optimal growth and differentiation of hematopoietic progenitor cells from various species (humans, non-human primates, mice, rats, and dogs). MethoCult™ is widely recognized as the “Gold Standard” for detection and quantification of hematopoietic progenitor cells in the CFU assay.

Key features of the CFU assay:

- Individual HPCs proliferate and differentiate to produce colonies of mature blood cells during culture in MethoCult™ medium containing growth factors and supplements.
- Progenitor cells of different lineages and stages of maturation produce colonies that differ in their size, morphology, and cellular composition.
- Each colony is derived from a single progenitor cell or CFU. The number of colonies provides a measure of the number of viable and functional CFUs in the cell sample being tested (i.e. 1 colony = 1 CFU).
- Enumerates all classes of myeloid and/or erythroid progenitor cells: erythroid (BFU-E and CFU-E), granulocyte-macrophage (CFU-GM, CFU-G, and CFU-M), and multi-potential progenitor cells (CFU-GEMM).

Why Use MethoCult™?

ROBUST. Ensure consistency and reproducibility in your experiments with standardized and rigorously tested media.

CONVENIENT. Identify and count total CFUs, erythroid (CFU-E and BFU-E), granulocyte/macrophage (CFU-GM, CFU-G, and CFU-M), and multi-lineage (CFU-GEMM) progenitor cells with ready-to-use formulations.

FLEXIBLE. Customize your cell-culture conditions by adding components of choice to suit your specific research needs.

BFU-E: burst-forming unit – erythroid; CFU-E: colony-forming unit – erythroid; CFU-GM: colony-forming unit – granulocyte/macrophage; CFU-G: colony-forming unit – granulocyte; CFU-M: colony-forming unit – macrophage; CFU-GEMM: colony-forming unit – granulocyte, erythrocyte, macrophage, megakaryocyte.

Applications

- Quantitation and characterization of human hematopoietic progenitor cells from cord blood, mobilized peripheral blood, and bone marrow⁴⁵
- Quantitation and characterization of hematopoietic progenitor cells from mouse bone marrow and other cell samples
- Quantitation of primitive hematopoietic progenitor cells from human and mouse long-term culture-initiating cell (LTC-IC) assays^{19,21}
- Evaluation of hematopoietic cell differentiation from human pluripotent stem cells (hPSCs)
- Quality control of cryopreservation, cell processing, and ex vivo manipulation procedures⁴⁶⁻⁵⁴
- Support of patient diagnosis, prognosis, and treatment in a clinical hematology lab⁵⁵⁻⁶⁰
- Support of the evaluation of donor samples, including cord blood, for stem cell transplants⁶¹⁻⁶⁶
- Study of the effects of cytokines, growth factors, hormones, or mimetics on hematopoietic progenitor cells⁶⁷⁻⁷¹
- Toxicity testing or drug screening assays⁷²⁻⁷⁵
- Optimization of gene transfer protocols and performing replating assays to study myeloid progenitor cell proliferation and self-renewal following genetic manipulation⁷⁶
- Quantitation of hematopoietic progenitor cells following ex vivo expansion^{27,68}



MethoCult™ H4034 Optimum Media in 24 x 3 mL and 100 mL Formats

The Colony Forming Unit (CFU) Assay

CFU Assay Workflow

Two of the most commonly performed assays for measuring the functionality of HSPCs are the in vivo mouse model engraftment study and the in vitro colony-forming unit (CFU) assay. Routine in vivo studies for potency assays are time-consuming, expensive, and difficult to validate. The CFU assay is significantly faster and more cost-effective compared to the in vivo mouse engraftment assay and has been shown to correlate well with in vivo engraftment in patients^{62, 65-66}; it is thus gaining popularity as the method of choice for evaluating the quality and potency of HSPCs used in cell therapy. There are currently no specific recommendations from regulatory agencies regarding the potency testing of hematopoietic cell therapy products. However, several cord blood bank accrediting bodies (including AABB and FACT) now require pre- and post-processing CFU data from cord blood collection as evidence of sample stability and to confirm the minimal loss of HSPC functionality during processing and cryopreservation. These requirements highlight the growing importance of incorporating the CFU assay in clinical workflows.

How Does the CFU Assay Work?




Figure 25. Colony-Forming Unit (CFU) Assay Procedure




Figure 26. Images of Human and Mouse Colonies Visualized on an Inverted Microscope



On-Demand Training

Access Free Instruction on CFU Set Up, Colony Identification, and Scoring
www.stemcell.com/hsc-training

MethoCult™ Media for Human Cells

| MethoCult™ Product | Catalog # | Unit Size | Components | | | | | Applications |
|--|-------------|---------------------|------------|-----|-----|-----------------------|--|--|
| | | | MC | FBS | BSA | Insulin + Transferrin | Growth Factors | |
| H4034 Optimum (GF H4034) | 04034 04044 | 100 mL 24 x 3 mL | • | • | • | | rhSCF, rhIL-3, rhG-CSF, rhEPO, rhGM-CSF | <ul style="list-style-type: none"> Detection of CFU-E, BFU-E, CFU-GM, CFU-GEMM in BM, MPB, PB, and CB Compatible with STEMvision™ |
| H4035 Optimum Without EPO (GF H4035) | 04035 04045 | 100 mL 24 x 3 mL | • | • | • | | rhSCF, rhIL-3, rhG-CSF, rhGM-CSF; no rhEPO | <ul style="list-style-type: none"> Detection of CFU-GM in BM, MPB, PB, and CB Compatible with STEMvision™ |
| H4434 Classic (GF H4434) | 04434 04444 | 100 mL 24 x 3 mL | • | • | • | | rhSCF, rhIL-3, rhEPO, rhGM-CSF | Detection of CFU-E, BFU-E, CFU-GM, CFU-GEMM in BM, MPB, PB, and CB |
| H4534 Classic Without EPO (GF H4534) | 04534 04544 | 100 mL 24 x 3 mL | • | • | • | | rhSCF, rhIL-3, rhGM-CSF; no rhEPO | Detection of CFU-GM in BM, MPB, PB, and CB |
| H4435 Enriched (GF+ H4435) | 04435 04445 | 100 mL 24 x 3 mL | • | • | • | | rhSCF, rhIL-3, rhIL-6, rhEPO, rhG-CSF, rhGM-CSF | <ul style="list-style-type: none"> Detection of CFU-E, BFU-E, CFU-GM, CFU-GEMM in BM, MPB, PB, and CB Recommended for CD34+-enriched cells and cells isolated by other purification methods |
| H4535 Enriched Without EPO (GF+ H4535) | 04535 04545 | 100 mL 24 x 3 mL | • | • | • | | rhSCF, rhIL-3, rhIL-6, rhG-CSF, rhGM-CSF; no rhEPO | <ul style="list-style-type: none"> Detection of CFU-GM in BM, MPB, PB, and CB Recommended for CD34+-enriched cells and cells isolated by other purification methods |
| SF H4436 | 04436 | 100 mL | • | | • | • | rhSCF, rhIL-3, rhIL-6, rhEPO, rhG-CSF, rhGM-CSF | Detection of CFU-E, BFU-E, CFU-GM, CFU-GEMM in BM, MPB, PB, and CB where a medium of defined composition is required |
| SF H4536 | 04536 | 100 mL | • | | • | • | rhSCF, rhIL-3, rhIL-6, rhG-CSF, rhGM-CSF; no rhEPO | Detection of CFU-GM in BM, MPB, PB, and CB where a medium of defined composition is required |
| SF H4636 | 04636 | 100 mL | • | | • | • | rhSCF, rhIL-3, rhG-CSF, rhGM-CSF, rhEPO | <ul style="list-style-type: none"> Detection of CFU-E, BFU-E, CFU-GM, and CFU-GEMM Recommended for the culture of human hPSC-derived hematopoietic progenitor cells in defined serum-free conditions Use in CFU assays of human primary hematopoietic progenitor cells isolated from BM, MPB, PB, and CB |
| SFH4636 Without EPO | 100-0945 | 100 mL | • | | • | • | rhSCF, rhIL-3, rhG-CSF, rhGM-CSF; no rhEPO | <ul style="list-style-type: none"> Detection of CFU-GM Use in CFU assays of human primary hematopoietic progenitor cells isolated from BM, MPB, PB, and CB |
| Express | 04437 04447 | 100 mL 24 x 3 mL | • | • | • | | rhEPO and other cytokines | <ul style="list-style-type: none"> Rapid CFU assays of human CB cells. Total colonies can be counted as early as 7 days after plating, without identification of colony type. If cultures are maintained for 14 - 16 days, BFU-E, CFU-GM, CFU-G, CFU-M, and CFU-GEMM colonies can be counted. Compatible with STEMvision™ |
| H4431 | 04431 | 100 mL | • | • | • | | Agar-LCM, rhEPO | <ul style="list-style-type: none"> Detection of CFU-E, BFU-E, CFU-GM, CFU-GEMM in BM and PB Suitable as a control medium for the detection of "EPO-independent" erythroid progenitor cells using MethoCult™ H4531 |
| H4531 | 04531 | 100 mL | • | • | • | | Agar-LCM; no rhEPO | <ul style="list-style-type: none"> Detection of CFU-GM in BM and PB Suitable for detection of "EPO-independent" erythroid progenitor cells |
| H4330 | 04330 | 90 mL | • | • | • | | contains serum, rhEPO, no other cytokines | Allows researchers to add cytokines of their choice for applications including: |
| H4230 | 04230 | 80 mL | • | • | • | | contains serum, no cytokines | <ul style="list-style-type: none"> Drug toxicity testing in vitro Detection of specific hematopoietic progenitor cells Investigation of the action of novel factors |
| SF H4236 | 04236 | 80 mL | • | | • | • | serum-free, contains serum substitute, no cytokines | <ul style="list-style-type: none"> Hematopoietic colony assays in other species Detection of genetically modified hematopoietic progenitor cells |
| H4100 | 04100 | 40 mL | • | | | | base methylcellulose, no serum, serum substitutes or cytokines | <ul style="list-style-type: none"> Cloning and selection of non-adherent cell lines |

For a full listing of products for analysis of HSPCs, visit www.stemcell.com/HSPCworkflow, under the "Analysis" tab.

CFU-E: colony-forming unit – erythroid; BFU-E: burst-forming unit – erythroid; CFU-GM: colony-forming unit – granulocyte/macrophage; CFU-GEMM: colony-forming unit – granulocyte, erythrocyte, macrophage, megakaryocyte; BM: bone marrow; MPB: mobilized peripheral blood; PB: mobilized peripheral blood; CB: cord blood; hPSC: human pluripotent stem cell; EPO: erythropoietin; MC: methylcellulose; FBS: fetal bovine serum; BSA: bovine serum albumin.

MethoCult™ Media for Mouse Cells

| MethoCult™ Product | Catalog # | Unit Size | Components | | | | | Applications |
|--------------------|----------------|---------------------|------------|-----|-----|-----------------------|--|--|
| | | | MC | FBS | BSA | Insulin + Transferrin | Growth Factors | |
| GF M3434 | 03434 03444 | 100 mL 24 x 3 mL | • | • | • | • | rmSCF, rmlL-3, rhIL-6, rhEPO | <ul style="list-style-type: none"> Assays of mouse hematopoietic progenitor cells (BFU-E, CFU-GM, CFU-G, CFU-M, CFU-GEMM) in BM, PB, spleen, and fetal liver Compatible with STEMvision™ |
| GF M3534 | 03534 | 100 mL | • | • | • | • | rmSCF, rmlL-3, rhIL-6; no rhEPO | <ul style="list-style-type: none"> Assays of mouse hematopoietic progenitor cells (CFU-GM, CFU-G, CFU-M) in BM, PB, spleen, and fetal liver Compatible with STEMvision™ |
| M3630 | 03630 | 100 mL | • | • | | | rhIL-7 | Assays of mouse pre-B clonogenic progenitor cells from BM and Whitlock-Witte long-term cultures |
| SF M3236 | 03236 | 80 mL | • | | • | • | no cytokines | Assays of mouse hematopoietic progenitor cells from BM, PB, spleen, and fetal liver where a medium of defined composition is required |
| SF M3436 | 03436 | 100 mL | • | | • | • | rhEPO and other cytokines | <ul style="list-style-type: none"> Assays of erythroid progenitor cell (BFU-E) derived colonies from BM and other tissues Compatible with STEMvision™ |
| M3334 | 03334 | 90 mL | • | • | • | • | contains serum, rhEPO, no additional cytokines | Detection of CFU-E and mature BFU-E from mouse BM, spleen, and fetal liver |
| M3234 | 03234 | 80 mL | • | • | • | • | contains serum, no cytokines | Base medium for CFU assays, allowing addition of growth factors of choice |
| M3231 | 03231 | 80 mL | • | • | • | | contains serum, no insulin or transferrin, no cytokines | <ul style="list-style-type: none"> Base medium for CFU assays, allowing addition of growth factors of choice Cloning of cell lines |
| M3134 | 03134 | 40 mL | • | | | | base methylcellulose, no serum, serum substitutes or cytokines | Base medium for CFU assays, allowing researchers flexibility in addition of desired components |

MC: methylcellulose; FBS: fetal bovine serum; BSA: bovine serum albumin; BM: bone marrow; MPB: mobilized peripheral blood; CB: cord blood; LCM: leukocyte-conditioned medium; BFU-E: burst-forming unit – erythroid; CFU-E: colony-forming unit – erythroid; CFU-G: colony-forming unit – granulocyte; CFU-GM: colony-forming unit – granulocyte/macrophage; CFU-M: colony-forming unit – macrophage; CFU-GEMM: colony-forming unit – granulocyte, erythrocyte, macrophage, megakaryocyte.

STEMvision™

Automated and Standardized Counting of CFU Assays

STEMvision™ is an instrument and software system designed for imaging and counting hematopoietic colonies in the colony-forming unit (CFU) assay.

Instead of manually counting colonies, users simply load a SmartDish™ (Catalog #27370) plate containing cells plated in MethoCult™ medium into STEMvision™. Digital images of individual 35 mm wells of the 6-well SmartDish™ plate containing human or mouse cells are acquired in approximately 1 minute or 30 seconds, respectively. Analysis requires approximately 1 minute for assays of human cells and 30 seconds for assays of mouse cells, and can be performed at a later time or overnight.

For human cultures, three STEMvision™ software packages have been designed for scoring and counting colonies produced by erythroid, myeloid, and multi-potential progenitor cells (BFU-E, CFU-E, CFU-GM/G/M, and CFU-GEMM) in 14-day CFU assays of cells isolated from cord blood (CB), bone marrow, and mobilized peripheral blood, cultured in MethoCult™ Optimum. A fourth software package is available to count the total number of colonies (without classification of CFU subtypes) in 7-day CFU assays of CB cells in MethoCult™ Express.

For mouse bone marrow cultures, three STEMvision™ software packages have been designed to count the total number of colonies produced by all combined myeloid and erythroid progenitor cells (BFU-E and CFU-GM/G/M) in MethoCult™ GF M3434, all myeloid progenitor cells (CFU-GM/G/M) in MethoCult™ GF M3534, and all erythroid progenitor cells (BFU-E) in MethoCult™ SF M3436.

SYSTEM IS SUPPLIED WITH:

- STEMvision™ instrument
- Computer and monitor
- Software for image acquisition, analysis, and review

REQUIRED REAGENTS:

- SmartDish™ meniscus-free cultureware
- MethoCult™ medium
- Red blood cell depletion reagent

CAPACITY:

- One 6-well SmartDish™ per run



STEMvision™ Instrument

Visit www.STEMvision.com for more information.

Why Use STEMvision™?

ACCURACY. Achieve faster and more accurate colony counts than with manual counting.

STANDARDIZE. Eliminate the subjectivity associated with manual colony counting.

CONVENIENT. Easily save and share your results. Images are automatically stored for record keeping.

Software for Human Assays

PRODUCT: STEMvision™ Human Cord Blood 7-Day CFU Analysis Package
CATALOG #: 22001

PRODUCT: STEMvision™ Human Cord Blood 14-Day CFU Analysis Package
CATALOG #: 22005

PRODUCT: STEMvision™ Human Bone Marrow 14-Day CFU Analysis Package
CATALOG #: 22006

PRODUCT: STEMvision™ Human Mobilized Peripheral Blood 14-Day CFU Analysis Package
CATALOG #: 22007

Software for Mouse Assays

PRODUCT: STEMvision™ Mouse Total CFU Analysis Package
CATALOG #: 22008

PRODUCT: STEMvision™ Mouse Myeloid CFU Analysis Package
CATALOG #: 22009

PRODUCT: STEMvision™ Mouse Erythroid CFU Analysis Package
CATALOG #: 22011

Software for Quality Management

PRODUCT: STEMvision™ 21 CFR Part 11 Compliance Software Add-On
CATALOG #: 500-0110

Human Cord Blood CFU Assays




Figure 27. Representative STEMvision™ Images Showing Colonies Derived from Human Cord Blood Progenitor Cells After 7 Days of Culture in MethoCult™ Express or After 14 Days of Culture in MethoCult™ Optimum

The images have been analyzed with the STEMvision™ Human Cord Blood (A) 7-Day (Catalog #22001) and (B) 14-Day (Catalog #22005) Analysis Packages. Green circles identify individual colonies in the 7-day CB CFU assay that counts the total number of CFUs. Red circles identify erythroid colonies (produced by BFU-E), yellow circles identify myeloid colonies (produced by CFU-GM, CFU-G or CFU-M), and blue circles identify mixed colonies (produced by CFU-GEMM) in the 14-day CB CFU assay. Erythroid and mixed colonies that contain hemoglobinized cells are shown in true red color.




Figure 28. STEMvision™ Automated Counts of Total, Erythroid (BFU-E), and Myeloid (CFU-G/M/GM) Colonies Are Highly Correlated to Manual Counts of 14-Day CB CFU Assays

Cryopreserved CB samples were thawed, plated in MethoCult™ Optimum, cultured for 14 days, and scored both manually using an inverted microscope and automatically using STEMvision™. The results show a strong correlation between automated counts using STEMvision™ and manual counts. Gray dashed lines represent a perfect linear correlation between manual and automated counts. Orange solid lines represent the actual linear correlation between manual and automated counts (n = 130 CFU assays).




Figure 29. STEMvision™ Automated Colony Counting of 7-Day and 14-Day CB CFU Assays Is More Reproducible Than Manual Counting

The coefficients of variation (CV) for total colony counts in (A) 7-day and (B) 14-day CFU assays of CB cells were determined by counting the same culture wells manually by 3 to 5 different people (gray diamonds) and automatically using 3 to 5 separate STEMvision™ instruments (orange squares). The average CVs for 7-day and 14-day total colony counts produced manually were 11% and 13%, respectively. CVs for 7-day and 14-day colony counts produced by STEMvision™ were 5%.

Mouse Bone Marrow CFU Assays

MethoCult™ GF M3434




Figure 30. Representative STEMvision™ Images Showing Colonies Derived from Mouse Bone Marrow Progenitor Cells After 12 Days of Culture in MethoCult™ GF M3434 Medium

Images of mouse bone marrow cells cultured in MethoCult™ GF M3434 medium were acquired using STEMvision™. The Total CFU Analysis Package (Catalog #22008) was used to analyze the image. Red circles identify the smallest colonies; size class 1, yellow circles; size class 2, blue circles; size class 3, and orange circles identify the largest colonies; size class 4.




Figure 31. STEMvision™ Automated Total Colony Counts (Myeloid Plus Erythroid) on Days 7, 10, and 12 Are Highly Correlated to Manual Counts of CFU Assays Using Mouse Bone Marrow Cells

Bone marrow cells were plated in MethoCult™ GF M3434. Colonies were counted on days (A) 7, (B) 10, and (C) 12 both manually using an inverted microscope, and automatically using STEMvision™ equipped with the Mouse Total CFU Analysis Package (Catalog #22008). We recommend counting CFU assays of mouse progenitor cells plated in M3434 between 7 and 12 days. Gray dashed lines represent a theoretical perfect linear correlation between manual and automated counts. Orange solid lines represent the actual linear correlation between manual and automated counts. Coefficients of determination (R^2) that describe each data set ($n = 104$ CFU assays in A, $n = 38$ in B, and $n = 99$ in C) are shown in orange.




Figure 32. STEMvision™ Automated Colony Counting of CFU Assays with Mouse Cells Is More Reproducible Than Manual Counting

The coefficients of variation (CVs) for total colony counts in CFU assays of mouse BM in MethoCult™ GF M3434 were determined by counting the same culture dishes manually by 3 to 5 different people (grey diamonds) and automatically using 3 to 5 separate STEMvision™ instruments (orange squares). The average CV for total colony counts produced manually was 16%. The average CV for colony counts produced by STEMvision™ was 4%.

SmartDish™ and STEMgrid™-6

Meniscus-Free Cultureware for More Accurate Counting of Hematopoietic Colonies

When the CFU assay is performed using traditional cultureware, the medium forms a meniscus where it meets the dish. The greater medium depth in the meniscus results in a greater proportion of colonies forming around the periphery (Figure 33A) where optical distortion can make it challenging to identify colonies (Figures 34A and 34C). This can reduce the accuracy of colony counting (i.e. result in undercounting of CFUs).

SmartDish™ 6-well culture plates have been designed to enable accurate and reproducible colony counting by preventing the formation of a meniscus. This allows for more uniform distribution of culture medium (Figure 33B), which results in a more even distribution of colonies throughout the entire 35 mm well. Additionally, the absence of the meniscus reduces optical distortion so that colonies present at the edge of the dish can be more easily counted (Figures 34B and 34D).

SmartDish™ has been designed to work with standard inverted microscopes when performing manual colony counts (Figure 34B), as well as with the STEMvision™ instrument for automated counting of hematopoietic CFU assays (Figure 34D). For manual counting of hematopoietic CFU assays, SmartDish™ is used with STEMgrid™-6, a detachable counting grid that facilitates navigation throughout the culture and also divides it into four quadrants for partial counting if desired.

SmartDish™ Meniscus-Free Cultureware

| | |
|-------------------|----------------------------|
| PRODUCT: | SmartDish™ (6-well plates) |
| CATALOG #: | 27370 (5/pack) |
| | 27371 (50/pack) |

RECOMMENDED FOR:

- Easier and more accurate colony counting
- Automated colony counting with STEMvision™ (required)

STEMgrid™-6 Counting Grid

| | |
|-------------------|----------------|
| PRODUCT: | STEMgrid™-6 |
| CATALOG #: | 27000 (1/pack) |

RECOMMENDED FOR:

- Manual colony counting in SmartDish™ cultureware under an inverted microscope





Figure 33. Schematic Illustration of Medium and Colony Distribution in Standard and SmartDish™ Cultureware

(A) The formation of a meniscus in standard cultureware results in more colonies forming around the dish edge where the culture medium is deeper. (B) The absence of a meniscus in SmartDish™ cultureware ensures more uniform distribution of culture medium and colonies throughout the entire well.

Manual Counting of Hematopoietic CFU Assays



Automated Counting of Hematopoietic CFU Assays




Figure 34. 14-Day Human Cord Blood CFU Assays Performed in Standard Non-Treated and SmartDish™ 6-Well Culture Plates

Shown are representative STEMvision™ images of 35 mm wells from either a (A, C) non-treated 6-well plate or (B, D) SmartDish™. The formation of a meniscus in standard cultureware causes more colonies to form around the periphery of the dish where the culture medium is deeper (A, C). Optical distortion obscures these colonies and makes them more difficult to count. Colonies are easier to count at the edge of the SmartDish™, which has been treated to eliminate the meniscus, allowing a more equal distribution of colonies (B, D).

MegaCult™ Collagen-Based Media

For the Detection of Megakaryocyte Progenitor Cells and for Permanent Records of Hematopoietic Progenitor Assays

Why Use Collagen-Based Media?

Cultures in collagen-based media can be dehydrated and fixed and have been shown to support the proliferation of hematopoietic progenitor cells.⁷⁷ This is beneficial when quantifying and detecting megakaryocyte colonies, which cannot be distinguished morphologically from macrophage colonies. Their identification requires the use of staining procedures to identify the expression of megakaryocyte-specific cell surface markers or enzymatic activity. Cellular and molecular analysis of cells cultured in methylcellulose-based media requires that the colonies be plucked and processed, which is time-consuming. In contrast, dehydrating and fixing a collagen-based culture prior to staining only requires about 30 minutes. Collagen-based cultures that have been dehydrated, fixed, and stained can be maintained as long-term records. This is in contrast to methylcellulose-based cultures that can only be maintained for approximately one week after the culture period.

Why Use MegaCult™?

CONVENIENT. Culture and stain your cells on one slide and evaluate immediately or store for examination at a later time.

OPTIMIZED. Quantify human megakaryocytic progenitor cells in serum-free conditions.

FLEXIBLE. Customize your cell-culture conditions by adding desired cytokines for the assay of human or mouse megakaryocytic progenitor cells.

For more information, visit www.MegaCult.com.

Human and Mouse MegaCult™ Media and Staining Kits

| Product Name | Catalog # | Components | Applications |
|---|-----------|---|--|
| MegaCult™-C Complete Kit with Cytokines | 04971 | <ul style="list-style-type: none"> Serum-free medium with recombinant cytokines (24 x 2 mL) Collagen solution (35 mL) CFU-Mk staining kit Double chamber slides (48) | Detection and staining of human megakaryocyte progenitor cells in light density or CD34+-enriched BM, MPB, and CB cells |
| MegaCult™-C Complete Kit Without Cytokines | 04970 | <ul style="list-style-type: none"> Serum-free medium without cytokines (24 x 1.7 mL) Collagen solution (35 mL) CFU-Mk staining kit Double chamber slides (48) | Detection and staining of human megakaryocyte progenitor cells in light density or CD34+-enriched BM, MPB, and CB cells (requires addition of appropriate cytokines) |
| MegaCult™-C Collagen and Medium with Cytokines | 04961 | <ul style="list-style-type: none"> Serum-free medium with recombinant cytokines (24 x 2 mL) Collagen solution (35 mL) | Detection of human megakaryocyte progenitor cells in light density or CD34+-enriched BM, MPB, and CB cells |
| MegaCult™-C Collagen and Medium Without Cytokines | 04960 | <ul style="list-style-type: none"> Serum-free medium without cytokines (24 x 1.7 mL) Collagen solution (35 mL) | Detection of human or mouse megakaryocyte or other progenitor cells (requires addition of appropriate cytokines) |
| MegaCult™-C Collagen and Medium with Lipids | 04974 | <ul style="list-style-type: none"> Serum-free medium with lipids, without cytokines (50 mL) Collagen solution (35 mL) | Detection of human or mouse megakaryocyte or other progenitor cells (requires addition of appropriate cytokines) |
| MegaCult™-C Staining Kit for CFU-Mk | 04962 | <ul style="list-style-type: none"> Anti-CD41 primary antibody Anti-TNP control antibody Biotin-conjugated secondary antibody Alkaline phosphatase detection system Human serum for dilutions and BSA for blocking Evans Blue counterstain | Immunocytochemical staining for detection of human megakaryocytes and platelets in CFU-Mk and BFU-E/Mk grown in MegaCult™-C |

CFU-Mk: colony-forming unit – megakaryocyte; BFU-E/Mk: burst-forming unit – erythroid/megakaryocyte; BM: bone marrow; MPB: mobilized peripheral blood; CB: cord blood.

Antibodies and ELISA Kits for HSPC Research

Analyze cells with antibodies that are verified to work with STEMCELL Technologies' cell isolation and cell culture reagents for select applications. These high-quality antibodies ensure consistent results for your downstream cell analysis, such as phenotyping and purity assessments of primitive cells and expanded progenitors.

For a complete product listing, including secondary antibodies and isotype controls, visit www.stemcell.com/antibodies. Alternatively, to view available ELISA kits for HSPC research, including the Erythropoietin (EPO) ELISA Kit (Catalog #01630), visit www.stemcell.com/ELISA.

Antibodies for Human Hematopoietic Cell Research

| Antigen | Clone | Isotype | Conjugate | Catalog # |
|-----------------------------|--------|---------------------------------|---|-----------|
| CD3 | SK7 | Mouse IgG ₁ , kappa | APC, FITC, PE, PerCP | 60127 |
| | UCHT1 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP, PerCP-Cy5.5 | 60011 |
| CD11b | ICRF44 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60040 |
| | M1/70 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, Pacific Blue™, PE, PerCP-Cy5.5 | 60001 |
| CD14 | MoP9 | Mouse IgG _{2b} , kappa | FITC | 60124 |
| | M5E2 | Mouse IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60004 |
| CD16 | 3G8 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60041 |
| CD19 | HIB19 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60005 |
| | 6D5 | Rat IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, Pacific Blue™, PE | 60006 |
| CD20 | 2H7 | Mouse IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60008 |
| CD34 | 8G12 | Mouse IgG ₁ , kappa | APC, FITC, PE | 60121 |
| | 563 | Mouse IgG ₁ , kappa | PE | 60119 |
| | 581 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60013 |
| CD38 | AT-1 | Mouse IgG ₁ , kappa | FITC | 60131 |
| | HIT2 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60014 |
| CD45 | 2D1 | Mouse IgG ₁ , kappa | FITC | 60123 |
| | HI30 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, Pacific Blue™, PE, PerCP-Cy5.5 | 60018 |
| CD45RA | HI100 | Mouse IgG _{2b} , kappa | Unconjugated, APC-Cyanine7, Biotin, PE | 100-0318 |
| CD71 (Transferrin Receptor) | OKT9 | Mouse IgG ₁ , kappa | Unconjugated, APC, Biotin, FITC, PE | 60106 |
| CD90 (Thy-1) | 5E10 | Mouse IgG ₁ , kappa | Unconjugated, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60045 |
| CD105 | 43A3 | Mouse IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60039 |
| CD117 (cKit) | 104D2 | Mouse IgG ₁ , kappa | Unconjugated, APC, Biotin, PE | 60087 |
| CD123 (IL-3Rα) | 6H6 | Mouse IgG ₁ , kappa | Unconjugated, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60110 |
| CD235a (Glycophorin A) | 2B7 | Mouse IgG ₁ , kappa | FITC | 60152 |
| CD235ab (Glycophorin A/B) | HIR2 | Mouse IgG _{2b} , kappa | Unconjugated, APC, Biotin, FITC, PE | 60111 |

Antibodies for Mouse Hematopoietic Cell Research

| Antigen | Clone | Isotype | Conjugate | Catalog # |
|--------------------|--------------|--|---|-----------|
| CD3 | 17A2 | Rat IgG _{2b} , kappa | Unconjugated, APC, FITC, PerCP-Cy5.5 | 100-1606 |
| CD3e | 145-2C11 | Hamster (Armenian) IgG ₁ , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60015 |
| | GK1.5 | Rat IgG _{2b} , kappa | Unconjugated, APC, FITC, PE, PerCP-Cy5.5 | 100-1601 |
| CD4 | RM4-5 | Rat IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60017 |
| | RM4-4 | Rat IgG _{2b} , kappa | APC, Biotin, FITC, PE | 60029 |
| CD8a | 2.43 | Rat IgG _{2b} , kappa | Unconjugated, APC, FITC | 100-1633 |
| | 53-6.7 | Rat IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60023 |
| CD11b | M1/70 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, Pacific Blue™, PE, PerCP-Cy5.5 | 60001 |
| CD19 | 6D5 | Rat IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, Pacific Blue™, PE | 60006 |
| CD43 | R2/60 | Rat IgM, kappa | Biotin | 60042 |
| CD45.1 | A20 | Mouse IgG _{2a} , kappa | APC, Biotin, FITC, PE | 60117 |
| CD45.2 | 104 | Mouse IgG _{2a} , kappa | APC, Biotin, FITC, PE | 60118 |
| CD45 | 30-F11 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60030 |
| CD48 (SLAMF2) | HM48-1 | Hamster (Armenian) IgG ₁ , lambda | APC | 60162 |
| CD45R (B220) | RA3-6B2 | Rat IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, Pacific Blue™, PE, PerCP-Cy5.5 | 60019 |
| CD80 (B7-1) | 16-10A1 | Hamster (Armenian) IgG | APC, FITC, PE | 100-1622 |
| CD86 (B7-2) | GL-1 | Rat IgG _{2a} , kappa | APC, FITC, PE | 100-1625 |
| CD117 (cKit) | 2B8 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60025 |
| | ACK2 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, FITC, PE | 60034 |
| CD150 (SLAM) | TC15-12F12.2 | Rat IgG _{2a} , lambda | Unconjugated, Alexa Fluor® 488, APC, Biotin, PE, PE-Cyanine7 | 60036 |
| EPCR (CD201) | RMEPCR1560 | Rat IgG _{2b} , kappa | Unconjugated, PE | 60038 |
| Gr-1 (Ly-6G/Ly-6C) | RB6-8C5 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60028 |
| Ly-6C | HK1.4 | Rat IgG _{2c} , kappa | APC, PE, PerCP-Cy5.5 | 100-1636 |
| Ly-6G | 1A8 | Rat IgG _{2a} , kappa | Unconjugated, APC, Biotin, FITC, Pacific Blue™, PE, PerCP-Cy5.5 | 60031 |
| Sca1 (Ly-6A/E) | E13-161.7 | Rat IgG _{2a} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE | 60032 |
| | T3G6 | Mouse IgG ₁ , kappa | Unconjugated | 60132 |
| Rhesus | T4G6 | Mouse IgG ₁ , kappa | | 60133 |
| TCR Beta | H57-597 | Hamster (Armenian) IgG | Unconjugated, APC, FITC, PE, PerCP-Cy5.5 | 100-1610 |
| TER119 | TER-119 | Rat IgG _{2b} , kappa | Unconjugated, Alexa Fluor® 488, APC, Biotin, FITC, PE, PerCP-Cy5.5 | 60033 |

Detection Assays for Cord Blood Potency Testing

ALDH^{br} Assay Kit



The ALDH^{br} Assay Kit (Catalog #01711) is optimized for the identification and quantitation of viable CD34⁺ cells that express high levels of the enzyme aldehyde dehydrogenase (ALDH) in human cord blood samples.

Product Information

| Products | Catalog # | Components |
|-------------------------------|-----------|---|
| ALDH ^{br} Assay Kit* | 01711 | ALDEFLUOR™ Kit <ul style="list-style-type: none">• ALDEFLUOR™ Reagent• ALDEFLUOR™ DEAB Reagent• ALDEFLUOR™ Assay Buffer• 2N HCl• DMSO Viability Dye and Antibodies <ul style="list-style-type: none">• 7-AAD Viability Dye• APC CD34 Antibody (Clone 581)• PE CD45 Antibody (Clone HI30)• PE-Cyanine5 CD235ab (GlyAB) Antibody (Clone HIR2)™ |

*Suitable for labeling up to 6 x 10⁷ cord blood cells

For more information, including data and a full protocol, visit www.stemcell.com/ALDHbrKit.

Tissue Culture Reagents and Supplies

A variety of support products are available to accompany STEMCELL Technologies' complete array of cell separation and specialized media products. For more details and a full list of tissue culture reagents and supplies, visit www.stemcell.com.

Tissue Culture Media

| Product Name | Catalog # | Unit Size |
|-----------------------------------|-----------|-----------|
| Agar Leukocyte Conditioned Medium | 02300 | 25 mL |
| DMEM with 4500 mg/L D-glucose | 36250 | 500 mL |
| DMEM with 1000 mg/L D-glucose | 36253 | 500 mL |
| DMEM/F-12 with 15mM HEPES | 36254 | 500 mL |
| Iscove's MDM (IMDM) | 36150 | 500 mL |
| IMDM with 2% FBS | 07700 | 100 mL |
| McCoy's 5A Medium | 36350 | 500 mL |
| MEM Alpha with Nucleosides | 36450 | 500 mL |
| MEM Alpha without Nucleosides | 36453 | 500 mL |
| MEM with Earle's & NEAA | 36550 | 500 mL |
| ReproTeSR™ (2-Component) | 05926 | 500mL |

Balanced Salt Solution

| Product Name | Catalog # | Unit Size |
|--------------------------|-----------|-----------|
| D-PBS | 37350 | 500 mL |
| D-PBS, 10X | 37354 | 500 mL |
| D-PBS with 2% FBS | 07905 | 500 mL |
| HBSS, Ca++ & Mg++ free | 37250 | 500 mL |
| HBSS, Without Phenol Red | 37150 | 500 mL |

Supplements

| Product Name | Catalog # | Unit Size |
|--------------------------------------|-----------|-----------|
| Heparin Solution | 07980 | 2 mL |
| Hydrocortisone Powder | 07904 | 100 mg |
| L-glutamine, 200 mM | 07100 | 100 mL |
| MEM Non-Essential Amino Acids, 100X | 07600 | 100 mL |
| Human Low-Density Lipoproteins (LDL) | 02698 | 5 mg |

Enzymes

| Product Name | Catalog # | Unit Size |
|-----------------------------------|-----------|-----------|
| Collagenase Type I (0.25%) | 07902 | 5 mL |
| DNase I, 1 mg/mL | 07900 | 1 mL |
| Trypsin in Citrate Saline (0.25%) | 07400 | 100 mL |
| Trypsin-EDTA (0.25%) | 07901 | 500 mL |
| Trypsin-EDTA (0.05%) | 07910 | 500 mL |

Serum Substitutes

| Product Name | Catalog # | Unit Size |
|--|-----------|-----------|
| 10% Bovine Serum Albumin (BSA) in Iscove's MDM | 09300 | 100 mL |
| BIT 9500 Serum Substitute | 09500 | 100 mL |

Miscellaneous Tissue Culture Reagents

| Product Name | Catalog # | Unit Size |
|------------------------------------|----------------|------------------|
| 3% Acetic Acid with Methylene Blue | 07060 | 100 mL |
| 7-AAD (7-Aminoactinomycin D) | 75001 | 500 tests |
| Ammonium Chloride Solution | 07800 07850 | 100 mL 500 mL |
| DAPI (Hydrochloride) | 75004 | 10 mg |
| Hydrocortisone | 74142 74144 | 100 mg 1 g |
| Propidium Iodide | 75002 | 10 mg |
| Trypan Blue | 07050 | 100 mL |

Tissue Culture Dishes & Slides

| Product Name | Catalog # | Unit Size |
|---|-------------------|--|
| Outer Dishes for CFU Assays | | |
| Corning® 245 mm Square Dish, Non-Treated | 38020 | 4/pack 16/case |
| Non-Adherent Culture Dishes *Recommended for CFU Assays using MethoCult™ | | |
| 35 mm Diameter Dishes | 27100 27150 | 10 Dishes 500 Dishes |
| SmartDish™ | 27370 27371 | 5/pack 50/pack |
| Other Dishes, Flasks, and Slides | | |
| 245 mm Square Dish, Tissue Culture-Treated | 38039 100-0084 | 16 Plates 4 Plates |
| Corning® 60 mm Gridded Scoring Dish | 38068 100-0085 | 500 Dishes 20 Dishes |
| Corning® Disposable Erlenmeyer Flask, 500 mL | 38013 | 25 Flasks |
| Corning® Disposable Erlenmeyer Flask, 250 mL | 38012 | 50 Flasks |
| Corning® Disposable Erlenmeyer Flask, 125 mL | 38011 | 50 Flasks |
| Double Chamber Slide Kit | 04963 | 48 double chamber slides, filter cards and spacers |
| Falcon® 24-Well Flat-Bottom Plate, Tissue Culture-Treated | 38021 | 50 Plates |
| Falcon® 96-Well Flat-Bottom Microplate, Tissue Culture-Treated | 38022 | 50 Plates |
| Falcon® Round-Bottom Tubes, 14 mL | 38008 | 500 Tubes |
| Tissue Culture-Treated Dishes, 100 mm | 38046 100-0082 | 500 Dishes 20 Dishes |
| Tissue Culture-Treated Dishes, 100 mm | 27125 27127 | 10 Dishes 240 Dishes |

Miscellaneous Tissue Culture Supplies

| Product Name | Catalog # | Unit Size |
|----------------------------|----------------|--------------------------------|
| 3 cc Syringes | 28240 | 100 Syringes/bag |
| Blunt-End Needles | 28110 28120 | 100/pack 2000/case |
| Reversible Strainer, 37 µm | 27215 27250 | Small, 50/box Large, 25/box |
| Hypoxia Chamber | 27310 | 1 chamber |

Proficiency Testing Programs and Quality Control Kits

STEMCELL Technologies is committed to helping you get the most out of your experiments. Our standardization tools include proficiency testing programs, quality control kits, training courses, and instructional materials to improve your competency in performing hematopoietic colony-forming unit (CFU) assays.

Proficiency Testing Programs

Evaluate your ability to perform all steps of the CFU assay, from thawing samples to plating cells in methylcellulose medium, in our Proficiency Testing programs. Participants are provided with a cell sample, MethoCult™ medium, additional reagents and supplies, and detailed instructions required to perform the CFU assay. The results from all participants are analyzed and compared to the cohort mean (Figure 35) according to guidelines outlined in ISO 13528.




Figure 35. Example Participant Data

Quality Control Kits

STEMCELL Quality Control Kits are recommended for technologists in cell processing laboratories wanting to monitor their ability to reproducibly set up and score colonies in hematopoietic CFU assays. The use of a standardized medium and identical aliquots of cells from a single source minimizes potential variability associated with these components of the assay, allowing a more accurate assessment of an individual's technique. Each kit includes sufficient supplies to perform monthly CFU assays over a 1 year period.



Proficiency Testing Supplies

| Proficiency Testing Program | | Catalog # |
|-------------------------------|----------------|----------------------|
| Frozen Human Bone Marrow (BM) | Spring Fall | 100-0926 100-0928 |
| Fresh Cord Blood (CBH) | Spring | 00606 |
| Frozen Cord Blood (CBZ) | Spring Fall | 100-0950 100-0952 |

For information on upcoming program dates and how to register, visit www.ProficiencyTesting.com.



Quality Control Kit Supplies

| Quality Control Kit | Catalog # | Application |
|------------------------|-----------|--|
| Human Bone Marrow (BM) | 00650 | Monitoring reproducibility in performing BM CFU assays |
| Human Cord Blood (CB) | 00651 | Monitoring reproducibility in performing CB CFU assays |


Training Courses and Instructional Materials

Our hematopoietic assay training courses are comprehensive and provide both theoretical and hands-on training. Our instructors share their knowledge and expertise to help overcome challenges in assay design, set up, and evaluation.


Live Training Course for Hematopoietic Progenitor Assays

| Course Name | Catalog # | Course Description |
|---|-----------|---|
| Standardization of the Hematopoietic Progenitor Assay | 00215 | <p>This 2-day course focuses on the standardization of the colony-forming unit (CFU) assay for human samples</p> <ul style="list-style-type: none"> Gain hands-on experience in the assay set-up, identification, and enumeration of hematopoietic progenitor cells using MethoCult™ |


STEMCELL Technologies also offers customized training courses, which can be highly beneficial to address specific research challenges or to standardize technical processes. Custom training packages provide your team with personalized instruction from a technical specialist at your facility or at our training laboratory. For more information about courses or customized training, visit www.stemcell.com/training or contact education@stemcell.com.



Catalog #29940



Catalog #28700



Catalog #28760

Instructional Materials

| Product Name | Catalog # | Description |
|--|-----------|---|
| Cord Blood Colony Atlas | 29940 | Detailed color images of hematopoietic colonies derived from human umbilical cord blood progenitor cells grown in methylcellulose-based media. Request a free hard copy at www.stemcell.com/cbc-atlas . |
| Human Hematopoietic Colony Atlas | 28700 | A guide for the process of identifying and evaluating colonies derived from human hematopoietic progenitor cells grown in methylcellulose-based media. Request a free hard copy at www.stemcell.com/human-atlas . |
| Human Hematopoietic Colonies in Health and Disease | 28760 | Colony-morphology from healthy individuals and patients with hematological disorders are illustrated in 230 detailed, color photographs. |

Technical Manuals for Hematopoietic Progenitor Assays

Technical Manuals are available to support hematopoietic media products and assay systems. These manuals provide detailed reagent information, step-by-step instructions for use, and valuable tips and hints. Print copies are available free of charge upon request and can also be downloaded at www.stemcell.com/technical-resources/product-information.html.

Free, On-Demand Training



Mouse CFU Assay

Learn to Set Up and Perform the CFU Assay for Mouse HSPCs
www.stemcell.com/hsc-mouse-training



Human CFU Assay


Get Step-by-Step Instruction on Set Up, Colony Identification, and Scoring
www.stemcell.com/hsc-training



CFU Assay Automation

Learn to Use STEMvision™ to Standardize CFU Scoring
www.stemcell.com/stemvision-training

Contract Assay Services



Contract Assay Services (CAS) is a contract research organization (CRO) established within STEMCELL Technologies that performs assay services based on in vitro and in vivo primary stem cell assays. Primary cells are thought to be more representative of in vivo functionality than cell lines, and can increase the biological relevance of in vitro assays. CAS specializes in providing the CFU assay for a variety of applications (see below), in addition to their portfolio of characterized assays and custom solutions for your individual needs. Since 2000, CAS has performed such studies for over 120 pharmaceutical, biotechnology, government, and academic life science organizations worldwide, and provides exceptional service through frequent communication, quality products, and unparalleled expertise.


Stem and Progenitor Cell Assays

Our in vitro assays can help you determine a compound's inhibitory or stimulatory effects on hematopoietic progenitor cells:


- CFU assays have been shown to be important in the evaluation of the potential inhibitory or stimulatory effects of a variety of compounds on hematopoietic and mesenchymal progenitor cells^{72,78-86}
- CFU assays can be used to assess proliferation and differentiation of hematopoietic progenitor cells to determine IC_{50} and IC_{90} values
- CFU assays for myeloid progenitor cells have been validated for the determination of maximum tolerated dose by the European Centre for the Validation of Alternative Methods (ECVAM)
- Assays for CFU-GM and CFU-Mk have been shown to be predictive of clinical outcomes such as neutropenia and thrombocytopenia^{72,81,82}

Our in vivo assays help you examine the effects of a compound on hematopoietic stem cells:


- Assess hematopoietic stem cell mobilization into the peripheral blood
- Determine kinetics of hematopoietic progenitor cell recovery following induction of myelosuppression




CFU-GM



CFU-GM with Inhibitory Compound



BFU-E



BFU-E with Inhibitory Compound

Figure 36. Effect of the Addition of an Inhibitory Compound on the Morphology of Human Bone Marrow CFU-GM and BFU-E Colonies



Interested in CAS?

Learn more about Contract Assay Services at STEMCELL Technologies
www.contractassay.com




Book a Consultation

Contact us to find how we can help you meet your goals
www.stemcell.com/contact-cas

Services for Cell Therapy Program

If you are looking to use STEMCELL Technologies products as ancillary materials (AMs) or raw materials in the manufacture of hematopoietic cellular therapies, we can work with you to develop solutions for your specific needs, such as custom product development, higher compliance product manufacture, regulatory support, and customer documentation. Under the Services for Cell Therapy Program, STEMCELL may be able to support the use of our products as AMs (as defined under USP Chapter <1043> or Ph. Eur. General Chapter 5.2.12) under an approved Investigational New Drug (IND) application, Clinical Trial Authorization/Application (CTA), or equivalent regulatory filing. We have successfully supported over 100 clinical trials around the globe.

Overview of STEMCELL's Custom Services for Cell Therapies Process



Why Choose the Services for Cell Therapy Program?

GLOBAL EXPERIENCE. We have worked with customers around the world to meet their country- or region-specific regulatory requirements.

RELIABLE. Our stringent supplier qualification and management processes ensure consistent quality and a dependable supply of materials and services.

CUTTING-EDGE. We can support your workflows from start to finish with innovative, specialized reagents.

HIGH-QUALITY. In-house manufacturing and testing of select products under relevant cGMPs with the ability to support AM qualification under USP <1043> or Ph. Eur. 5.2.12.

ACCESSIBLE. Our network of qualified distribution centers enables worldwide delivery of our products.

Taking Your Research to the Clinic?

STEMCELL's Services for Cell Therapy program has a team of experts who can help support your regulatory filing by providing custom solutions such as quality documentation, additional product testing, and customized product manufacturing. To learn more about how we can support your preclinical and clinical research needs, visit us at www.stemcell.com/services/cell-therapy.html.

References

1. Delaney C et al. (2010) *Nat Med* 16(2): 232–6.
2. Cutler C et al. (2013) *Blood* 122(17): 3074–81.
3. Leberbauer C et al. (2005) *Blood* 105(1): 85–94.
4. Cantù C et al. (2011) *Blood* 117(13): 3669–79.
5. Flygare J et al. (2011) *Blood* 117(12): 3435–44.
6. Heckl D et al. (2011) *Blood* 117(14): 3737–47.
7. Satchwell TJ et al. (2011) *Blood* 118(1): 182–91.
8. Gaikwad A et al. (2007) *Exp Hematol* 35(4): 587–95.
9. Kumkhaek C et al. (2013) *Blood* 121(16): 3216–27.
10. Ohmine S et al. (2011) *Stem Cell Res Ther* 2(6): 46.
11. Lechman ER et al. (2012) *Cell Stem Cell* 11(6): 799–811.
12. Chin JY et al. (2013) *Mol Ther* 21(3): 580–7.
13. Boitano AE et al. (2010) *Science* 329(5997): 1345–8.
14. Fares I et al. (2014) *Science* 345(6203): 1509–12.
15. Pabst S et al. (2014) *Nature Meth* 11(4): 436–42.
16. Dexter TM et al. (1977) *J Cell Physiol* 91(3): 335–44.
17. Gartner S & Kaplan HS (1980) *Proc Natl Acad Sci USA* 77(8): 4756–9.
18. Sutherland HJ et al. (1990) *Proc Natl Acad Sci USA* 87(9): 358–8.
19. Lemieux ME et al. (1995) *Blood* 86(4): 1339–74.
20. Conneally E et al. (1997) *Proc Natl Acad Sci USA* 94(18): 9836–41.
21. Hogge DE et al. (1996) *Blood* 88(10): 3765–73.
22. Sutherland HJ et al. (1991) *Blood* 78(3): 666–72.
23. Ploemacher RE et al. (1991) *Blood* 78(10): 2527–33.
24. Sutherland HJ et al. (1989) *Blood* 74(5): 1563–70.
25. Prosper F et al. (1997) *Blood* 89(11): 3991–7.
26. Ponchio L et al. (1995) *Blood* 86(9): 3314–21.
27. Petzer AL et al. (1996) *Proc Natl Acad Sci USA* 93(4): 1470–4.
28. Petzer AL et al. (1996) *J Exp Med* 183(6): 2551–8.
29. Zandstra PW et al. (1997) *Proc Natl Acad Sci USA* 94(9): 4698–703.
30. Fraser CC et al. (1990) *Blood* 76(6): 1071–6.
31. Stewart AK et al. (1997) *Cancer Gene Ther* 4(3): 148–56.
32. Zandstra PW et al. (1994) *Biotechnology* 12(9): 909–14.
33. Kogler G et al. (1998) *Bone Marrow Transplant Suppl* 3: S48–53.
34. Cashman JD et al. (1990) *Blood* 75(1): 96–101.
35. Eaves CJ et al. (1991) *Blood* 78(1): 110–7.
36. Verfaillie CM (1993) *Blood* 82(7): 2045–53.
37. Ghaffari S et al. (1997) *Br J Haematol* 97(1): 22–8.
38. Udomsakdi C et al. (1992) *Proc Natl Acad Sci USA* 89(13): 6192–6.
39. Eaves CJ et al. (1993) *Proc Natl Acad Sci USA* 90(24): 12015–9.
40. Petzer AL et al. (1996) *Blood* 88(6): 2162–1.
41. Ailles LE et al. (1997) *Blood* 90(7): 2555–64.
42. Maciejewski JP et al. (1996) *Blood* 88(6): 1983–91.
43. Cavazzana-Calvo M et al. (1996) *Blood* 88(10): 3901–9.
44. Gong JH et al. (1994) *Leukemia* 8(4): 652–8.
45. Eaves CJ (1995) (E Beutler, MA Lichtman, BS Coller, TJ Kipps, eds.) *Williams Hematology*, Vol. 5: L22-6 McGraw-Hill, Inc.
46. Broxmeyer HE et al. (2003) *Proc Natl Acad Sci USA* 100(2): 645–50.
47. Guimond M et al. (2000) *Blood* 100(2): 375–82.
48. Itoh T et al. (2003) *Transfusion* 43(9): 1303–8.
49. Koliakos G et al. (2007) *Cytotherapy* 9(7): 654–9.
50. Rubinstein P et al. (1995) *Proc Natl Acad Sci USA* 92(22): 10119–22.
51. Slaper-Coretenbach IC et al. (1999) *Rheumatology* 38(8): 751–4.
52. Timeus F et al. (2003) *Haematologica* 88(1): 74–9.
53. Balducci E et al. (2003) *Stem Cells* 21(1): 33–40.
54. Ito CY et al. (2010) *Blood* 115(2): 257–60.
55. Douay L et al. (1986) *Exp Hematol* 14(5): 358–65.
56. Haas R et al. (1995) *Blood* 85(12): 3754–61.
57. Jagannath S et al. (1992) *Blood* 80(7): 1666–72.
58. Marit G et al. (1998) *Leukemia* 12(9): 1447–56.
59. Sagaster V et al. (2003) *Haematologica* 88(11): 1204–12.
60. Alonso JM 3rd et al. (2001) *Cytotherapy* 3(6): 429–33.
61. Spitzer G et al. (1980) *Blood* 55(2): 317–23.
62. Migliaccio AR et al. (2000) *Blood* 96(8): 2717–22.
63. Iori AP et al. (2004) *Bone Marrow Transplant* 33(11): 1097–105.
64. Yoo KH et al. (2007) *Bone Marrow Transplant* 39(9): 515–21.
65. Prasad VK et al. (2008) *Blood* 112(7): 2979–89.
66. Page KM et al. (2011) *Biol Blood Marrow Transplant* 17(9): 1362–74.
67. Frostad S et al. (1998) *Stem Cells* 16(5): 334–42.
68. Mayani H et al. (1993) *Blood* 81(12): 3252–8.
69. Qureshi SA et al. (1999) *Proc Natl Acad Sci USA* 96(21): 12156–61.
70. Schwartz GN et al. (1996) *Stem Cells* 14(3): 337–50.
71. Gribaldo L et al. (1999) *Exp Hematol* 27(11): 1593–8.
72. Pessina A et al. (2003) *Toxicol Sci* 75(2): 355–67.
73. Pessina A et al. (2001) *Toxicol In Vitro* 15(6): 729–40.
74. Volpe DA & Warren MK (2003) *Toxicol In Vitro* 17(3): 271–7.
75. Quintas-Cardama A et al. (2010) *Blood* 115(15): 3109–17.

76. Conneally E et al. (1996) *Blood* 87(2): 456–64.
77. Dobo I et al. (1995) *J Hematother* 4(4): 281–7.
78. Casati S et al. (2003) *Toxicol In Vitro* 17(1): 69–75.
79. Pessina A et al. (2009) *Toxicol In Vitro* 23(1): 194–200.
80. Parent-Massin D (2001) *Cell Biol Toxicol* 17(2): 87–94.
81. Van den Heuvel RL et al. (2001) *Cell Biol Toxicol* 17(2): 107–16.
82. Pyatt DW et al. (2000) *Mol Pharmacol* 57(3): 512–8.
83. Freund YR et al. (2002) *Toxicol Appl Pharmacol* 181(1): 16–26.
84. Froquet R et al. (2001) *Toxicol In Vitro* 15(6): 691–9.
85. Gribaldo L et al. (2000) *Toxicol Sci* 58(1): 96–101.
86. Gonzalez-Cid M et al. (2000) *Cell Biol Toxicol* 16(4): 235–41.

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, ALDEFLUOR, HetaSep, StemSpan, MyeloCult, MethoCult, STEMvision, SmartDish, STEMgrid, MegaCult, ImmunoCult, EasySep, RoboSep, RosetteSep, RapidSpheres, L-Calc, ErythroClear, ReproTeSR, STEMdiff, StemSep, and SepMate are trademarks of STEMCELL Technologies Canada Inc. Lymphoprep is a trademark of Serumwerk Bernburg AG. The products sold under the Lymphoprep brand name are also manufactured by Serumwerk Bernburg AG. Mozobil is a registered trademark of Genzyme Corporation. Neuopogen is a registered trademark of Amgen Inc. CryoStor, HypoThermosol and BloodStor are registered trademarks of BioLife Solutions. X-VIVO and HPGM are trademarks of Lonza Group Ltd. Alexa Fluor, Pacific Blue, and StemPro are trademarks of Life Technologies Inc. CellGro is a registered trademark of CellGenix GmbH. Stemline is a registered trademark of Sigma-Aldrich Co. LLC. Ficoll-Paque is a registered trademark of Cytiva. TeSR, E8, and mTeSR are trademarks of WARE. Corning and Falcon are registered trademarks of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.

UNLESS OTHERWISE STATED, PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES. FOR PRODUCT-SPECIFIC COMPLIANCE AND INTENDED USE INFORMATION, REFER TO THE PRODUCT INFORMATION SHEET. GENERAL INFORMATION ON QUALITY AT STEMCELL MAY BE FOUND AT WWW.STEMCELL.COM/COMPLIANCE.

HEMATOPOIETIC STEM AND PROGENITOR CELLS

Products for Your Research



TOLL FREE PHONE 1 800 667 0322

PHONE +1 604 877 0713

INFO@STEMCELL.COM

TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT WWW.STEMCELL.COM

DOCUMENT #29054 VERSION 8.0.0 APRIL 2025