

CELL ISOLATION PRODUCTS FOR
REGULATORY T CELLS

Cell Isolation Products For Regulatory T Cells





Table of Contents


3	Your First Step Towards Cutting-Edge Treg Research
4	Regulatory T Cells
5	The STEMCELL Advantage
6	Treg Isolation Procedure
7	Significant Time Savings
8	Classic Tregs
9	Highest Purity Tregs
10	Untouched Tregs
11	Isolate Tregs for Functional Analysis
12	Isolation of Mouse Tregs
13	Flexible Cell Isolation
14	References
15	Complete Product Listing

Your First Step Towards

Cutting-Edge Treg Research



Regulatory T cells (Tregs) are a specialized subset of T cells that play an essential role in the prevention of autoimmunity and the maintenance of immune system homeostasis and tolerance to self antigens. Harnessing the suppressive activity of Tregs for therapeutic purposes is a major area of interest as they hold great potential for the treatment of autoimmune diseases, prevention of transplant rejection and cancer immunotherapy.



The first step towards cutting-edge Treg research is the isolation of highly purified, functional Tregs. With unique cell isolation platforms, a team of dedicated and knowledgeable scientists, and a collaborative network of academic researchers, STEMCELL Technologies has developed a full range of products for the rapid and efficient isolation of highly functional Tregs. Choose the best Treg kit based on your desired phenotype, your downstream application and your starting sample, or let us design a kit that is tailor-made for your research. Start your Treg research right with high purity Tregs and STEMCELL Technologies.

Regulatory T Cells

Regulatory T cells (Tregs) have the ability to suppress T cell responses and play a critical role in peripheral tolerance and immune regulation.¹ In humans, natural Tregs are phenotypically CD4⁺CD25^{high}, and express high levels of the intracellular transcription factor FoxP3.^{2,3} Unlike their mouse counterparts, human Tregs appear to display a heterogeneous cell surface phenotype and regulatory capacity, leading to the identification and classification of various Treg subtypes within the CD4⁺CD25^{high}FoxP3⁺ T cell population.⁴ Human peripheral blood Tregs are rare, comprising only 2 – 3% of the CD4⁺ T cell population.² For these reasons, the identification and isolation of human Tregs can be challenging and time-consuming.

To overcome these technical difficulties, researchers have exploited the expression of other surface markers, either alone or in combination, in order to enrich for distinct Treg populations. The expression of CD127 (IL-7 receptor alpha chain) inversely correlates with FoxP3 expression,^{6,7} while CD49d, the α -chain of the integrin VLA-4 ($\alpha_4\beta_1$) is expressed on a majority of pro-inflammatory effector cells.⁸ Removal of CD127^{high} cells within the CD4⁺CD25^{high} population enriches for a subset of Tregs with high levels of FoxP3 expression and high expansion potential,⁷ whereas the use of CD127 and CD49d antibodies in combination provides a negative selection (e.g. depletion) approach to isolate purified untouched Tregs with suppressive activity.^{7,8}

More recently, new subpopulations of Tregs have been identified based on the expression of CD45RA.^{4,9,10} The heterogeneous expression of CD45RA allows for the identification and isolation of Treg subtypes that are functionally distinct and have different suppression and expansion capabilities. These subtypes include resting Tregs, activated Tregs, and cytokine-secreting Tregs.

In order to support and facilitate research in this exciting field, STEMCELL Technologies has developed cutting-edge products for the isolation and assessment of Tregs. Our Treg isolation kits provide the fastest and easiest way to isolate various subtypes of Tregs directly from human whole blood, buffy coat, PBMC or from mouse spleen with high purity and recovery. Purified cells are highly functional and ready to use in any downstream Treg assay.

The STEMCELL Advantage



FAST

- Save hours on a complete separation
- Isolate cells directly from whole blood, buffy coat or leukapheresis samples

COST-EFFECTIVE

- No expensive columns
- Save time on costly flow cytometry sorting
- Kits are stable for up to two years

VERSATILE

- Isolate different subtypes of Tregs
- Choose the best kit for your downstream application or starting sample
- Customize kits to get the Tregs you want

Did you know?

EasySep™ cell isolation can be fully automated with RoboSep™, which can isolate from four samples simultaneously.

Treg Isolation Procedure

Isolate Tregs directly from whole blood or buffy coat using our **Complete Kits for Regulatory T Cells**. Tregs are first pre-enriched during a Ficoll-Paque™ PLUS spin using RosetteSep™ (Step 1). The CD25^{high} cells are then positively selected using EasySep™ immunomagnetic separation system (Step 2).

For isolation of CD4⁺CD25^{high} Tregs from PBMCs, Tregs are pre-enriched with EasySep™ instead of RosetteSep™, followed by the EasySep™ CD25^{high} positive selection procedure (Step 2).

1 Pre-enrichment of Cells





FIGURE 1. RosetteSep™ Tetrameric Antibody Complexes (TAC)

2 Selection of CD25^{high} Tregs





FIGURE 2. EasySep™ Tetrameric Antibody Complexes (TAC)

FIGURE 1. Unwanted cells are crosslinked to red blood cells (RBC) using bispecific Tetrameric Antibody Complexes (TACs). These labeled cells pellet with the RBC when centrifuged over density medium, leaving the enriched CD4⁺ T cells at the plasma: density medium interface.

FIGURE 2. TACs selectively bind to desired cells (positive selection) or unwanted cells (negative selection) and crosslink them to magnetic particles. The labeled cells are then separated from unlabeled cells using the EasySep™ magnet.

Significant Time Savings

Protocol Comparison Between Treg Isolation Kits by STEMCELL Technologies and a Competitor

	STEMCELL TECHNOLOGIES	COMPETITOR
Time	Complete Kits for Human Regulatory T Cell Isolation	Column-based Human Regulatory T Cell Isolation Kits
0:00	RosetteSep™ CD4 ⁺ or CD4 ⁺ subset cocktail incubation	
0:15		Ficoll-Paque™ PLUS spin
0:30	Ficoll-Paque™ PLUS spin	
0:45		Wash cells x 2
1:00	Wash cells x 2	
1:15		CD4 ⁺ biotin cocktail incubation
1:30	Resuspend in buffer	Microbead incubation
1:45	CD25 cocktail incubation	Centrifuge
2:00	Magnetic particle incubation	Load and run column
2:15	Magnetic separation	Wash column x 2
2:30		
2:45	Wash cells	Wash cells
3:00	Tregs ready for use	CD25 microbead incubation
3:15		Wash cells
3:30		Load and run column
3:45		Wash column x 3
4:00		Load and run column
4:15		Wash column x 3
4:30		Wash cells
4:45		Tregs ready for use

3 hr
No columns
4 centrifugation steps

4 hr 45 min
3 columns
9 centrifugation steps

Classic Tregs

Initial phenotypic analysis of human Tregs showed that in vitro regulatory activity is restricted to CD4⁺ T cells that express high levels of CD25.^{2,3} Isolation of CD4⁺CD25^{high} T cells allows researchers to isolate and study the broadest phenotypic range of Tregs.


STEMCELL Technologies' kits provide the **fastest** isolation of highly purified CD4⁺CD25^{high} T cells. These cells can be used either directly in downstream assays, or as a pre-enriched population for reducing flow sorting time.

Complete Kit for Human CD4⁺CD25^{high} T Cells (CATALOG #15862)

Starting with whole blood, isolate CD4⁺CD25^{high} T cells quickly and easily in two simple steps.

HOW IT WORKS:

CD4⁺ T cells are enriched by negative selection using RosetteSep™ followed by positive selection of CD25^{high} cells using EasySep™ or RoboSep™ (see page 6).




Purity and phenotype of CD3⁺CD4⁺ cells following pre-enrichment with RosetteSep™ ("Pre-enriched") and positive selection with EasySep™ ("Enriched").

EasySep™ Human CD4⁺CD25^{high} T Cell Isolation Kit (CATALOG #18062)

Starting with PBMC, isolate CD4⁺CD25^{high} T cells without columns or washes.

HOW IT WORKS:

CD4⁺ T cells are enriched by negative selection using EasySep™ or RoboSep™ followed by positive selection of CD25^{high} cells using EasySep™ or RoboSep™.



Purity and phenotype of CD3⁺CD4⁺ cells following pre-enrichment with EasySep™ ("Pre-enriched") and positive selection with EasySep™ ("Enriched").

Highest Purity Tregs


CD127 and CD49d are useful cell surface markers for removing effector cells from CD4⁺CD25⁺ T cell populations allowing for further enrichment of highly suppressive Tregs.^{6,7,8} STEMCELL has developed kits using these markers to give you the highest purity of FoxP3⁺ cells for downstream analysis without the need for additional flow sorting.

Complete Kit for Human CD4⁺CD127^{low}CD25^{high} Regulatory T Cells (CATALOG #15861)

Get high purity and recovery of FoxP3⁺ Tregs directly from whole blood.

HOW IT WORKS:

CD4⁺CD127^{low} T cells are enriched by negative selection using RosetteSep™ followed by positive selection of CD25^{high} cells using EasySep™ or RoboSep™ (see page 6).




Purity and phenotype of CD3⁺CD4⁺ cells following pre-enrichment with RosetteSep™ ("Pre-enriched") and positive selection with EasySep™ ("Enriched").

Complete Kit for Human CD4⁺CD127^{low}CD49d⁺CD25^{high} Regulatory T Cells (CATALOG #15864)

Get the highest purity of FoxP3⁺ activated/memory Tregs directly from whole blood.

HOW IT WORKS:

CD4⁺CD127^{low}CD49d⁺ T cells are enriched by negative selection using RosetteSep™ followed by positive selection of CD25^{high} cells using EasySep™ or RoboSep™ (see page 6).



Purity and phenotype of CD3⁺CD4⁺ cells following pre-enrichment with RosetteSep™ ("Pre-enriched") and positive selection with EasySep™ ("Enriched").

Untouched Tregs

For researchers who prefer a gentle method of isolating untouched Tregs, STEMCELL has developed negative selection kits that enrich Tregs without the need for positive selection. Our EasySep™ isolation kits give highly functional, untouched Tregs in a single step from PBMCs in one hour.


EasySep™ Human CD4⁺CD127^{low}CD49d⁻ Regulatory T Cells (CATALOG #19232)

Isolate untouched FoxP3⁺ T cells in a single step from PBMC.


HOW IT WORKS:

CD4⁺CD127^{low}CD49d⁻ T cells are enriched by EasySep™ negative selection, which depletes unwanted cells and leaves desired cells untouched.

Start:
Gate: CD3⁺CD4⁺ cells



Enriched:
Gate: CD3⁺CD4⁺ cells



Purity and phenotype of CD3⁺CD4⁺ cells before ("Start") and after ("Enriched") EasySep™.


EasySep™ Human CD4⁺CD127^{low} T Cells (CATALOG #19231)

Isolate untouched CD4⁺CD127^{low} T cells in a single step from PBMC and save time on flow sorting.


HOW IT WORKS:

CD4⁺CD127^{low} T cells are enriched by EasySep™ negative selection, which depletes unwanted cells and leaves desired cells untouched.

Start:
Gate: CD3⁺CD4⁺CD127^{low} cells



Enriched:
Gate: CD3⁺CD4⁺CD127^{low} cells




Purity and phenotype of CD3⁺CD4⁺ cells before ("Start") and after ("Enriched") EasySep™.

Isolate Tregs for Functional Analysis

The suppression activity of isolated CD4⁺CD25^{high} T cells was assessed by measuring their impact on CD4⁺CD25^{neg} proliferative responses to stimulation with anti-CD3/CD28 coated beads in the absence of IL-2 (Figure 3). Mixing purified CD4⁺CD25^{high} T cells with CD4⁺CD25^{neg} T cells at a ratio of 1:1 almost completely eliminates detectable CD4⁺CD25^{neg} T cell proliferation responses, and suppression of proliferation continues to be detected when cells are mixed in a ratio of 1:10.

Next, Tregs isolated with RosetteSep™ and EasySep™ were first expanded ex vivo for 14 days in the presence of rapamycin, and then co-cultured with autologous CD4⁺CD25^{depleted} responder T cells labelled with CFSE at a 1 Treg : 2 responders ratio (Figure 4). Whereas, CD4⁺CD25^{depleted} T cells expanded in the absence of rapamycin have no suppressive capacity, all of the expanded Tregs effectively suppressed proliferation of the autologous CD4⁺CD25^{depleted} responder cells.

FIGURE 3. Isolated CD4⁺CD25^{high} T cells are anergic and suppress proliferation of CD4⁺CD25^{neg} T cells




Purified CD4⁺CD25^{high} T cell fractions were assessed for anergy by measuring their proliferation response to anti-CD3/CD28 coated beads (in grey). The suppression activity of purified CD4⁺CD25⁺ T cells was assessed by measuring their ability to reduce the proliferative response of CD4⁺CD25^{neg} T cells to CD3/CD28 beads (in black). T cell proliferation was quantified by measuring dilution of the fluorescent dye CFSE with flow cytometry. Results are expressed as proliferation detected relative to control CD4⁺CD25^{neg} cell populations stimulated with CD3/CD28 beads for 7 days in RPMI medium containing 5% human AB serum. The bars are the means of 5 experiments with error bars indicating standard deviations.

Acknowledgements:

We would like to thank Dr. Rajendra Pahwa (Diabetes Research Institute, University of Miami Miller School of Medicine) for samples and technical assistance during the development of some RoboSep™ protocols.

FIGURE 4. Expanded Tregs Suppress Proliferation of CD4⁺CD25^{depleted} Responder Cells



Left: Freshly isolated Tregs were stimulated with anti-CD3/CD28 coated beads in the presence of exogenous IL-2 and 100 nm rapamycin. On day 7 the cells were restimulated and the cultures continued until day 14. The CD4⁺CD25⁻ depleted cells in grey, were cultured in the absence of rapamycin for the duration of culture.

Right: Following the 14 days of culture, the expanded Tregs were co-cultured with autologous CD4⁺CD25⁻ depleted responder T cells labelled with CFSE at a 1 Treg : 2 responders ratio.

STEMCELL's CD4⁺CD25^{high} Regulatory T Cell Isolation Kit and related products have been isolated for:

- Activation and expansion^{1,2,3,4}
- Suppression of allogeneic mixed lymphocyte response (MLR)^{1,5}
- Co-culture with dendritic cells^{2,5,6}
- Adoptive transfer^{4,7}
- RNA extraction^{3,8}
- Gene expression profile⁸
- Real-time PCR analysis^{8,5}
- Chromatin immunoprecipitation–DNA microarray (ChIP-chip)⁸

Isolation of Mouse Tregs

Mouse regulatory T cells can be isolated from single cell suspensions of spleen or other tissues using EasySep™ Mouse CD4⁺CD25⁺ Regulatory T Cell Isolation Kit (Catalog #19782). Mouse CD4⁺ T cells are enriched using EasySep™ Mouse CD4⁺ T Cell Pre-Enrichment Kit (Catalog #19772), which depletes cells expressing CD8, CD11b, CD19, CD24, CD45R, and CD49b. CD25⁺ cells are labeled using EasySep™ anti-CD25 PE Labeling Reagent, PE Selection Cocktail, and EasySep™ Magnetic Particles, and separated from unlabeled cells using an EasySep™ magnet.

FIGURE 5. Typical FACS Analysis of EasySep™ Mouse CD4⁺CD25⁺ Regulatory T Cell Isolation Kit






FIGURE 6. EasySep™ isolated CD4⁺CD25⁺ T cell fractions are highly enriched for FoxP3-expressing cells.



Intracellular FoxP3 expression was detected in 4% paraformaldehyde-fixed CD4⁺CD25⁺ T cells using a FITC-labeled anti-FoxP3 antibody. Similar results were obtained in 4 separate experiments.

FIGURE 7. CD4⁺CD25⁺ T cells isolated using EasySep™ are able to suppress T cell proliferation responses.



EasySep™ isolated CD4⁺CD25⁺ T cells were co-cultured at a ratio of 1:1 and 1:3 with stimulated CD4⁺CD25^{neg} T cells to determine their ability to suppress T cell proliferation. Results are expressed relative to the level of proliferation detected in stimulated CD4⁺CD25^{neg} cell cultures. The bars are the means of 4 experiments with error bars indicating standard deviations.

Flexible Cell Isolation

Pre-Enrichment Kits

When purities of over 90% are required, STEMCELL's Pre-Enrichment Kits are an ideal way to save hours on flow sorting while getting purities you need for that crucial experiment. RosetteSep™ or EasySep™ Pre-Enrichment Kits provide excellent yields of pre-enriched Tregs from whole blood, buffy coat or PBMC, and allow you to go from sample to sorting in less than 90 minutes.

STARTING SAMPLE	PRODUCT	CATALOG #
Whole blood	RosetteSep™ Human CD4 ⁺ T Cell Enrichment Cocktail	15062
Whole blood	RosetteSep™ Human CD4 ⁺ CD127 ^{low} T Cell Enrichment Cocktail	15361
PBMC	EasySep™ Human CD4 ⁺ CD127 ^{low} T Cell Enrichment Kit	19231
PBMC	EasySep™ Human CD4 ⁺ T Cell Enrichment Kit	19052
Single cell suspensions	EasySep™ Mouse CD4 ⁺ T Cell Pre-Enrichment Kit	19772

Customizable Cell Isolation

STEMCELL's Treg isolation kits are easily customizable. If you can't find the kit you need, we can help you find the antibodies and method to isolate the Treg subtype that you want. We have a team of technical experts who will work with you to help you develop the right cell isolation solution to meet your Treg research needs.

Did you know?

STEMCELL has made-to-order kits for the isolation of Tregs from rat and non-human primates. Contact techsupport@stemcell.com for more information.

References

1. Sakaguchi S. Naturally arising Foxp3-expressing CD25⁺CD4⁺ regulatory T cells in immunological tolerance to self and non-self. *Nat Immunol.* 2005; 6(4): 345-352.
2. Kang SM, Tang Q and Bluestone JA. CD4⁺CD25⁺ Regulatory T Cells in Transplantation: Progress, Challenges and Prospects. *American Journal of Transplantation.* 2007; 7: (6); 1457–1463.
3. Baecher-Allan C, Wolf E, and Hafler DA. Functional analysis of highly defined, FACS-isolated populations of human regulatory CD4⁺CD25⁺ T cells. *Clin Immunol.* 2005; 115(1): 10-18.
4. Baecher-Allan C, Brown JA, Freeman GJ, and Hafler DA. CD4⁺CD25^{high} Regulatory cells in human peripheral blood. *J Immunol;* 2001 1: 167(3): 1245-1253.
5. Miyara M, Yoshiokay, Kitoh A et al. Functional delineation and differentiation dynamics of human CD4⁺ T cells expressing the FoxP3 transcription factor. *Immunity.* 2009; 30(6): 899-911.
6. Liu W, Putnam AL, Xu-Yu Z et al. CD127 expression inversely correlates with FoxP3 and suppressive function of human CD4⁺ T reg cells. *J Exp Med.* 2006; 203: 1701–1711.
7. Seddiki N, Santner-Nanan B, Martinson J et al. Expression of interleukin (IL)-2 and IL-7 receptors discriminates between human regulatory and activated T cells. *J Exp Med.* 2006; 203(7): 1693-1700.
8. Kleinewietfeld M, Starke M, Di Mitri D et al. CD49d provides access to "untouched" human Foxp3⁺ Treg free of contaminating effector cells. *Blood.* 2009; 113(4): 827-836.
9. Seddiki N, Santner-Nanan B, Tangye SG et al. Persistence of naive CD45RA⁺ regulatory T cells in adult life. *Blood.* 2006; 107(7), 2830–2838.
10. Hoffman P, Eder R, Boeld TJ, et al. Only the CD45RA⁺ subpopulation of CD4⁺CD25^{high} T cells gives rise to homogeneous regulatory T-cell lines. *Blood.* 2006; 108(13): 4260-4267.

Selected Product References

1. Wang X-N, Haniffa MA, Holtick U, et al. Regulatory T-Cell Suppression of CD8⁺ T-Cell-Mediated Graft-Versus-Host Reaction Requires Their Presence During Priming. *Transplantation.* 88(2): 188-197.
2. Nguyen KD, Vanichsarn C, and Nadeau KC. Impaired IL-10-dependent Induction of Tolerogenic Dendritic Cells by CD4⁺CD25^{hi}CD127^{lo} Natural Regulatory T Cells in Human Allergic Asthma. *J Respir And Crit Care Med.* 2009; 180: 823-833.
3. Venken K, Hellings N, Broekmans T et al. Natural Naïve CD4⁺CD25⁺CD127^{low} Regulatory T Cell (Treg) Development and Function Are Disturbed in Multiple Sclerosis Patients: Recovery of Memory Treg Homeostasis during Disease Progression. *J Immunol* 2008; 180: 6411-6420.
4. Trzonkowski P, Bieniaszewska M, Juciska et al. First-in-man clinical results of the treatment of patients with graft versus host disease with human ex vivo expanded CD4⁺CD25⁺CD127⁺ T regulatory cells. *Clinical Immunology.* 2009; 133(1): 22-26.
5. Shalev I, Liu H, Kosciuk C, et al. Targeted Deletion of fgl2 Leads to Impaired Regulatory T Cell Activity and Development of Autoimmune Glomerulonephritis. *J Immunol* 2008; 180: 249-260.
6. Muthuswamy R, Urban J, Lee J, Reinhart TA, Bartlett D, and Kalinski P. Ability of Mature Dendritic Cells to Interact with Regulatory T Cells Is Imprinted during Maturation. *Cancer Res.* 2008; 68(14): 5972-5978.
7. Richer MJ, Straka N, Fang D, Shanina I and Horwitz MS. Regulatory T-Cells Protect From Type 1 Diabetes After Induction by Coxsackievirus Infection in the Context of Transforming Growth Factor. *Diabetes.* 2008; 57(5): 1302-1311.
8. Liu w, Putnam AL, Xu-yu Z, et al. CD127 expression inversely correlates with FoxP3 and suppressive function of human CD4⁺ T reg cells. *J Exp Med* 2006; 203 (7): 1701-1711.

Complete Product Listing

Cell Isolation Kits

SPECIES	STARTING SAMPLE	PHENOTYPE OF CELLS	PRODUCT	CATALOG #	PAGE
Human	Whole blood	CD4 ⁺ CD25 ⁺ T cells	Complete Kit for Human CD4 ⁺ CD25 ⁺ T Cells	15862	8
		CD4 ⁺ CD127 ^{low} CD25 ⁺ T cells	Complete Kit for Human CD4 ⁺ CD127 ^{low} CD25 ⁺ Regulatory T Cells	15861	9
		CD4 ⁺ CD127 ^{low} CD49d ⁺ CD25 ⁺ T cells	Complete Kit for Human CD4 ⁺ CD127 ^{low} CD49d ⁺ CD25 ⁺ Regulatory T Cells	15864	9
	PBMC	CD4 ⁺ CD25 ⁺ T cells	EasySep™/RoboSep™ Human CD4 ⁺ CD25 ⁺ T Cell Isolation Kit	18062	8
		CD4 ⁺ CD127 ^{low} T cells	EasySep™/RoboSep™ Human CD4 ⁺ CD127 ^{low} T Cell Enrichment Kit	19231	10
		CD4 ⁺ CD127 ^{low} CD49d ⁺ T cells	EasySep™/RoboSep™ Human CD4 ⁺ CD127 ^{low} CD49d ⁺ Regulatory T Cell Enrichment Kit	19232	10
Mouse	Spleen or other tissues	CD4 ⁺ CD25 ⁺ T cells	EasySep™/RoboSep™ Mouse CD4 ⁺ CD25 ⁺ T Cell Isolation Kit	19782	12

Pre-Enrichment Kits

SPECIES	STARTING SAMPLE	PRODUCT	CATALOG #	PAGE
Human	Whole blood	RosetteSep™ Human CD4 ⁺ T Cell Enrichment Cocktail	15062	13
		RosetteSep™ Human CD4 ⁺ CD127 ^{low} T Cell Enrichment Cocktail	15361	13
	PBMC	EasySep™ Human CD4 ⁺ CD127 ^{low} T Cell Enrichment Kit	19231	13
		EasySep™ Human CD4 ⁺ T Cell Enrichment Kit	19052	13
Mouse	Spleen or other tissues	EasySep™ Mouse CD4 ⁺ T Cell Enrichment Kit	19772	13

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL-FREE PHONE 1 800 667 0322 • PHONE 1 604 877 0713

TOLL-FREE FAX 1 800 567 2899 • FAX 1 604 877 0704

ORDERS@STEMCELL.COM • INFO@STEMCELL.COM

FOR FULL CONTACT DETAILS WORLDWIDE VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT FOR THERAPEUTIC OR DIAGNOSTIC USE. CATALOG #29976 VERSION 1.1.0 NOVEMBER 2011

