

BIT 9500 Serum Substitute

Supplement for serum-free cell culture media

Catalog #09500

100 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

BIT 9500 Serum Substitute has been developed for use in applications where a serum-free culture medium of defined composition is required. This product contains bovine serum albumin (BSA), insulin, and transferrin (BIT) in Iscove's Modified Dulbecco's Medium (IMDM). It contains pre-screened batches of BSA that have been selected to support the optimal growth of human hematopoietic progenitor cells in serum-free media formulations. It is also suitable for the culture of mouse hematopoietic progenitor cells in serum-free conditions.

Properties

Storage: Store at -20°C.

Shelf Life: Stable until expiry date (EXP) on label.

Contains:

- Bovine serum albumin
- Recombinant human insulin
- Human transferrin (iron-saturated)
- Iscove's MDM

This product contains material derived from human plasma. Donors have been tested and found negative for HIV-1 and -2, hepatitis B, and hepatitis C prior to donation. However, this product should be considered potentially infectious and treated in accordance with universal handling precautions.

Directions for Use

1. Thaw at room temperature (15 - 25°C) or overnight at 2 - 8°C.
2. Swirl the bottle to mix.
3. Store at 2 - 8°C for up to 1 month. Alternatively, aliquot and store at -20°C. After thawing aliquots, store at 2 - 8°C for up to 1 month; do not re-freeze. Do not exceed the shelf life of the supplement.

NOTE: Use BIT 9500 Serum Substitute at a final concentration of 20%. For example, add 20 mL of BIT 9500 Serum Substitute to 80 mL of medium plus supplements, for a final volume of 100 mL.

Related Products

For related products, including specialized culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com/HSPCworkflow. For available fresh and cryopreserved peripheral blood, cord blood, and bone marrow products, visit www.stemcell.com/primarycells.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2026 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.