

# STEMprep™ Sample Tubes



Scientists Helping Scientists™ | [www.stemcell.com](http://www.stemcell.com)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #200-0800

25 Units

## Product Description

STEMprep™ Sample Tubes are specially designed for use with the STEMprep™ Tissue Dissociator instrument. The instrument drives an internal rotor in each STEMprep™ Sample Tube, enabling robust and efficient dissociation of a wide variety of tissue types.


## Properties

<b>Stability and Storage:</b>	Store at 15 - 25°C. Stable until expiry date (EXP) on label.
<b>Usage:</b>	Single use only
<b>Tissue Size:</b>	20 - 4000 mg
<b>Min/Max Volume:</b>	1 - 10 mL NOTE: The total combined sample volume should not exceed the "MAX" line indicated on the tube.
<b>Maximum Rotational Speed:</b>	3000 RPM in forward direction; 40 RPM in reverse direction
<b>Maximum Total Revolutions:</b>	4000 revolutions

## Directions for Use

STEMprep™ Sample Tubes come fully assembled. Each tube consists of a lid, a rotor, and a reservoir with blades (Figure 1). To open, turn the lid counterclockwise to separate the lid and rotor from the reservoir. Secure the lid by turning it clockwise, past the initial resistance, until it stops firmly and cannot be tightened any further.

For complete instructions on using STEMprep™ Sample Tubes, refer to the Technical Manual: STEMprep™ Tissue Dissociator (Document #10000030598), available at [www.stemcell.com](http://www.stemcell.com), or contact us to request a copy.



- A. Lid
- B. Rotor
- C. Reservoir
- D. Blades

**Figure 1. STEMprep™ Sample Tube**

## Notes and Tips

The maximum amount of sample size is dependent on the STEMprep™ protocol. For optimal results and to avoid damage to the tubes, it is recommended to pre-cut the tissue into pieces (< 500 mg) before loading it into the Sample Tubes.

The tubes may be damaged if hard material, such as bone, is used.

If the lid and rotor become separated, the rotor can be reinserted into the lid. However, disassembling the lid from rotor is not recommended, as it may compromise sealing performance and could result in sample loss.

## Related Products

For more information about STEMprep™ kits and protocols, visit [www.stemcell.com/stemprep](http://www.stemcell.com/stemprep), contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO [WWW.STEMCELL.COM/COMPLIANCE](http://WWW.STEMCELL.COM/COMPLIANCE).

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and STEMprep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.