

BloodStor® 55-5

Optimized biopreservation reagent for hematopoietic cells and tissues

Catalog # 07943 6 x 8 mL Syringes
07937 16 x 7.2 mL
07950 1000 mL Bag



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

BloodStor® is a series of generic cGMP freeze media products used to cryopreserve stem cells and other cells isolated from umbilical cord blood, peripheral blood, bone marrow, and other biologics. BloodStor® 55-5 is pre-formulated with 55% dimethyl sulfoxide (DMSO) USP, 5% Dextran-40 USP, and water for injection (WFI) quality water.

- Ready-to-use
- Animal component-free
- Sterile, USP vial
- cGMP manufactured with USP grade / highest quality components
- Sterility, endotoxin, and cell-based quality control testing

Product Information

PRODUCT NAME	CATALOG #	SIZE	STORAGE	SHELF LIFE	CONTAINS
Blood Stor® 55-5 Syringe Bundle	07943	6 x 8 mL Syringes	Store at 2 - 8°C.	Stable for 24 months from date of manufacture (MFG) on label. Protect product from prolonged exposure to light.	55% DMSO and 5% Dextran-40
Blood Stor® 55-5	07937	16 x 7.2 mL	Store at 2 - 8°C.	Stable until expiry date (EXP) on label. Protect product from prolonged exposure to light.	55% DMSO and 5% Dextran-40
Blood Stor® 55-5	07950	1000 mL	Store at 2 - 8°C.	Stable for 24 months from date of manufacture (MFG) on label. Protect product from prolonged exposure to light.	55% DMSO and 5% Dextran-40

Please refer to the Safety Data Sheet (SDS) for hazard information.

Product may be shipped at room temperature (15 - 25°C) and should be refrigerated upon receipt.

Handling / Directions for Use

FREEZING

1. Wipe down the outside of the BloodStor® 55-5 container with 70% ethanol or isopropanol before opening.
2. Obtain a cell suspension using a red cell depletion or reduction protocol.
3. Add BloodStor® 55-5 at a ratio of 1 part BloodStor® 55-5 to 5 parts sample (e.g. 2 mL BloodStor® 55-5 to 10 mL sample).
4. Mix thoroughly and transfer the suspension to a cryovial.
5. Freeze cells using a standard slow rate-controlled cooling protocol (approximately -1°C/minute) or an isopropanol freezing container and store at liquid nitrogen temperature (-135°C).

NOTE: Long-term storage at -80°C is not recommended.

THAWING

1. Quickly thaw cells in a 37°C water bath by gently shaking the vial. Do not submerge the vial. Remove the vial when only a small frozen cell pellet remains. Do not vortex cells.
2. Wipe the outside of the vial with 70% ethanol or isopropanol.
3. Dilute in cell culture medium of choice at a ratio of 1 part sample in 10 parts medium.
4. Plate cells immediately according to desired protocol. If a smaller volume is required, centrifuge the cells appropriately and resuspend the cells in a volume of culture medium appropriate for the desired plating vessel.

THIS PRODUCT IS MANUFACTURED UNDER A cGMP QUALITY MANAGEMENT SYSTEM COMPLIANT TO 21 CFR 820.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. BloodStor is a registered trademark of BioLife Solutions. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.