

0.2% Heparin Sodium Salt in PBS

Cell Culture Supplement

Catalog # 07980

2 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Heparin is a mucopolysaccharide with anticoagulant properties. Heparin supports the binding of Fibroblast Growth Factor (FGF) to its receptor and increases the stability of FGF. Heparin is used together with Epidermal Growth Factor (EGF) and FGF in the culture of embryonic rat or adult human neural stem and progenitor cells, and adult mouse neural stem and progenitor cells. Heparin and FGF are not required for embryonic mouse neural stem and progenitor cell culture.

Properties

Storage: Store at 2 - 8°C.

Shelf Life: Stable until expiry date (EXP) on label.

Contains: 0.2% (2 mg/mL; 1000 IU/mL) heparin sodium salt in phosphate-buffered saline (PBS)

Handling / Directions For Use

Use 0.2% Heparin Sodium Salt in PBS with the NeuroCult™ NS-A Proliferation Kits (Human; Catalog #05751 and Rat; Catalog #05771) for the growth and expansion of neural stem and progenitor cells. This product is added to the culture medium together with recombinant human Basic Fibroblast Growth Factor (rh bFGF; Catalog #02634). Add 10 µL of 0.2% Heparin Sodium Salt in PBS to 10 mL of medium to achieve a final concentration of 2 µg/mL (0.0002%).

For complete instructions refer to the Technical Manual: In Vitro Proliferation and Differentiation of Neural Stem and Progenitor Cells Using NeuroCult™ (Human; Document #28724 or Rat; Document #28725), available on our website at www.stemcell.com or contact us to request a copy.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, NeuroCult and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.