

StemSpan™ SFEM



Serum-free medium for culture and expansion of hematopoietic cells

Catalog # 09600 100 mL
 09650 500 mL

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

StemSpan™ Serum-Free Expansion Medium (SFEM) has been developed and tested for the in vitro culture and expansion of human hematopoietic cells. This medium contains pre-tested bovine serum albumin, insulin and transferrin, and supplements in Iscove's MDM. Recombinant hematopoietic growth factors, required for the optimal growth and expansion of hematopoietic cells, have not been added to StemSpan™ SFEM. This allows users the flexibility to prepare medium that meets their requirements.

In addition to the culture and expansion of human hematopoietic cells, StemSpan™ SFEM has been used for the culture and expansion of hematopoietic stem and progenitor cells from other species, including mouse, non-human primate, and dog (Miller & Eaves; Zhang & Lodish; Sandrin et al.; Bauer et al.). StemSpan™ SFEM has also been used for the generation, culture, and assay of human and mouse dendritic cells, and for the culture of various other hematopoietic and non-hematopoietic cell types. In most applications, addition of specific hematopoietic growth factors, cytokines, and/or other compounds is required for optimal growth.

Properties

- Storage:** Store at -20°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- Iscove's MDM
 - Bovine serum albumin
 - Recombinant human insulin
 - Human transferrin (iron-saturated)
 - 2-Mercaptoethanol
 - Supplements

Donors have been tested and found negative for hepatitis B surface antigen (HBsAg) and HIV-1 antibodies and/or HIV-1 antigen. However, this product should be considered potentially infectious and treated in accordance with universal handling precautions.

Handling / Directions For Use

1. Thaw StemSpan™ SFEM at room temperature (15 - 25°C) or overnight at 2 - 8°C. Mix well.
NOTE: If not used immediately, aliquot into tubes and store at -20°C. Once aliquots are thawed, do not re-freeze.
2. Add desired cytokines, growth factors, and other components to StemSpan™ SFEM and mix well.
NOTE: Added components and cells in sterile cell culture medium (e.g. Iscove's MDM or DMEM) should not exceed ~10% of total volume.
3. Add cells, mix well, and set up cultures as desired.

ASSESSMENT OF HEMATOPOIETIC CELLS

Assessment of CD34+ cells before and after culture may be performed by flow cytometry using the following fluorochrome-conjugated antibody clones:

- Anti-Human CD34 Antibody, Clone 581 (Catalog #60013) or Clone 563 (Catalog #60119) or Clone 8G12 (Catalog #60121)
- Anti-Human CD45 Antibody, Clone HI30 (Catalog #60018) or Clone 2D1 (Catalog #60123)

Notes and Tips

STEMCELL Technologies Inc. recommends the use of Human LDL (Catalog #02698) as a culture supplement. It has been pre-screened for the culture, expansion, and colony assay of human hematopoietic and non-hematopoietic cells in serum-free culture media. It promotes the proliferation and survival of human hematopoietic and other progenitor cells in culture, resulting in increased cell output in expansion cultures and increased colony numbers and/or colony size in colony assays.

Selection of an optimal cytokine combination is dependent upon the source and type of cells and the experimental objectives of the researcher. StemSpan™ expansion supplements, described below, are suitable for use with StemSpan™ SFEM.

- StemSpan™ CD34+ Expansion Supplement (10X) (Catalog #02691)
 - Culture and expansion of large numbers of human CD34+ progenitor cells
 - Contains: rh SCF, rh TPO, rh IL-3, rh IL-6, rh Flt3 ligand, other additives
- StemSpan™ CC100 (Catalog #02690)
 - Culture and expansion of human hematopoietic cells
 - Contains: rh Flt3 ligand, rh SCF, rh IL-3, rh IL-6
- StemSpan™ CC110 (Catalog #02697)
 - Culture and expansion of human hematopoietic cells
 - Contains: rh Flt3 ligand, rh SCF, rh TPO
- StemSpan™ Erythroid Expansion Supplement (100X) (Catalog #02692)
 - Expansion and lineage-specific differentiation of human CD34+ cells into erythroid progenitor cells
 - Contains: rh SCF, rh IL-3, rh EPO
- StemSpan™ Megakaryocyte Expansion Supplement (100X) (Catalog #02696)
 - Expansion and lineage-specific differentiation of human CD34+ cells into megakaryocyte progenitor cells
 - Contains: rh SCF, rh TPO, rh IL-6, rh IL-9
- StemSpan™ Myeloid Expansion Supplement (100X) (Catalog #02693)
 - Expansion and lineage-specific differentiation of human CD34+ cells into myeloid progenitor cells
 - Contains: rh SCF, rh TPO, rh G-CSF, rh GM-CSF

SCF = stem cell factor; EPO = erythropoietin; TPO = thrombopoietin; rh = recombinant human; IL = interleukin; Flt = fms-like tyrosine kinase

RELATED PRODUCTS

For related products, including specialized culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com/HSPCworkflow or contact us at techsupport@stemcell.com. For available fresh and cryopreserved peripheral blood, cord blood, and bone marrow products in your region, visit www.stemcell.com/primarycells.

References

- Abdelwahab SF et al. (2003) HIV-1-suppressive factors are secreted by CD4+ T cells during primary immune responses. *Proc Natl Acad Sci U S A* 100(25): 15006–10.
- Agosti V et al. (2004) Critical role for Kit-mediated Src kinase but not PI 3-kinase signaling in pro T and pro B cell development. *J Exp Med* 199(6): 867–78.
- Baksh D et al. (2005) Soluble factor cross-talk between human bone marrow-derived hematopoietic and mesenchymal cells enhances in vitro CFU-F and CFU-O growth and reveals heterogeneity in the mesenchymal progenitor cell compartment. *Blood* 106(9): 3012–9.
- Bauer TR et al. (2006) Correction of the disease phenotype in canine leukocyte adhesion deficiency using ex vivo hematopoietic stem cell gene therapy. *Blood* 108(10): 3313–20.
- Gotze K et al. (2001) gp130-Stimulating designer cytokine hyper-interleukin-6 synergizes with murine stroma for long-term survival of primitive human hematopoietic progenitor cells. *Exp Hematol* 29(7): 822–32.
- Lansdorp PM & Dragowska W. (1992) Long-term erythropoiesis from constant numbers of CD34+ cells in serum-free cultures initiated with highly purified progenitor cells from human bone marrow. *J Exp Med* 175(6): 1501–9.
- Miller CL & Eaves CJ. (1997) Expansion in vitro of adult murine hematopoietic stem cells with transplantable lympho-myeloid reconstituting ability. *Proc Natl Acad Sci U S A* 94(25): 13648–53.
- Petzer AL et al. (1996) Self-renewal of primitive human hematopoietic cells (long-term-culture-initiating cells) in vitro and their expansion in defined medium. *Proc Natl Acad Sci U S A* 93(4): 1470–4.
- Petzer AL et al. (1996) Differential cytokine effects on primitive (CD34+CD38-) human hematopoietic cells: novel responses to Flt3-ligand and thrombopoietin. *J Exp Med* 183(6): 2551–8.

Sandrin V et al. (2002) Lentiviral vectors pseudotyped with a modified RD114 envelope glycoprotein show increased stability in sera and augmented transduction of primary lymphocytes and CD34+ cells derived from human and nonhuman primates. *Blood* 100(3): 823–32.

Zandstra PW et al. (1998) Ontogeny-associated changes in the cytokine responses of primitive human haemopoietic cells. *Br J Haematol* 101(4): 770–8.

Zhang CC & Lodish HF. (2005) Murine hematopoietic stem cells change their surface phenotype during ex vivo expansion. *Blood* 105(11): 4314–20.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and StemSpan are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.