RosetteSep™ HLA Lymphoid Cell Enrichment Kit

REF 15271HLA

For 20 tests (10 mL of whole blood per test)

Components:

RosetteSep™ HLA Lymphoid Cell Enrichment Cocktail 5 x 2 mL
RosetteSep™ DM-L Density Medium 2 x 100 mL

ENGLISH

INTENDED USE

RosetteSep™ cell enrichment cocktails are designed for the in vitro enrichment of specific cell subsets from human whole blood.

PRODUCT DESCRIPTION

The RosetteSep™ antibody cocktail crosslinks unwanted cells in human whole blood to multiple red blood cells (RBCs), forming immunorosettes. These dense immunorosettes pellet along with the free RBCs when centrifuged over a density gradient medium such as RosetteSep™ DM-L (Catalog #15705). Desired cells are never labeled with antibody and are easily collected as a highly enriched population at the interface between the plasma and the density gradient medium.

RosetteSep™ HLA Lymphoid Cell Enrichment REF #15271HC.1 Cocktail

This cocktail contains a combination of mouse and rat monoclonal antibodies. These antibodies are bound in bispecific Tetrameric Antibody Complexes (TACs) which are directed against cell surface antigens on unwanted cells (CD16, CD19, CD20, CD36, CD56) and glycophorin A on red blood cells. The mouse monoclonal antibody subclass is IgG₁.

RosetteSep™ DM-L Density Medium REF #15705

RosetteSep™ DM-L is a density separation medium designed specifically for use with RosetteSep™ cocktails for the enrichment of human lymphocytes.

Density: 1.081 g/mL

QUALITY CONTROL

RosetteSep™ HLA Lymphoid Cell Enrichment REF #15271HC.1 Cocktail

RosetteSep™ cell enrichment cocktails are manufactured using aseptic technique and tightly controlled processes.

Each lot of RosetteSep™ cell enrichment cocktail is sterility tested according to USP methods and Quality Control performance tested in cell separation assays using human whole blood.

RosetteSep™ DM-L Density Medium REF #15705

RosetteSep™ DM-L Density Medium is manufactured using aseptic technique and tightly controlled processes.

Each lot of RosetteSep™ DM-L Density Medium is sterility tested according to USP methods.

STORAGE AND STABILITY

RosetteSep™ HLA Lymphoid Cell Enrichment REF #15271HC.1 Cocktail

Store at 2 - 8°C. This product may be shipped at 15 - 25°C, but should be refrigerated upon receipt. Do not freeze. Product stable at 2 - 8°C until expiry date as indicated on label.

RosetteSep™ DM-L Density Medium REF #15705

Store at 15 - 25°C. Storage at 2 - 8°C is acceptable, but ensure that the medium equilibrates to 15 - 25°C and invert bottle to mix contents before use. Keep protected from direct light. Product stable at 15 - 25°C until expiry date as indicated on label.

WARNINGS AND PRECAUTIONS

1. For professional users only.
2. This product is for in vitro diagnostic use.
3. Do not use cocktail or density medium if vial contents have leaked. Unused cocktail or density medium may be disposed of according to standard laboratory procedures for non-hazardous liquids.
4. This product should be handled by trained personnel observing good laboratory practices. Once this product is added to human cells, treat the suspension as potentially biohazardous. Handling of reagents and disposal of wastes should observe all local, state, or national regulations.
5. This product is a potential irritant to eyes, respiratory system, and skin. This product may also be harmful if ingested. Avoid exposure through skin, eye contact, inhalation, and ingestion.

SPECIAL MATERIALS REQUIRED BUT NOT PROVIDED

Recommended Medium

Phosphate buffered saline with 2% fetal bovine serum (PBS + 2% FBS, Catalog #07905).


STEMCELL Technologies Inc | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | www.stemcell.com

For Technical Assistance

Tel: +1.604.877.0713

European Toll-Free Number: 00800 7836 2355

e-mail: techsupport@stemcell.com

MDSS GmbH

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30175 Hannover • Germany

Document #29682

Version 3.2.1

2015

HANDLING AND DIRECTIONS FOR USE

Ensure that blood sample, recommended medium (see Special Materials Required but not Provided), RosetteSep™ DM-L Density Medium and centrifuge are all at room temperature (15 - 25°C).

1. Aliquot 10 mL of whole blood into a 50 mL tube. If desired, retain a small aliquot of blood (500 µL) for flow cytometric analysis of the start sample.
2. Add RosetteSep™ HLA Cocktail at **50 µL/mL** of whole blood (e.g. for 10 mL of whole blood, add 500 µL of cocktail). Mix well.
3. Incubate **20 minutes** at room temperature (15 - 25°C).
4. Dilute sample with an equal volume of PBS + 2% FBS and mix gently.
5. Layer the diluted sample on top of 10 mL of RosetteSep™ DM-L Density Medium (Catalog #15705).

Be careful to minimize mixing of density gradient medium and sample.


6. Centrifuge for **25 minutes** at 330 x g (see Notes) at room temperature (15 - 25°C) with the brake off.
7. Remove the enriched cells from the RosetteSep™ DM-L Density Medium : plasma interface.

Note: It is sometimes difficult to see the cells at the interface, especially when very rare cells are enriched. It is advisable to remove some of the density gradient medium along with the enriched cells in order to ensure optimal recovery.

8. Wash enriched cells with PBS + 2% FBS.
9. Use enriched cells as desired. If you wish to evaluate the cell purity by flow cytometry, we recommend lysing enriched samples with ammonium chloride to remove residual RBCs (this can be done as the wash step).

ROSETTESEP™ PROCEDURE

Numbers refer to steps in Handling and Directions for Use.


STEMCELL Technologies Inc | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | www.stemcell.com

For Technical Assistance

Tel: +1.604.877.0713

European Toll-Free Number: 00800 7836 2355

e-mail: techsupport@stemcell.com

MDSS GmbH

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30175 Hannover • Germany

Document #29682

Version 3.2.1

2015

NOTES

Density Medium

RosetteSep™ DM-L Density Medium has been formulated to optimize lymphoid cell recovery. Using a different density medium may cause cell loss.

Conversion of g to RPM

To convert g to rpm, use the following formula:

$$RPM = \sqrt{\frac{RCF}{(1.118 \times 10^{-5}) \times (\text{Radius})}}$$

Where: RPM = centrifuge speed in revolutions per minute
RCF = relative centrifugal force (g)
Radius = radius of rotor (cm)

Assessing Purity

Purity of lymphoid (CD3⁺) cells can be measured by flow cytometry after staining with a fluorochrome-conjugated anti-CD3 antibody (e.g. Anti-Human CD3, Clone UCHT1, Catalog #60011), or a combination of other lymphoid cell specific antibodies.

To reduce non-specific antibody binding, add normal human serum to all flow cytometry samples (start and enriched) prior to the addition of the antibody stain, at a concentration of 2 µL human serum/100 µL cells.

Typical Results

These results are for illustrative purposes only. They were obtained using samples from normal, healthy adults. Results from individual patient samples may vary.

CATALOG #	CELL TYPE ENRICHED	PURITY
15271HLA	Lymphoid Cells (CD3 ⁺)	> 85%

 REF Catalog or reference number	 LOT Batch code	 Use by: YYYY-MM
 Caution, consult accompanying documents	 IVD In Vitro Diagnostic Medical Device	 For storage within temperature limits
 CE Mark	 Manufacturers identification (name & address)	 EC REP Authorized EC representative in the European Community

TECHNICAL ASSISTANCE

For technical support please contact us by email at techsupport@stemcell.com or call either +1.604.877.0713 or the European Toll-Free number 00800 7836 2355. For more information please visit www.stemcell.com.

If you require a printed copy or a translated version of this document in a certain language please contact technical support.


 MDSS GmbH

Schiffgraben 41
30175 Hannover, Germany

 STEMCELL Technologies Inc | 1618 Station Street, Vancouver, BC | V6A 1B6 | Canada | www.stemcell.com

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