

ArciTect™ tracrRNA Kit

Trans-activating crRNA for guide RNA generation in CRISPR-Cas9 genome editing

Catalog #	76016	5 nmol Kit
	76017	10 nmol Kit
	76018	20 nmol Kit



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

ArciTect™ tracrRNA is a trans-activating CRISPR RNA (crRNA), one of two RNA components required to make a guide RNA (gRNA) for CRISPR-Cas9 genome editing. ArciTect™ tracrRNA contains a complementary sequence to the linker sequence in ArciTect™ crRNA, allowing fast and efficient gRNA duplex formation. ArciTect™ tracrRNA is recognized by all ArciTect™ Cas9 nucleases and is essential for ribonucleoprotein (RNP) complex formation.

The ArciTect™ tracrRNA Kit includes tracrRNA as well as Annealing Buffer (5X), which is required for forming a gRNA duplex between ArciTect™ tracrRNA and ArciTect™ crRNA. ArciTect™ crRNA is available by custom order for user-specified target sequences, in 2 nmol (Catalog #76010), 10 nmol (Catalog #76011), and 20 nmol (Catalog #76012) sizes.

Ordering Information

PRODUCT NAME	CATALOG #	SIZE	COMPONENTS
ArciTect™ tracrRNA Kit	76016	5 nmol Kit	- ArciTect™ tracrRNA (5 nmol), Component #76016A - Annealing Buffer (5X) (100 µL), Component #76019
ArciTect™ tracrRNA Kit	76017	10 nmol Kit	- ArciTect™ tracrRNA (10 nmol), Component #76017A - Annealing Buffer (5X) (100 µL), Component #76019
ArciTect™ tracrRNA Kit	76018	20 nmol Kit	- ArciTect™ tracrRNA (20 nmol), Component #76018A - Annealing Buffer (5X) (100 µL), Component #76019

Storage and Stability

The following components are sold as part of the ArciTect™ tracrRNA Kits (see Ordering Information) and are not available for individual sale.

COMPONENT NAME	COMPONENT #	STORAGE	SHELF LIFE
ArciTect™ tracrRNA	76016A 76017A 76018A	Store at -80°C. Alternatively, store at -20°C for up to 6 months.	Stable for 24 months from date of manufacture (MFG) on label.
Annealing Buffer (5X)	76019	Store at -20°C. Alternatively, store at 2 - 8°C for up to 6 months.	Stable for 24 months from date of manufacture (MFG) on label.

Materials Required But Not Included

- Nuclease-free water (e.g. Sigma Catalog # LSKNF0500)

Directions for Use

The following protocol is for preparation of a 200 μ M tracrRNA stock solution.

1. Briefly centrifuge the vial before opening.
2. Add nuclease-free water as outlined in Table 1.

Table 1. Preparation of 200 μ M* tracrRNA Stock Solution

ArciTect™ tracrRNA	VOLUME OF NUCLEASE-FREE WATER (μ L)
5 nmol (Component #76016A)	25
10 nmol (Component #76017A)	50
20 nmol (Component #76018A)	100

*200 μ M is equal to 200 pmol/ μ L

3. Mix thoroughly. If not used immediately, aliquot and store at -80°C for up to 1 month. After thawing the aliquots, use immediately. Do not re-freeze.

For complete instructions on CRISPR-Cas9 genome editing, refer to the Technical Bulletin: Human Pluripotent Stem Cell Genome Editing Using the ArciTect™ CRISPR-Cas9 System (Document #27084), available at www.stemcell.com or contact us to request a copy.

Related Products

For related products, including other genome editing tools, specialized cell culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com or contact us at techsupport@stemcell.com.

References

- Gundry MC et al. (2016) Highly efficient genome editing of murine and human hematopoietic progenitor cells by CRISPR/Cas9. *Cell Rep* 17(5): 1453–61.
- Hultquist JF et al. (2016) A Cas9 ribonucleoprotein platform for functional genetic studies of HIV-host interactions in primary human T cells. *Cell Rep* 17(5): 1438–52.
- Kim S et al. (2014) Highly efficient RNA-guided genome editing in human cells via delivery of purified Cas9 ribonucleoproteins. *Genome Res* 24(6): 1012–9.
- Liang X et al. (2015) Rapid and highly efficient mammalian cell engineering via Cas9 protein transfection. *J Biotechnol* 208: 44–53.
- Ran FA et al. (2013) Double nicking by RNA-guided CRISPR Cas9 for enhanced genome editing specificity. *Cell* 154(6): 1380–9.
- Rupp LJ et al. (2017) CRISPR/Cas9-mediated PD-1 disruption enhances anti-tumor efficacy of human chimeric antigen receptor T cells. *Sci Rep* 7(1): 737.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ArciTect are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. The purchase of the ArciTect™ products conveys to the purchaser the limited, non-transferable right to use, in accordance with all applicable laws and regulations, the ArciTect™ products and any related material solely for research purposes only, not for any clinical, human, agricultural, veterinary, livestock, or commercial purpose, with no warranty (express or implied) from the owners or assignees of the Patents or STEMCELL Technologies each of whom expressly disclaim any warranty regarding results obtained through the use of the ArciTect™ products, and without any exception the owners or assignees of the Patents or STEMCELL Technologies, or as applicable a director, trustee, officer, employee, agent, faculty, official investigator or student thereof, will not suffer or be exposed to any liability, damages, loss, or expense of any kind whatsoever arising from or related to the use of the ArciTect™ products by a purchaser of same. Distribution of ArciTect™ products by STEMCELL Technologies is covered under at least US 8,697,359, US 8,771,945, US 8,795,965, US 8,865,406, US 8,871,445, US 8,889,356, US 8,889,418, US 8,895,308, US 8,906,616 and foreign equivalents ("Patents"). While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.