

Dissociation Reagents

Collagenase Type IV

For digestion of native collagen fibrils



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 07426
07427

100 mg
1 g

Product Description

Collagenase is a protease consisting of a single polypeptide chain approximately 1,000 amino acid residues in length. Collagenase is capable of digesting native collagen fibrils commonly found in connective tissues and therefore is frequently used for tissue dissociation. Collagenase preparations contain the activity of several proteases, including collagenase, caseinase, clostripain, and trypsin (Kessler & Yaron). Collagenase Type IV contains low levels of tryptic activity and is frequently used for the dissociation of pancreatic islets (Taguchi et al.) and for applications where maintenance of receptor integrity is required.

Product Information

Alternative Names:	Clostridium histolyticum collagenase; Collagenase 4; Collagenase Type 4; Collagenase IV
Format:	Lyophilized powder
Storage:	Store at 2 - 8°C.
Stability:	Stable until expiry date (EXP) on label.
Reconstitution:	Dissociation reagents can be reconstituted in a balanced salt solution or buffer of choice.
Molecular Weight:	68 - 130 kDa
CAS Number:	9001-12-1
Optimum pH:	6.3 - 8.5
Cleavage Site:	-Pro-X- ↑ -Gly-Pro-Y- : X = neutral Y = nonspecific

Specifications

Source:	Clostridium histolyticum
Activity:	Collagenase: ≥ 160 CDU/mg dry weight (mgdw); Caseinase: ≥ 100 u/mgdw; Clostripain: ≤ 3.0 u/mgdw; Trypsin: ≤ 0.1 u/mgdw. See Notes for further information.

Related Products

For a complete list of dissociation reagents, as well as related products available from STEMCELL Technologies, please visit our website at www.stemcell.com or contact us at techsupport@stemcell.com.

Notes

ACTIVITY UNITS

Collagenase: 1 collagenase digestion unit (CDU) equals 1 μ mol of L-leucine equivalents released from collagen in 5 hours at 37°C, pH 7.5.

Caseinase: 1 unit equals 1 μ mol of L-leucine equivalents released from 25 mg vitamin-free casein in 5 hours at 37°C, pH 7.5. Measures non-specific proteolytic activity.

Clostripain: 1 unit hydrolyzes 1 μ mol of N α -benzoyl-L-arginine ethyl ester (BAEE)/minute at 25°C at pH 7.6, after activation in 2.5 mM dithiothreitol (DTT).

Trypsin: 1 unit hydrolyzes 1 μ mol of BAEE/minute at 25°C at pH 7.6.

References

- Baruch K et al. (2015) Breaking immune tolerance by targeting Foxp3(+) regulatory T cells mitigates Alzheimer's disease pathology. *Nat Commun* 6: 7967.
- Choudhery MS et al. (2014) Effect of mild heat stress on the proliferative and differentiative ability of human mesenchymal stromal cells. *Cytotherapy* 17: 359–68.
- Dietrich J et al. (2005) EIF2B5 mutations compromise GFAP+ astrocyte generation in vanishing white matter leukodystrophy. *Nat Med* 11(3): 277–83.
- Fahlgren A et al. (2003) Increased expression of antimicrobial peptides and lysozyme in colonic epithelial cells of patients with ulcerative colitis. *Clin Exp Immunol* 131(1): 90–101.
- Imoukhuede PI & Popel AS. (2012) Expression of VEGF Receptors on Endothelial Cells in Mouse Skeletal Muscle. *PLoS One* 7(9).
- Kessler E & Yaron A. (1973) A novel aminopeptidase from clostridium histolyticum. *Biochem Biophys Res Commun* 50(2): 405–12.
- Presnell SC et al. (2011) Isolation, characterization, and expansion methods for defined primary renal cell populations from rodent, canine, and human normal and diseased kidneys. *Tissue Eng C, Methods* 17(3): 261–73.
- Renn CN et al. (2006) TLR Activation of Langerhans Cell-Like Dendritic Cells Triggers an Antiviral Immune Response. *J Immunol* 177(1): 298–305.
- Taguchi Y et al. (2010) Impaired insulin secretion from the pancreatic islets of hypothyroidal growth-retarded mice. *J Endocrinol* 206(2): 195–204.
- Vatakis DN et al. (2012) Using the BLT humanized mouse as a stem cell based gene therapy tumor model. *J Vis Exp* (70): e4181.
- Wang J et al. (2015) Interleukin 18 function in atherosclerosis is mediated by the interleukin 18 receptor and the Na-Cl co-transporter. *Nat Med* 21(7): 820–6.
- Welte Y et al. (2013) Patient derived cell culture and isolation of CD133+ putative cancer stem cells from melanoma. *J Vis Exp* (73): e50200.
- Xue X & Shah YM. (2013) In vitro organoid culture of primary mouse colon tumors. *J Vis Exp* (75): e50210.
- You Y & Brody SL. (2013) Culture and differentiation of mouse tracheal epithelial cells. *Methods Mol Biol* 945: 123–43.
- Zhou YY et al. (2000) Culture and adenoviral infection of adult mouse cardiac myocytes: methods for cellular genetic physiology. *Am J Physiol Hear Circ Physiol* 279(1): H429–36.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.