

ArciTect™ T7 Endonuclease I Kit

For estimation of CRISPR-Cas9 genome editing efficiency

Catalog # 76021 1 Kit 25 Reactions
76022 1 Kit 125 Reactions



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

ArciTect™ T7 Endonuclease I is the preferred enzyme for detecting genome editing such as insertions or deletions (INDELs) generated by CRISPR-Cas9. ArciTect™ T7 Endonuclease I Kit is comprised of ArciTect™ T7 Endonuclease I and ArciTect™ T7 Endonuclease I Buffer (10X), which have been tested and validated for use with the ArciTect™ CRISPR-Cas9 genome editing system. ArciTect™ T7 Endonuclease I recognizes and cleaves mismatched DNA, cruciform DNA structures, Holliday structures or junctions, heteroduplex DNA, and, less efficiently, nicked double-stranded DNA. Since the cleavage efficiency is proportional to the number of INDELs created at a specific DNA target, ArciTect™ T7 Endonuclease I Kit is used to estimate gene-editing efficiency in a rapid and cost effective manner.

Ordering Information

PRODUCT NAME	CATALOG #	SIZE	COMPONENTS
ArciTect™ T7 Endonuclease I Kit	76021	1 Kit - 25 Reactions	- ArciTect™ T7 Endonuclease I (25 Reactions) - ArciTect™ T7 Endonuclease I Buffer (1 mL)*
ArciTect™ T7 Endonuclease I Kit	76022	1 Kit - 125 Reactions	- ArciTect™ T7 Endonuclease I (125 Reactions) - ArciTect™ T7 Endonuclease I Buffer (1 mL)*

*1 mL of buffer is sufficient for up to 500 reactions.

Component Storage and Stability

The following components are sold as part of ArciTect™ T7 Endonuclease I Kits (see Ordering Information) and are not available for individual sale.

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE
ArciTect™ T7 Endonuclease I	76023	25 Reactions	Store at -20°C.	Stable for 24 months from date of manufacture (MFG) on label.
ArciTect™ T7 Endonuclease I	76024	125 Reactions	Store at -20°C.	Stable for 24 months from date of manufacture (MFG) on label.
ArciTect™ T7 Endonuclease I Buffer (10X)	76025	1 mL	Store at -20°C.	Stable for 24 months from date of manufacture (MFG) on label.

Specifications

Source: E. coli
Concentration: 10 Units (U)/μL
Activity: Each Unit converts > 90% of 1 μg of supercoiled cruciform pUC (AT) plasmid DNA into linear DNA in 1 hour at 37°C.
Purity: > 95%
Cleavage Site: Single strand cleavage at the first, second, or third phosphodiester bond 5' to a DNA mismatch

Materials Required But Not Included

PRODUCT NAME	CATALOG #
Genomic DNA isolation kit	e.g. Norgen Biotek 24700
PCR tubes	e.g. Corning PCR-02-C
Forward and reverse primers*	---
Nuclease-free water	e.g. Sigma LSKNF0500
High-fidelity DNA polymerase master mix	e.g. New England BioLabs M0494
Thermocycler	---
PCR purification kit	e.g. QIAGEN 28104
Microvolume spectrophotometer	---
Proteinase K	e.g. New England BioLabs P8107S
DNA loading buffer	e.g. New England BioLabs B7024S
DNA ladder	e.g. New England BioLabs N0550S
Agarose gel apparatus and reagents	---

*Forward and reverse primers should be designed to amplify a ~1000 base pair region surrounding the target site, offset so that T7 digestion will result in two fragments of distinctly different sizes for resolution on an agarose gel (e.g. 700 and 300 base pair fragments). For further information, refer to the Technical Bulletin: Design Considerations for the ArciTect™ CRISPR-Cas9 Genome Editing System (Document #27083), available at www.stemcell.com or contact us to request a copy.

Directions for Use

- Isolate genomic DNA (gDNA) from edited cells.
- Prepare Reagent Mix for PCR amplification of target region from 100 ng of gDNA as indicated in Table 1.
NOTE: Indicated reaction volumes are for Q5® Hot Start High-Fidelity 2X Master Mix (New England BioLabs Catalog #M0494). For other DNA polymerases, adjust component concentrations as required.

Table 1. Reagent Mix for PCR Amplification of Target Region

COMPONENT	FINAL CONCENTRATION OR AMOUNT	VOLUME PER REACTION
DNA polymerase master mix	1X	25 µL
10 µM Forward primer	0.5 µM	2.5 µL
10 µM Reverse primer	0.5 µM	2.5 µL
Template DNA	100 ng	Concentration-dependent
Nuclease-free water	---	As needed
Total volume	---	50 µL

- Amplify the target region by PCR, using the conditions indicated in Table 2.

Table 2. PCR Cycling Conditions for Amplification of Target Region

STEP	TEMPERATURE	TIME
Initial denaturation	98°C	30 seconds
35 cycles	98°C	10 seconds
	50 - 72°C* (annealing temperature)	15 seconds
	72°C	30 seconds
Final extension	72°C	2 minutes
Hold	4°C	Infinite

*The annealing temperature should be ~5°C below the lowest melting temperature of the primers. When using ArciTect™ Human HPRT Positive Control Kit (Catalog #76013), the recommended annealing temperature is 67°C.

4. Extract PCR product using a PCR purification kit, then measure the concentration using a microvolume spectrophotometer.
5. Transfer 200 ng of PCR product to a new PCR tube, and adjust to a final volume of 17 µL with nuclease-free water. Add 2 µL of ArciTect™ T7 Endonuclease I Buffer (10X).

NOTE: If an uncut control is desired, prepare an identical second tube. This control will demonstrate that PCR generates a band of the expected size, which is only cleaved in the presence of T7 Endonuclease I.

6. Denature and anneal the PCR product in a thermocycler using the conditions indicated in Table 3.

NOTE: If a thermocycler with the indicated ramp rates is unavailable, heat to 95°C for 10 minutes and let samples cool to room temperature (15 - 25°C).

Table 3. Denaturing and Annealing Conditions for PCR Products

STEP	TEMPERATURE	RAMP RATE	TIME
Initial denaturation	95°C	N/A	5 minutes
Annealing	95 - 85°C	-2°C/second	---
	85 - 25°C	-0.1°C/second	---
Hold	4°C	N/A	Infinite

N/A: Not applicable.

7. Add 1 µL of ArciTect™ T7 Endonuclease I. Incubate at 37°C for 15 minutes.

NOTE: Do not add ArciTect™ T7 Endonuclease I to the uncut control tube.

NOTE: Avoid an incubation temperature above 42°C, as this may cause an increase in non-specific nuclease activity.

8. Terminate the ArciTect™ T7 Endonuclease I reaction by adding 1 µL of proteinase K and incubate at 37°C for 5 minutes.

NOTE: If not used immediately, store at -20°C. Thaw at room temperature (15 - 25°C) prior to use.

9. Perform gel electrophoresis by loading each reaction product with DNA loading buffer into individual wells of a 1% agarose gel.

10. Using a gel imager, quantify cut (edited) and uncut (non-edited) band intensities for the digested condition(s). Using these values, estimate cleavage efficiency for each edited sample as follows:

Cleavage efficiency (%) = (cut band intensity ÷ sum of cut and uncut band intensities) x 100

NOTE: Band intensity, due to dye intercalation, is proportional not only to the number of DNA copies, but also to DNA length. As a result, the smaller cleavage product may be difficult to quantify and can be excluded from this estimate. Moreover, since ArciTect™ T7 Endonuclease I does not recognize single base pair mismatches or homozygous mutations, the actual efficiency of genome editing may be higher than estimated.

See Figure 1 for representative examples of INDEL detection by T7 Endonuclease I assay on agarose gels.




Figure 1. INDEL Detection by T7 Endonuclease I Assay. Human embryonic stem (ES) and induced pluripotent stem (iPS) cells (A) and T cells (B) were edited and INDEL formation was assessed using ArciTect™ Human HPRT Positive Control Kit (Catalog #76013) and ArciTect™ T7 Endonuclease I Kit. Following CRISPR-mediated editing at the HPRT locus, genomic DNA was isolated and a 1 kb region surrounding the target site was amplified by PCR using ArciTect™ CRISPR Human HPRT Primer Mix (included with Catalog #76013). PCR products were purified, then denatured, re-annealed, and cut with ArciTect™ T7 Endonuclease I. Samples were resolved on a 1% agarose gel alongside Quick-Load® Purple 2-Log DNA Ladder (0.1 - 10.0 kb; New England BioLabs Catalog #N0550S), and band intensities were determined using a ChemiDoc™ MP Imaging System (Bio-Rad). Relative intensities of the 1 kb and 0.8 kb bands were used to calculate the INDEL efficiency (%). Control: Uncut PCR product (no T7 added); Test, Donor 1, and Donor 2: T7-digested PCR product.

Related Products

For related products, including other genome editing tools, specialized cell culture and storage media, supplements, antibodies, cytokines, and small molecules, visit www.stemcell.com or contact us at techsupport@stemcell.com.

References

- Mean RJ et al. (2004) Modification of the enzyme mismatch cleavage method using T7 endonuclease I and silver staining. *Biotechniques* 36(5): 758–60.
- Vouillot L et al. (2015) Comparison of T7E1 and surveyor mismatch cleavage assays to detect mutations triggered by engineered nucleases. *G3 (Bethesda)* 5(3): 407–15.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ArciTect are trademarks of STEMCELL Technologies Canada Inc. Q5 and Quick-Load are registered trademarks of New England Biolabs, Inc. ChemiDoc is a trademark of Bio-Rad Laboratories, Inc. All other trademarks are the property of their respective holders. The purchase of the ArciTect™ products conveys to the purchaser the limited, non-transferable right to use, in accordance with applicable laws and regulations, such purchased products and any related material solely for research purposes only, not for any commercial purposes, with no warranty from the owner or assignee of the Patents, and with no liability to the owner or assignee of the Patents, or a director, trustee, officer, employee, agent, faculty, official investigator or student thereof. Distribution of ArciTect™ products is covered under at least US 8,697,359, US 8,771,945, US 8,795,965, US 8,865,406, US 8,871,445, US 8,889,356, US 8,889,418, US 8,895,308, US 8,906,616 and foreign equivalents ("Patents"). While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.