

mTeSR™1 Without Phenol Red

Defined, feeder-free maintenance medium for human ES and iPS cells

Catalog #05876

1 Kit



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

mTeSR™1 medium is a complete, serum-free, defined formulation designed for the feeder-free maintenance and expansion of human embryonic stem (ES) cells^{1,2} and human induced pluripotent stem (iPS) cells³⁻⁵ in the undifferentiated state. mTeSR™1 Without Phenol Red is for applications in which the absence of phenol red may be desirable, such as live fluorescence imaging of cell cultures. Complete mTeSR™1 Without Phenol Red medium (Basal Medium + 5X Supplement) contains recombinant human basic fibroblast growth factor (rh bFGF) and recombinant human transforming growth factor β (rh TGF β). Addition of further growth factors is not required.

mTeSR™1 Without Phenol Red may be used with either Corning® Matrigel® hESC-Qualified Matrix (Corning Catalog #354277) or Vitronectin XF™ (Catalog #07180, a matrix developed and manufactured by Primorigen Biosciences) as the culture matrix.

Each lot of mTeSR™1 5X Supplement is used to prepare complete mTeSR™1 medium and then performance tested in a culture assay using human pluripotent stem cells.

Component Storage and Stability

The following components are sold as part of the mTeSR™1 Without Phenol Red kit (Catalog #05876) and are not available for individual sale.

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE
mTeSR™1 Without Phenol Red Basal Medium	05877	400 mL	Store at 2 - 8°C. Do not freeze.	Stable for 12 months from date of manufacture (MFG) on label.
mTeSR™1 5X Supplement	05852	100 mL	Store at -20°C.	Stable until expiry date (EXP) on label.

Preparation of Complete mTeSR™1 Without Phenol Red Medium

Use sterile techniques to prepare complete mTeSR™1 Without Phenol Red medium (Basal Medium + 5X Supplement). The following example is for preparing 500 mL of complete medium. If preparing other volumes, adjust accordingly.

NOTE: Thaw supplement or complete medium at room temperature (15 - 25°C) or overnight at 2 - 8°C. Do not thaw in a 37°C water bath.

1. Thaw mTeSR™1 5X Supplement and mix thoroughly.

NOTE: Once thawed, use supplement immediately or aliquot and store at -20°C for up to 3 months. Do not exceed the shelf life of the supplement. After thawing the aliquoted supplement, use immediately. Do not re-freeze.

2. Add 100 mL of mTeSR™1 5X Supplement to 400 mL of mTeSR™1 Without Phenol Red Basal Medium. Mix thoroughly.

NOTE: If not used immediately, store complete mTeSR™1 Without Phenol Red medium at 2 - 8°C for up to 2 weeks. Alternatively, aliquot and store at -20°C for up to 6 months. Do not exceed the shelf life of the individual components. After thawing the aliquoted complete medium, use immediately or store at 2 - 8°C for up to 2 weeks. Do not re-freeze.

If prepared aseptically, complete mTeSR™1 Without Phenol Red medium is ready for use. If desired, the medium can be filtered using a 0.2 - 0.22 μ m low protein binding polyethersulfone (PES) filter unit (e.g. Fisher 09-741-04 [0.2 μ m, 250 mL]; Fisher SCGP00525 [0.22 μ m, 50 mL]).

Directions for Use

For complete instructions on how to maintain human ES and iPS cells in mTeSR™1, refer to the Technical Manual: Maintenance of Human Pluripotent Stem Cells in mTeSR™1 (Document #28315) available at www.stemcell.com or contact us to request a copy.

Assessment of hPSCs

The following antibodies can be used to characterize hPSCs by flow cytometry or immunocytochemistry:

- Anti-Human SSEA-4 Antibody, Clone MC-813-70 (Catalog #60062)
- Anti-Human TRA-1-60 Antibody, Clone TRA-1-60R (Catalog #60064)
- Anti-Human OCT4 (OCT3) Antibody, Clone 3A2A20 (Catalog #60093)

For complete flow cytometry protocols and antibodies that can be used, refer to the Technical Manual: Maintenance of Human Pluripotent Stem Cells in mTeSR™1 (Document #28315), available at www.stemcell.com or contact us to request a copy.

Related Products

For related products, including specialized cell culture and storage media, matrices, antibodies, cytokines, and small molecules, visit www.stemcell.com/hPSCworkflow or contact us at techsupport@stemcell.com.

References

1. Ludwig TE et al. (2006) Derivation of human embryonic stem cells in defined conditions. *Nat Biotechnol* 24(2): 185–7.
2. Ludwig TE et al. (2006) Feeder-independent culture of human embryonic stem cells. *Nat Methods* 3(8): 637–46.
3. Yu J et al. (2007) Induced pluripotent stem cell lines derived from human somatic cells. *Science* 318(5858): 1917–20.
4. Masaki H et al. (2007) Heterogeneity of pluripotent marker gene expression in colonies generated in human iPS cell induction culture. *Stem Cell Res* 1(2): 105–15.
5. Sun N et al. (2009) Feeder-free derivation of induced pluripotent stem cells from adult human adipose stem cells. *Proc Natl Acad Sci USA* 106(37): 15720–5.



This product was developed under license to intellectual property owned by WiCell™ Research Institute.

This product is sold for research use only (whether the buyer is an academic or for-profit entity) under a non-transferable, limited-use license. Purchase of this product does not include the right to sell, use or otherwise transfer this product for commercial purposes (i.e., any activity undertaken for consideration, such as use of this product for manufacturing, or resale of this product or any materials made using this product, or use of this product or any materials made using this product to provide services) or clinical use (i.e., administration of this product or any material using this product to humans) or the right to implant any material made using this product into an animal by, or in collaboration with, a for-profit entity, for purposes other than basic pre-clinical research applications (including without limitation teratoma assays) to validate the function of the cells. Purchasers wishing to use the product for purposes other than research use should contact Asterias Biotherapeutics, Inc. legal department at (650) 433-2900 or legal@asteriasbio.com. Purchasers who do not agree to the terms and conditions set forth above should return the product in acceptable conditions to the seller for a refund.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. mTeSR is a trademark of WARF. Corning and Matrigel are registered trademarks of Corning Incorporated. Vitronectin XF is developed and manufactured by Primorigen Biosciences, Inc. and Vitronectin XF is a trademark of Primorigen Biosciences, Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.