

Alpha MEM without Nucleosides

Minimum Essential Medium Eagle - alpha modification (Alpha MEM) without nucleosides

Catalog #36453

500 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Minimum Essential Medium Eagle - Alpha Modification (Alpha MEM) without Nucleosides is recommended for a wide variety of cell culture applications. Selection of a suitable nutrient medium is dependent on type of cell, conditions of culture, and degree of chemical definition. MEM has typically been used for the cultivation of cells grown in monolayers; however, additional supplementation to Alpha MEM without Nucleosides has expanded the applications of this medium.

Alpha MEM without Nucleosides contains L-glutamine and is supplemented with sodium bicarbonate.

See Formulation on page 2 for a full list of components.

Properties

Storage: Store at 2 - 8°C.

Shelf Life: Stable for 12 months from date of manufacture (MFG) on label.

Directions for Use

This product contains labile amino acid L-glutamine, which has a half-life of approximately 12 months when stored at 2 - 8°C. For certain cell culture applications, medium stored for greater than 12 months following date of manufacture should be supplemented with additional L-glutamine. Add 5 mL of 200 mM L-Glutamine (Catalog #07100) to 500 mL of Alpha MEM without Nucleosides to achieve a final concentration of 2 mM.

NOTE: This product does not contain antibiotics. If desired, add antibiotics and use within 1 week.

Formulation

INORGANIC SALTS	g/L
Calcium Chloride•2H ₂ O	0.2
Magnesium Sulfate (anhydrous)	0.09767
Potassium Chloride	0.4
Sodium Chloride	6.8
Sodium Phosphate Dibasic	0.14
AMINO ACIDS	g/L
L-Alanine	0.025
L-Arginine	0.105
L-Asparagine•H ₂ O	0.05
L-Aspartic Acid	0.03
L-Cysteine HCl•H ₂ O	0.1
L-Cystine•2HCl	0.024
L-Glutamic Acid	0.075
L-Glutamine	0.292
Glycine	0.05
L-Histidine	0.031
L-Isoleucine	0.0524
L-Leucine	0.0524
L-Lysine	0.058
L-Methionine	0.015
L-Phenylalanine	0.032
L-Proline	0.04
L-Serine	0.025
L-Threonine	0.048
L-Tryptophan	0.01
L-Tyrosine•2Na•2H ₂ O	0.052
L-Valine	0.046
VITAMINS	g/L
Ascorbic Acid	0.05
Biotin	0.0001
Choline Chloride	0.001
Folic Acid	0.001
myo-Inositol	0.002
Niacinamide	0.001
D-Calcium Pantothenate	0.001
Pyridoxal•HCl	0.001
Riboflavin	0.0001
Thiamine•HCl	0.001
Vitamin B12	0.00136
OTHER	g/L
D-Glucose	1
Phenol Red	0.01
Sodium Bicarbonate	2.2
Sodium Pyruvate	0.11
Thioctic Acid	0.0002

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.