

LY294002

PI3K/AKT pathway inhibitor;
Inhibits PI3K


Catalog #72152 5 mg

Catalog #72154 25 mg

Product Description

LY294002 is a PI3K inhibitor that has greater potency and selectivity than quercetin, the structure on which it is based. LY294002 inhibits PI3K ($IC_{50} = 1.4 \mu M$) and also shows activity against CK2, but not PI4K, EGFR, PDGFR, MAPK, PKA, or PKC. (Davies et al.; Vlahos et al.)

| | |
|---------------------|--|
| Alternative Names: | Not applicable |
| CAS Number (Model): | 154447-36-6 |
| Chemical Formula: | $C_{19}H_{17}NO_3$ |
| Molecular Weight: | 307.3 g/mol |
| Purity: | $\geq 98\%$ |
| Chemical Name: | 2-(4-morpholinyl)-8-phenyl-4H-1-benzopyran-4-one |
| Structure: | |



Properties

| | |
|-------------------------------|---|
| Product Format: | A crystalline solid |
| Stability and Storage: | Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt. |
| Preparation: | <p>Solubility:</p> <ul style="list-style-type: none">· DMSO \leq 50 mM· Absolute ethanol \leq 50 mM· PBS (pH 7.2) \leq 160 μM <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 5 mg in 1.63 mL of fresh DMSO.</p> <p>Prepare stock solution fresh before use. Stock solutions in DMSO or ethanol are stable for up to 6 months if stored at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p> |

Published Applications

MAINTENANCE AND SELF-RENEWAL

- Suppresses proliferation and self-renewal of mouse embryonic stem (ES) cells (Lianguzova et al.; Paling et al.).

DIFFERENTIATION

- Promotes differentiation to insulin-producing cells from mouse ES cells (Hori et al.).
- Inhibits myotube formation from myoblasts (Coolican et al.; Jiang et al.).

References

- Coolican SA et al. (1997) The mitogenic and myogenic actions of insulin-like growth factors utilize distinct signaling pathways. *J Biol Chem* 272 (10): 6653-62.
- Davies SP et al. (2000) Specificity and mechanism of action of some commonly used protein kinase inhibitors. *Biochem J* 351(Pt 1): 95-105.
- Hori Y et al. (2002) Growth inhibitors promote differentiation of insulin-producing tissue from embryonic stem cells. *Proc Natl Acad Sci USA* 99 (25): 16105-10.
- Jiang BH et al. (1998) An essential role of phosphatidylinositol 3-kinase in myogenic differentiation. *Proc Natl Acad Sci USA* 95(24): 14179-83.
- Lianguzova MS et al. (2007) Phosphoinositide 3-kinase inhibitor LY294002 but not serum withdrawal suppresses proliferation of murine embryonic stem cells. *Cell Biol Int* 31(4): 330-7.
- Paling NRD et al. (2004) Regulation of embryonic stem cell self-renewal by phosphoinositide 3-kinase-dependent signaling. *J Biol Chem* 279(46): 48063-70.
- Vlahos CJ et al. (1994) A specific inhibitor of phosphatidylinositol 3-kinase, 2-(4-morpholinyl)-8-phenyl-4H-1-benzopyran-4-one (LY294002). *J Biol Chem* 269(7): 5241-8.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

Warning

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.