

Human Recombinant TGF-beta 2

Transforming growth factor beta 2

Catalog #78174 10 µg

Catalog #78174.1 50 µg

Product Description

Transforming growth factor (TGF) beta 2 is a member of the TGF- β superfamily and regulates diverse cellular phenotypes. Similar to TGF- β 1 & - β 3, TGF- β 2 signals via serine-threonine kinase type I and II receptors and activates signal transduction via SMAD family proteins, regulating a variety of functions, such as cell proliferation, differentiation, wound healing, apoptosis, and metabolism (de Caestecker; Massague; Zuniga et al.). TGF- β 2 is important in many developmental processes, for example, mice with TGF- β 2 deletions show defects in the development of cardiac, lung, craniofacial, limb, eye, ear, and urogenital systems (Dunker & Kreigstein).

Product Information

Alternative Names:	BSC-1 cell growth inhibitor, Cetermin, Glioblastoma-derived T-cell suppressor factor, G-TSF, Polyergin, TGFB2
Accession Number:	P61812
Amino Acid Sequence:	ALDAAYCFRN VQDNCCLRPL YIDFKRDLGW KWIHEPKGYN ANFCAGACPY LWSSDTQHSR VLSLYNTINP EASASPCCVS QDLEPLTILY YIGKTPKIEQ LSNMIVKSCK CS
Predicted Molecular Mass:	12.7 kDa
Species:	Human
Product Formulation:	Lyophilized from a filtered solution containing hydrochloric acid.
Source:	Human Cells
Purity:	≥ 95%

Specifications

Activity:	The specific activity is $\geq 5.5 \times 10^6$ units/mg ($EC_{50} \leq 0.18$ ng/mL) as determined by the ability to inhibit the IL-4-dependent proliferation of TF-1 human erythroleukemic cells.
Endotoxin Level:	Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is ≤ 1 EU/ μ g protein.

Preparation and Storage

Stability and Storage:	Store at -80°C. Stable as supplied for 12 months from date of receipt.
Preparation:	Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex.

OPTIONAL: After reconstitution, if product will not be used immediately, dilute with concentrated bovine serum albumin (BSA) to a final BSA concentration of 0.1%. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at 2 - 8°C for more than 1 week or at -20 to -80°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Data




Figure 1. Biological Activity and Molecular Mass of Human Recombinant TGF-beta 2

(A) The biological activity of Human Recombinant TGF-beta 2 was tested by its ability to inhibit the IL-4-dependent proliferation of TF-1 human erythroleukemic cells. The EC₅₀ is defined as the effective concentration of the growth factor at which inhibition is at 50% of maximum. The EC₅₀ in the example above is less than 0.18 ng/mL.

(B) 5 µg of Human Recombinant TGF-beta 2 was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant TGF-beta 2 has a predicted molecular mass of 12.7 kDa.

Related Products

For a complete list of cytokines or peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com/cytokines, or contact us at techsupport@stemcell.com.

References

- de Caestecker M. (2004) The transforming growth factor-beta superfamily of receptors. *Cytokine Growth Factor Rev* 15(1): 1-11.
- Dünker N & Kriegstein K. (2000) Targeted mutations of transforming growth factor-beta genes reveal important roles in mouse development and adult homeostasis. *Eur J Biochem* 267(24): 6982-8.
- Massagué J. (2000) How cells read TGF-beta signals. *Nat Rev Mol Cell Biol* 1(3): 169-78.
- Zúñiga JE et al. (2005) Assembly of betaRI:TbetaRII:TGFbeta ternary complex in vitro with receptor extracellular domains is cooperative and isoform-dependent. *J Mol Biol* 354(5): 1052-68.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2025 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.