

# Fetal Bovine Serum

## Fetal bovine serum for cell culture and research applications

Catalog #100-0179 (NZ Origin) 500 mL  
Catalog #100-0180 (AUS Origin) 500 mL



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

STEMCELL fetal bovine serum (FBS) can be used in a broad range of applications, including as a serum supplement for culture media or as a research assay reagent. STEMCELL FBS is sterile and has low endotoxin and hemoglobin levels.

## Properties

**Storage:** Store at -20°C.  
**Shelf Life:** Stable until expiry date (EXP) on label.  
**Physical Appearance:** Clear amber liquid

## Directions For Use

Thaw FBS at 2 - 8°C prior to use. Use as directed in the protocol of choice.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright ©2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.