

Dyes and Stains

Rhod-2 AM (Bromide)

A cell-permeable calcium dye that exhibits fluorescence upon increasing Ca²⁺ concentrations

Catalog #100-0992

1 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Rhod-2 AM (Bromide) is a cell-permeable acetoxymethyl (AM) ester dye which fluoresces in the presence of Ca²⁺. Cellular uptake is facilitated by the AM group and is removed by cytoplasmic esterases, resulting in intracellular accumulation of a cell-impermeable Ca²⁺ indicator, Rhod-2 (Smith et al.). Increasing Ca²⁺ concentrations significantly enhance Rhod-2 AM fluorescence, which makes the long wavelength of this dye valuable for applications involving cells with high levels of autofluorescence. Increased intracellular calcium concentrations are a critical indicator leading to the activation of apoptotic pathways (Mukherjee et al.). Use Rhod-2 AM (Bromide) dye to measure mitochondrial calcium influx by detection using flow cytometry, fluorescence microscopy, and microplate readers.

CAS Number:	145037-81-6
Chemical Formula:	C ₅₂ H ₅₉ N ₄ O ₁₉ • Br
Molecular Weight:	1124.0 g/mol
Excitation Wavelength:	553 nm
Emission Wavelength:	577 nm
Cutoff:	570 nm
Structure:	



Properties

Storage:	Store at -20°C.
Shelf Life:	Stable until expiry date (EXP) on label. Protect product from prolonged exposure to light.
Format:	Red solid

Directions for Use

Please read the entire protocol before proceeding. The following protocol is for staining adherent cells in a black-wall/clear-bottom plate.

Preparation of Rhod-2 AM (Bromide) Stock and Working Solution

1. To prepare a stock solution, dissolve Rhod-2 AM (Bromide) in dimethyl sulfoxide (DMSO) at 2 - 5 mM.
NOTE: If not used immediately, aliquot and store at -20°C. Do not exceed the expiry date as indicated on the label. Avoid repeated freeze-thaw cycles.
2. To prepare the Rhod-2 AM (Bromide) working solution, dilute the stock solution (prepared in step 2) to 2 - 20 µM with Hanks' Balanced Salt Solution with 20 mM HEPES (HHBS) or a buffer of choice. Use the working solution immediately; do not store.
NOTE: Although 4 - 5 µM is the recommended final concentration for most cell lines, titrate the dye for optimal performance.
3. (OPTIONAL): Additional reagents may be added to the working solution to improve performance as follows:
 - To increase the solubility of Rhod-2 AM (Bromide), add 0.04% Pluronic® F-127 (AAT Bioquest Catalog #20053).
 - To reduce leakage of de-esterified indicators in cells that contain organic anion transporters, add RadiUse™ probenecid (AAT Bioquest) at a final concentration of 0.5 - 1.0 mM.

Staining Samples

1. Plate cells in a black-wall/clear-bottom plate with culture medium and allow cells to adhere overnight.
2. Add the Rhod-2 AM (Bromide) working solution directly to the culture medium at the desired concentration.
NOTE: If serum in the culture medium interferes with your stimulant of choice, replace the culture medium with HHBS before adding the working solution.
3. Incubate at 37°C and 5% CO₂ for at least 30 minutes.
NOTE: The optimal incubation time should be determined for different cell lines.
4. Remove the Rhod-2 AM (Bromide) working solution and wash the cells with HHBS or a buffer of choice.
NOTE: For cells that contain organic anion transporters, use buffer containing 1 mM RadiUse™ probenecid.

Fluorescence Detection

Add a stimulant of choice and monitor the fluorescence intensity as follows:

- Fluorescence microscope: TRITC or appropriate filter set
- Fluorescence microplate reader containing a programmable liquid-handling system, such as a functional drug-screening system (FDSS), FLIPR®, or FlexStation®: Ex/Em = 540/590 nm (cutoff = 570 nm)

References

1. Mukherjee SB et al. (2002) Increase in cytosolic Ca²⁺ levels through the activation of non-selective cation channels induced by oxidative stress causes mitochondrial depolarization leading to apoptosis-like death in *Leishmania donovani* promastigotes. *J Biol Chem* 277(27): 24717–27.
2. Smith NA et al. (2018) Fluorescent Ca²⁺ indicators directly inhibit the Na,K-ATPase and disrupt cellular functions. *Sci Signal* 11(515): eaal2039.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2022 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Pluronic is a registered trademark of BASF Corporation. RadiUse is a trademark of AAT Bioquest Inc. FLIPR and FlexStation are registered trademarks of Molecular Devices, LLC. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.