

Costar® 12-Well Flat-Bottom Plate, Tissue Culture-Treated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

**Sterile, clear polystyrene, flat-bottom, tissue culture-treated
multiwell plate with lid**

Catalog #200-0624

50 Plates

Product Description

Costar®12-Well Flat-Bottom Plate, Tissue Culture-Treated is ideal for attachment-dependent cell and tissue culture in a multiwell plate format. These polystyrene plates are transparent and easy to stack and feature non-reversible lids with condensation rings to reduce cross-contamination. Each well is labeled alphanumerically for easy identification.

Catalog #200-0624, 50 Plates; individually wrapped (Corning Catalog #3513)

APPLICATIONS

- Cell culture
- Microscopy

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Stable until expiry date on manufacturer label.

Features

MATERIALS

- Plate and lid: clear and transparent polystyrene
- Tissue culture-treated for optimal cell attachment
- Hydrophilic and negatively charged
- Non-pyrogenic

DIMENSIONS

- Conforms to standard plate footprint and dimensions

CAPACITY

- Approximately 6.9 mL per well
- Recommended working volume of 0.76 - 1.14 mL per well
- Approximate cell growth area of 3.8 cm² per well

STERILITY

- Irradiation

Directions for Use

Use as directed in protocol of choice.

Notes and Tips

For chemical resistance or thermoplastic property information, contact us at techsupport@stemcell.com.

Related Products

For a complete list of related products, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Costar is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness.