

## Small Molecules

### BMN 673

Poly ADP ribose polymerase (PARP) inhibitor

Catalog #100-1129

10 mg



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

BMN 673 is a poly ADP ribose polymerase (PARP) inhibitor (PARP1 IC<sub>50</sub> = 0.57 nM; Shen et al.). BMN 673 is an inhibitor of PARP in cancer cells with BRCA1, BRCA2, or PTEN mutations (Shen et al.). BMN 673 binds to PARP and prevents PARP-mediated repair of single-strand DNA breaks, leading to DNA damage (Shen et al.).

Alternative Names:	Talazoparib
CAS Number:	1207456-01-6
Chemical Formula:	C <sub>19</sub> H <sub>14</sub> F <sub>2</sub> N <sub>6</sub> O
Molecular Weight:	380.4 g/mol
Purity:	≥ 98%
Chemical Name:	(8S,9R)-5-fluoro-8-(4-fluorophenyl)-2,7,8,9-tetrahydro-9-(1-methyl-1H-1,2,4-triazol-5-yl)-3H-pyrido [4,3,2-de]phthalazin-3-one

Structure:



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	<ul style="list-style-type: none"><li>• DMSO ≤ 50 mM</li><li>• Absolute ethanol ≤ 655 μM</li></ul> For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 2.63 mL of DMSO. Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application. Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.

## Published Applications

### CANCER RESEARCH

- Traps PARP1 at the DNA damage site and induces cytotoxic effects in breast cancer cells (Sethy & Kundu).
- Increases the number of CD8+ T cells and NK cells, as well as production of IFN- $\gamma$  and TNF- $\alpha$ , to inhibit growth of BRCA1-deficient murine epithelial ovarian cancer cell line BR5FVB1-Akt (Huang et al.).
- Inhibits BRCA2-associated pancreatic ductal adenocarcinoma (PDAC) growth in a PDAC murine xenograft model (Andrei et al.).

## References

Andrei AZ et al. (2015) Increased in vitro and in vivo sensitivity of BRCA2-associated pancreatic cancer to the poly (ADP-ribose) polymerase-1/2 inhibitor BMN 673. *Cancer Lett* 364(1): 8–16.

Huang J et al. (2015) The PARP1 inhibitor BMN 673 exhibits immunoregulatory effects in a BRCA1  $-/-$  murine model of ovarian cancer. *Biochem Biophys Res Commun* 463(4): 551–6.

Sethy C & Kundu CN. (2022) PARP inhibitor BMN-673 induced apoptosis by trapping PARP-1 and inhibiting base excision repair via modulation of pol- $\beta$  in chromatin of breast cancer cells. *Toxicol Appl Pharmacol* 436: 115860.

Shen Y et al. (2013) BMN 673, a novel and highly potent PARP1/2 inhibitor for the treatment of human cancers with DNA repair deficiency. *Clin Cancer Res* 19(18): 5003–15.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

**This product is hazardous. Please refer to the Safety Data Sheet (SDS).**

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.