

# EBV (EBNA-3A) Peptide Pool

## Epstein-Barr virus (EBNA-3A) peptide pool for immune cell activation

Catalog #100-1386

~25 µg (15 nmol)/peptide



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

EBV (EBNA-3A) Peptide Pool is a lyophilized mixture of 234 peptides from Epstein-Barr nuclear antigen 3 (EBNA-3A) of Epstein-Barr virus (EBV; strain B95-8). The pool consists of 15-mer peptides with 11-amino-acid overlaps that cover amino acids 1 - 944 on EBNA-3A. EBNA-3A contributes to B cell growth transformation (Hertle et al.; Tomkinson et al.). It has been found to interact with host cell proteins, where it activates (Young et al.) or represses (Cludts & Farrell; Hickabottom et al.) transcription of specific cellular genes. One unit of this product (i.e. ~25 µg/peptide) is sufficient for stimulating  $2.5 \times 10^8$  cells.

### APPLICATIONS

- Antigen-specific T cell stimulation
- Cellular immune response
- Immune monitoring
- T cell assays
- T cell expansion

## Product Information

Number of Peptides: 234

Source: Epstein-Barr virus (strain B95-8) (also known as human herpesvirus 4 [HHV-4])

Accession Number: P12977

Protein Name: Epstein-Barr nuclear antigen 3A (EBNA-3A)

Protein Sequence: MDKDRPGPPALDDNMEEEVPSTSVVQEQVSAGDWENVLIELSDSSSEKEAEDAHLPAQKGTKRKRVDHDAGG  
SAPARPMPLPPQPDLPGREAILRRFPLDLRTLLQAIGAAATRIDTRAIDQFFGSQISNTEMYIMYAMAIRQAIRDRRRN  
PASRRDQAKWRLQTLAAGWPMGYQAYSSWMYSYTDHQTTPTFVHLQATLGCTGGRRCHVTFSGATFKLPRCTP  
GDRQWLYVQSSVGNIVQSCNPRYSIFFDYMAIHRSLTKIWEEVLTPDQRVSFMEFLGFLQRTDLSYIKSFVSDALGT  
TSIQTPWIDNPNSTETAQAWNAGFLRGRAYGIDLLRTEGEHVEGATGETREESEDTESDGDDDELPCIVSRGGPKV  
KRPIFIRRLHRLLLMRAGKRTEQGKEVLEKARGSTYGTTPRPVVKPRPEVPQSDETATSHGSAQVPEPPTIHLAAQ  
GMAYPLHEQHGMAPCPVAQAPPTPLPPVSPGDQLPGVFSQDGRVACAPVPAPAGPIVRPWEPSLTQAAGQAFAPV  
RPQHMPVEVPVPVPTVALERPVPYKPVPRAPPKIAMQGPGETSGIRRARERWRPAPWTPNPPRSPSQMSVRDLAR  
LRAEAQVKQASVEVQPPQLTQVSPQQPMQPLVPEQQMFPGAPFSQVADVVRAPGVAMQPPQYFDLPLIQPISQ  
GAPVAPLRASMGPPVPVPATQPPYFDIPLTEPINQGASAAHFLPQQPMQPLVPEQQWMFPGAALSQSVRPGVAQS  
QYFDLPLTQPINHGAPAAHFLHQPMEGPWWPEQWMFQGAPPSQGTQVQHLQDALGYTLHGLNHPGVVPVSPA  
VNQYHLSQAAFGLPIDEDESGEGSDTSEPCEALDLSIHGRPCQAPWVPVQEEGGQDATEVLDSLHGRPRPRT  
EWPVQGEQQNVGTGP ETRRVVSAVVHMCQDDEFDPDLQDPPDEA

Gene Name: EBNA3

Purity: Average 70%

Formulation: Lyophilized as trifluoroacetate salts

## Preparation and Storage

Storage: Store at -20°C.

Stability: Stable as supplied until expiry date (EXP) on label.

Preparation: Warm to room temperature (15 - 25°C) before reconstitution. Add pure dimethyl sulfoxide (DMSO; ~40 µL) and dilute with water to the desired concentration. Final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system. If not used immediately, aliquot and store at -20°C. Protect from light. Avoid repeated freeze-thaw cycles.

## Related Products

For a complete list of peptide pools, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Cludts I & Farrell PJ. (1998) Multiple functions within the Epstein-Barr virus EBNA-3A protein. *J Virol* 72(3): 1862–9.
- Hertle ML et al. (2009) Differential gene expression patterns of EBV infected EBNA-3A positive and negative human B lymphocytes. *PLoS Pathog* 5(7): e1000506.
- Hickabottom M et al. (2002) Two nonconsensus sites in the Epstein-Barr virus oncoprotein EBNA3A cooperate to bind the co-repressor carboxyl-terminal-binding protein (CtBP). *J Biol Chem* 277(49): 47197–204.
- Tomkinson B et al. (1993) Epstein-Barr virus nuclear proteins EBNA-3A and EBNA-3C are essential for B-lymphocyte growth transformation. *J Virol* 67(4): 2014–25.
- Young P et al. (2008) Epstein-Barr virus nuclear antigen (EBNA) 3A induces the expression of and interacts with a subset of chaperones and co-chaperones. *J Gen Virol* 89(4):866–77.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.