

HPV16 (L1) Peptide Pool

Human papillomavirus 16 (L1) peptide pool for immune cell activation

Catalog #100-1392

~25 µg (15 nmol)/peptide



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

HPV16 (L1) Peptide Pool is a lyophilized mixture of 124 peptides from major capsid protein L1 of human papillomavirus type 16 (HPV16). The pool consists of 15-mer peptides with 11-amino-acid overlaps that cover amino acids 1 - 505 on L1. L1 enters the nuclei of host cells via a Kap $\alpha\beta$ 1-mediated pathway (Nelson et al.). It cooperates with E2, another viral protein, to enhance viral DNA replication (Siddiqi et al.). L1 pentamers can assemble into virus-like particles (Chen et al., 2000; Chen et al., 2001). One unit of this product (i.e. 25 µg/peptide) is sufficient for stimulating 2.5×10^8 cells.

APPLICATIONS

- Antigen-specific T cell stimulation
- Immune monitoring
- T cell assays
- T cell expansion

Product Information

Number of Peptides: 124

Source: Human papillomavirus type 16

Accession Number: Q9WLQ6

Protein Name: Major capsid protein L1

Protein Sequence: MSLWLPSEATVYLPVPVSKVVSTDEYVARTNIYYHAGTSRLLAVGHPYFPIKKPNNNKILVPKVSGLQYRVFRIHLP
DPNKFGFPDTSFYNPDTQRLWACVGVEVGRGQPLGVGISGHPLLNLDDTENASAYAANAGVDNRECISMDYK
QTQLCLIGCKPPIGEHWGKGSPTNVAVNPGDCPPLELINTVIQDGMVDTGFGAMDFTTLQANKSEVPLDICTSI
CKYPDYIKMVSEPYGDSLFFYLRRQMFVRHLFNAGAVGENVPDDLYIKGSGSTANLASSNYFPTPSGSMVTSDA
QIFNKPYWLQRAQGHNNIGICWGNQLFVTVDTRSTNMSLCAAISTSETTYKNTNFKEYLRHGEEYDLQFIFQLCKI
TLTADVMTYIHSMNSTILEDWNFGLQPPPGGTLEDYRFVTSQAIACQKHTPPAPKEDPLKKYTFWEVNLKEKFSAD
LDQFPLGRKFLLQAGLKAKPKFTLGKRRKATPTTSSTSTAKRKKRKL

Gene Name: L1

Purity: Average 70%

Formulation: Lyophilized as trifluoroacetate salts

Preparation and Storage

Storage: Store at -20°C.

Stability: Stable as supplied until expiry date (EXP) on label.

Preparation: Warm to room temperature (15 - 25°C) before reconstitution. Add pure dimethyl sulfoxide (DMSO; ~40 µL) and dilute with water to the desired concentration. Final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system. If not used immediately, aliquot and store at -20°C. Protect from light. Avoid repeated freeze-thaw cycles.

Related Products

For a complete list of peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

References

- Chen XS et al. (2000) Structure of small virus-like particles assembled from the L1 protein of human papillomavirus 16. *Mol Cell* 5(3): 557–67.
- Chen XS et al. (2001) Papillomavirus capsid protein expression in *Escherichia coli*: purification and assembly of HPV11 and HPV16 L1. *J Mol Biol* 307(1): 173–82.
- Nelson LM et al. (2002) Nuclear import strategies of high risk HPV16 L1 major capsid protein. *JBC* 277(26): 23958–64.
- Siddiq A et al. (2015) The human papillomavirus type 16 L1 protein directly interacts with E2 and enhances E2-dependent replication and transcription activation. *J Gen Virol* 96(Pt 8): 2274–85.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.