

HHV8 (K8.1) Peptide Pool

Human herpesvirus 8 (K8.1) peptide pool for immune cell activation

Catalog #100-1404

~25 µg (15 nmol)/peptide



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

HHV8 (K8.1) Peptide Pool is a lyophilized mixture of 55 peptides from glycoprotein K8.1 of human herpesvirus 8 (HHV8). The pool consists of 15-mer peptides with 11-amino-acid overlaps that cover amino acids 1 - 228 on K8.1. K8.1 is a viral envelope protein that is expressed in the late-lytic phase (Zoetewij et al.). K8.1 binds to cell surface heparin sulfate and facilitates HHV8 infection (Birkmann et al.; Wang et al.). One unit of this product (i.e. 25 µg/peptide) is sufficient for stimulating 2.5×10^8 cells.

APPLICATIONS

- Antigen-specific T cell stimulation
- Cellular immune response
- Immune monitoring
- T cell assays
- T cell expansion

Product Information

Number of Peptides:	55
Source:	Human herpesvirus 8 (also known as Kaposi's sarcoma-associated herpesvirus [KSHV])
Accession Number:	O36551
Protein Name:	Glycoprotein K8.1
Protein Sequence:	MSSTQIRTEIPVALLILCLCLVACHANCPTYRSHLGFQEGWSGQVYQDWLGRMNCSYENMTALEAVSLNGTRLA AGSPSSEYPNVSVSVEDTSASGSGEDAIDESGSGEEERPVTSHVTFMTQSVQATTELTDALISAFSGSYSSGEPST TRIRVSPAENGRNSGASNRVPFSATTTTTRGRDAHYNAEIRTHLYILWAVGLLLGLVLILYLCVPRCRRKKPYIV
Gene Name:	K8.1
Purity:	Average 70%
Formulation:	Lyophilized as trifluoroacetate salts

Preparation and Storage

Storage:	Store at -20°C.
Stability:	Stable as supplied until expiry date (EXP) on label.
Preparation:	Warm to room temperature (15 - 25°C) before reconstitution. Add pure dimethyl sulfoxide (DMSO; ~40 µL) and dilute with water to the desired concentration. Final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system. If not used immediately, aliquot and store at -20°C. Protect from light. Avoid repeated freeze-thaw cycles.

Related Products

For a complete list of peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com, or contact us at techsupport@stemcell.com.

References

Birkmann A et al. (2001) Cell surface heparan sulfate is a receptor for human herpesvirus 8 and interacts with envelope glycoprotein K8.1. J Virol 75(23): 11583–93.

Wang FZ et al. (2001) Human herpesvirus 8 envelope glycoprotein K8.1A interaction with the target cells involves heparan sulfate. J Virol 75(16): 7517–27.

Zoetewij JP et al. (1999) Identification and rapid quantification of early-and late-lytic human herpesvirus 8 infection in single cells by flow cytometric analysis: characterization of antiherpesvirus agents. J Virol 73(7): 5894–902.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.