

# Anti-Human IL-2 Antibody, Clone MT8G10, FITC

Mouse monoclonal antibody against human, cynomolgus IL-2, FITC-conjugated

Catalog #100-1464

100 Tests

5 µL/test

## Product Description

This mouse monoclonal antibody (clone MT8G10) reacts with human interleukin 2 (IL-2). IL-2 is an ~16 kDa cytokine composed of a four-alpha-helix bundle produced primarily by CD4+ T helper (Th) cells and to a lesser degree by CD8+ T cells, natural killer (NK) cells, and NK T cells. It binds to the IL-2 receptor  $\alpha$  (IL-2R $\alpha$ ; CD25), which complexes with IL-2R $\beta$  (CD122) and  $\gamma$ c (CD132) subunits and leads to the activation of JAK1/3, PI3K-AKT, and MAPK signaling pathways. IL-2 mediates T cell growth and proliferation, facilitates the differentiation of CD4+ T cells to Th1 and Th2 cells, promotes CD8+ T cell and NK cell cytolytic activity, and drives the development of CD4+FOXP3+ regulatory T cells. The MT8G10 antibody is suitable for the detection of intracellular IL-2 by flow cytometry.

<b>Target Antigen:</b>	IL-2
<b>Alternative Names:</b>	IL2, interleukin 2, interleukin-2, lymphokine, T cell growth factor, TCGF
<b>Gene ID:</b>	3558
<b>Species Reactivity:</b>	Human, Cynomolgus
<b>Host Species:</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Clone:</b>	MT8G10
<b>Isotype:</b>	IgG1
<b>Immunogen:</b>	Recombinant human IL-2
<b>Conjugate:</b>	FITC (Fluorescein isothiocyanate)

## Applications


Verified Applications: FC

Abbreviations: CellSep: Cell separation; CHIP: Chromatin immunoprecipitation; FA: Functional assay; FACS: Fluorescence-activated cell sorting; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IHC-P: Immunohistochemistry (paraffin-embedded); IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

## Properties

- Purification:** The antibody was purified by affinity chromatography and conjugated with FITC.
- Stability and Storage:** Product stable at 2 - 8°C when stored undiluted. Do not freeze. Protect product from prolonged exposure to light. Stable until expiry date (EXP) on label.
- Directions for Use:** For flow cytometry, the suggested use of this antibody is 5 µL per 1 x 10<sup>6</sup> cells in 50 µL. It is recommended that the antibody be titrated for optimal performance for each application.

## Data



**Figure 1. Flow Cytometry Analysis of Unstimulated and Stimulated Human T Cells Labeled With Anti-Human IL-2 Antibody, Clone MT8G10 FITC**

(A) Flow cytometry analysis of unstimulated human T cells labeled with anti-human IL-2 Antibody, Clone MT8G10, FITC and anti-human CD8a antibody, clone RPA-T8, APC (Catalog #60022AZ.1). (B) Flow cytometry analysis of PMA/ionomycin-stimulated human T cells labeled with anti-human IL-2 Antibody, Clone MT8G10, FITC and anti-human CD8a antibody, clone RPA-T8, APC (Catalog #60022AZ.1).

## Related Products

For a complete list of antibodies, including other conjugates, sizes, and clones, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/antibodies](http://www.stemcell.com/antibodies), or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

Arenas-Ramirez N et al. (2015) Interleukin-2: biology, design and application. *Trends Immunol* 36(12): 763–77.

Boyman O & Sprent J. (2012) The role of interleukin-2 during homeostasis and activation of the immune system. *Nat Rev Immunol* 12(3): 180–90.

Liao W et al. (2013) Interleukin-2 at the crossroads of effector responses, tolerance, and immunotherapy. *Immunity* 38(1): 13–25.

Malek TR & Castro I. (2010) Interleukin-2 receptor signaling: at the interface between tolerance and immunity. *Immunity* 33(2): 153–65.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.