

# Lestaurtinib

Epigenetic modifier; Inhibits FMS-like tyrosine kinase 3 (FLT3)

Catalog #100-1558

5 mg

## Product Description

Lestaurtinib is an FMS-like tyrosine kinase 3 (FLT3) inhibitor ( $IC_{50} = 2$  nM). FLT3 mutations cause increased activation via autophosphorylation, promoting cell growth, and inhibiting apoptosis. Lestaurtinib inhibits FLT3 phosphorylation and its downstream targets in acute myeloid leukemia (AML; Levis et al.). Other tyrosine kinases such as Janus kinase 2 (JAK2), which are mutated in myeloproliferative disorders, can also be inhibited by Lestaurtinib ( $IC_{50} = 1$  nM) in erythroid cells (Hexner et al.).

**Alternative Names:** CEP-701, KT-5555

**CAS Number:** 111358-88-4


**Chemical Formula:**  $C_{26}H_{21}N_3O_4$

**Molecular Weight:** 439.5 g/mol

**Purity:**  $\geq 98\%$

**Chemical Name:** Not applicable

**Structure:**



## Properties

<b>Product Format:</b>	A white to off-white powder
<b>Stability and Storage:</b>	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
<b>Preparation:</b>	<ul style="list-style-type: none"><li>• DMSO ≤ 225 mM</li><li>• Absolute ethanol ≤ 45 mM</li></ul> For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 228 µL of DMSO.  Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.  Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.

## Published Applications

### CANCER RESEARCH

- Induces a cytotoxic response in human acute myeloid leukemia cell lines and prolongs the survival in a leukemia mouse model (Levis et al.).
- Reduces tumor cell proliferation, increases DNA damage, and induces apoptosis in a medulloblastoma mouse model (Pallavicini et al.).
- Induces apoptosis in human glioma cell lines and suppresses tumor growth in glioma mouse xenograft models (Cao et al.).
- Reduces viability of human B cell lymphoma (BCL) cells and suppresses tumor growth in a BCL mouse xenograft model (Beck et al.).

## References

- Beck D et al. (2016) Synthetic lethal screen demonstrates that a JAK2 inhibitor suppresses a BCL6-dependent IL10RA/JAK2/STAT3 pathway in high grade b-cell lymphoma. *J Biol Chem* 291(32): 16686–98.
- Cao Y et al. (2020) Lestaurtinib potentiates TRAIL-induced apoptosis in glioma via CHOP-dependent DR5 induction. *J Cell Mol Med* 24(14): 7829–40.
- Hexner EO et al. (2008) Lestaurtinib (CEP701) is a JAK2 inhibitor that suppresses JAK2/STAT5 signaling and the proliferation of primary erythroid cells from patients with myeloproliferative disorders. *Blood* 111(12): 5663–71.
- Levis M et al. (2002) A FLT3-targeted tyrosine kinase inhibitor is cytotoxic to leukemia cells in vitro and in vivo. *Blood* 99(11): 3885–91.
- Pallavicini G et al. (2023) Lestaurtinib inhibits Citron kinase activity and medulloblastoma growth through induction of DNA damage, apoptosis and cytokinesis failure. *Front Oncol* 13: 1202585.
- Wu M et al. (2018) FLT3 inhibitors in acute myeloid leukemia. *J Hematol Oncol* 11(1): 133.

## Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.