

3-Deazaneplanocin A (Hydrochloride)

Epigenetic modifier; inhibits S-adenosyl-L-homocysteine hydrolase (SAHase)

Catalog #100-1640

1 mg

Product Description

3-Deazaneplanocin A (Hydrochloride) inhibits S-adenosyl-L-homocysteine hydrolase (SAHase; $IC_{50} = 0.08 - 0.24 \mu M$), which results in reduction of the histone methyltransferase enhancer of zeste homolog 2 (EZH2) protein levels (Tseng et al.). EZH2 is part of the polycomb repressive complex 2 (PRC2), an important regulator of tumor suppressor genes. EZH2 is associated with regulating cell proliferation, apoptosis, senescence, and DNA damage repair. It can be selectively targeted by 3-Deazaneplanocin A (Hydrochloride), making it an attractive target for cancer therapy (Duan et al.; Tan et al.). This product is supplied as the hydrochloride salt of the molecule.

Alternative Names: 2,3-DMCC, DZnep, NSC 617989

CAS Number: 120964-45-6


Chemical Formula: $C_{12}H_{14}N_4O_3 \cdot HCl$

Molecular Weight: 298.7 g/mol

Purity: $\geq 98\%$

Chemical Name: 5R-(4-amino-1H-imidazo[4,5-c]pyridin-1-yl)-3-(hydroxymethyl)-3-cyclopentene-1S,2R-diol, monohydrochloride

Structure:



Properties

Product Format:	A crystalline solid
Stability and Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Preparation:	<ul style="list-style-type: none">• Phosphate-buffered saline (PBS; pH 7.2) ≤ 16 mM• DMSO ≤ 10 mM For example, to prepare a 10 mM stock solution in PBS (pH 7.2), resuspend 1 mg in 335 µL of PBS (pH 7.2). Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

DIFFERENTIATION

- Supports differentiation of reprogrammed urine-derived cells into myotubes by promoting expression of myogenin (Takizawa et al.).
- Contributes to erythroid differentiation of human leukemia cells by promoting expression of erythroid-related genes (Fujiwara et al.).

CANCER RESEARCH

- Induces cell death and reduces cell migration in human chondrosarcoma cell lines (Girard et al.).
- Increases apoptosis in human gastric cancer cells with high expression of EZH2 (Pan et al.).

References

- Duan R et al. (2020) EZH2: A novel target for cancer treatment. *J Hematol Oncol* 13(1): 104.
- Fujiwara T et al. (2014) 3-Deazaneplanocin A (DZNep), an inhibitor of S-adenosylmethionine-dependent methyltransferase, promotes erythroid differentiation. *J Biol Chem* 289(12): 8121–34.
- Girard N et al. (2014) 3-Deazaneplanocin A (DZNep), an inhibitor of the histone methyltransferase EZH2, induces apoptosis and reduces cell migration in chondrosarcoma cells. *PLoS One* 9(5): e98176.
- Pan YM et al. (2016) STAT3 signaling drives EZH2 transcriptional activation and mediates poor prognosis in gastric cancer. *Mol Cancer* 15(1): 79.
- Takizawa H et al. (2019) Modelling Duchenne muscular dystrophy in MYOD1-converted urine-derived cells treated with 3-deazaneplanocin A hydrochloride. *Sci Rep* 9(1): 3807.
- Tan J et al. (2007) Pharmacologic disruption of polycomb-repressive complex 2-mediated gene repression selectively induces apoptosis in cancer cells. *Genes Dev* 21(9): 1050–63.
- Tseng CKH et al. (1989) Synthesis of 3-deazaneplanocin A, a powerful inhibitor of S-adenosylhomocysteine hydrolase with potent and selective in vitro and in vivo antiviral activities. *J Med Chem* 32(7): 1442–6.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.