

GSK-J4 (Hydrochloride)

Histone modifier; Inhibits Jumonji domain-containing protein D3 (JMJD3)

Catalog #100-1642

25 mg

Product Description

GSK-J4 (Hydrochloride) is an inhibitor of the Jumonji domain-containing protein D3 (JMJD3) histone demethylase ($IC_{50} = 9 \mu M$). JMJD3 demethylates the tri-methylated lysine 27 on histone H3 (H3K27me3) to activate transcription (Kruidenier et al.). Dysregulation of H3K27me3 marks have been associated with different cancers. By inhibiting JMJD3, GSK-J4 can impact key pathways such as cell proliferation, cell cycle, and apoptosis that are often targeted by anti-cancer agents (Dalpatraj et al.).

Alternative Names: Not applicable

CAS Number: 1797983-09-5


Chemical Formula: $C_{24}H_{27}N_5O_2 \cdot HCl$

Molecular Weight: 454 g/mol

Purity: $\geq 98\%$

Chemical Name: N-[2-(2-Pyridinyl)-6-(1,2,4,5-tetrahydro-3H-3-benzazepin-3-yl)-4-pyrimidinyl]- β -alanine, ethyl ester, hydrochloride

Structure:



Properties

Product Format:	A yellow powder
Stability and Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Preparation:	<ul style="list-style-type: none">• Water ≤ 20 mM• DMSO ≤ 75 mM <p>For example, to prepare a 10 mM stock solution in water, resuspend 10 mg in 2.2 mL of water. If not fully dissolved, warm the 10 mM stock solution in a 37°C water bath or incubator with periodic mixing until the solution is clear.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.</p>

Published Applications

CANCER RESEARCH

- Inhibits cell proliferation and migration and induces apoptosis in human glioma cell lines (Sui et al.).
- Induces apoptosis and cell cycle arrest in human acute myeloid leukemia cells and slows disease progression in a leukemia mouse model (Chu et al.; Li et al.).
- Reduces cell proliferation of castration-resistant human prostate cancer cells (Morozov et al.).
- Induces cell cycle arrest and inhibits cell proliferation of human and mouse medulloblastoma cell lines (Deng et al.).
- Increases differentiation, endoplasmic reticulum stress, and apoptosis in human neuroblastoma cell lines and reduces tumor growth in a mouse xenograft model (Lochmann et al.).

References

- Chu X et al. (2020) GSK-J4 induces cell cycle arrest and apoptosis via ER stress and the synergism between GSK-J4 and decitabine in acute myeloid leukemia KG-1a cells. *Cancer Cell Int* 20(1): 209.
- Dalpatraj N et al. (2023) GSK-J4: An H3K27 histone demethylase inhibitor, as a potential anti-cancer agent. *Int J Cancer* 153(6): 1130–8.
- Deng H et al. (2022) Targeting H3K27me3 demethylase to inhibit Shh signaling and cholesterol metabolism in medulloblastoma growth. *Front Oncol* 12: 1057147.
- Kruidenier L et al. (2012) A selective jumonji H3K27 demethylase inhibitor modulates the proinflammatory macrophage response. *Nature* 488 (7411): 404–8.
- Li Y et al. (2018) Therapeutic potential of GSK-J4, a histone demethylase KDM6B/JMJD3 inhibitor, for acute myeloid leukemia. *J Cancer Res Clin Oncol* 144(6): 1065–77.
- Lochmann TL et al. (2018) Targeted inhibition of histone H3K27 demethylation is effective in high-risk neuroblastoma. *Sci Transl Med* 10(441): eaoo4680.
- Morozov VM et al. (2017) Inhibitor of H3K27 demethylase JMJD3/UTX GSK-J4 is a potential therapeutic option for castration resistant prostate cancer. *Oncotarget* 8(37): 62131–42.
- Sui A et al. (2017) The pharmacological role of histone demethylase JMJD3 inhibitor GSK-J4 on glioma cells. *Oncotarget* 8(40): 68591–8.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

GSK-J4 (Hydrochloride)

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.