

# A-485

p300/CREB-binding protein (CBP) inhibitor

Catalog #100-1645

5 mg

## Product Description

A-485 is an inhibitor of p300 and CREB-binding protein (CBP;  $IC_{50} = 60$  nM). p300 and CBP are histone acetyltransferases (HATs) that regulate transcription through acetylating lysine on the histone tail or by acting as a transcriptional co-activator. A-485 selectively binds to the catalytic active site of p300 and competes with acetyl-coenzyme A (Acetyl-CoA) to inhibit p300 activity (Lasko et al.). The inhibitory effect of A-485 on histone acetylation is most pronounced on histone 3 lysine 18 (H3K18) and lysine 27 (H3K27) (Weinert et al.).

**Alternative Names:** Not applicable

**CAS Number:** 1889279-16-6


**Chemical Formula:**  $C_{25}H_{24}F_4N_4O_5$

**Molecular Weight:** 536.5 g/mol

**Purity:**  $\geq 98\%$

**Chemical Name:** (1R)-N-[(4-fluorophenyl)methyl]-2,3-dihydro-5-[[[(methylamino)carbonyl]amino]-2',4'-dioxo-N-[(1S)-2,2,2-trifluoro-1-methylethyl]-spiro[1H-indene-1,5'-oxazolidine]-3'-acetamide

**Structure:**



## Properties

<b>Product Format:</b>	A crystalline solid
<b>Stability and Storage:</b>	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
<b>Preparation:</b>	<ul style="list-style-type: none"> <li>• DMSO <math>\leq</math> 25 mM</li> <li>• Absolute ethanol <math>\leq</math> 35 mM</li> </ul> <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 186 <math>\mu</math>L of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p>

## Published Applications

### CANCER RESEARCH

- Reduces cell growth, growth hormone secretion and expression of key oncogenes in mouse pituitary adenoma cells in vitro and in a mouse xenograft model (Ji et al.).
- Inhibits human breast cancer cell growth by disrupting estrogen receptor function (Waddell et al.).
- Induces senescence through reactive oxygen species production to activate autophagy pathway, resulting in cell growth arrest of human lung cancer cells (Ansari et al.).
- Induces apoptosis and differentiation of mouse leukemia stem cells and delays leukemia progression in a mouse acute myeloid leukemia model (Pan et al.).

## References

- Ansari MSZ et al. (2023) Pharmacological targeting of CBP/p300 drives a redox/autophagy axis leading to senescence-induced growth arrest in non-small cell lung cancer cells. *Cancer Gene Ther* 30(1): 124–36.
- Ji C et al. (2022) The p300 inhibitor A-485 exerts antitumor activity in growth hormone pituitary adenoma. *J Clin Endocrinol Metab* 107(6): e2291–300.
- Lasko LM et al. (2017) Discovery of a selective catalytic p300/CBP inhibitor that targets lineage-specific tumours. *Nature* 550(7674): 128–32.
- Pan F et al. (2023) Enhancer remodeling drives MLL oncogene-dependent transcriptional dysregulation in leukemia stem cells. *Blood Adv* 7(11): 2504–19.
- Waddell A et al. (2021) Pharmacological inhibition of cbp/p300 blocks estrogen receptor alpha (E $\alpha$ ) function through suppressing enhancer h3k27 acetylation in luminal breast cancer. *Cancers (Basel)* 13(11): 2799.
- Weinert BT et al. (2018) Time-resolved analysis reveals rapid dynamics and broad scope of the CBP/p300 acetylome. *Cell* 174(1): 231–44.

## Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.