

Cyclosporin A

Immunosuppressive agent; Inhibits calcineurin

Catalog #100-1646

100 mg

Product Description

Cyclosporin A is an immunosuppressive agent that inhibits calcineurin ($IC_{50} = 7 \text{ nM}$) and has also been shown to inhibit cytochrome P-450 3A (CYP3A) in human liver microsomes ($K_i = 0.98 \text{ }\mu\text{M}$; Amundsen et al.). Cyclosporin A binds to a cytosolic protein called cyclophilin to form a complex that inhibits calcineurin (Handschumacher et al.). Calcineurin is a Ca^{2+} - and calmodulin-dependent serine/threonine protein phosphatase involved in interleukin-2 transcription and T cell activation (Hermann-Kleiter & Baier).

Alternative Names: CsA, Cyclosporine, Sandimmune

CAS Number: 59865-13-3


Chemical Formula: $\text{C}_{62}\text{H}_{111}\text{N}_{11}\text{O}_{12}$

Molecular Weight: 1202.6 g/mol

Purity: $\geq 98\%$

Chemical Name: Not applicable

Structure:



Properties

Product Format:	A white powder
Stability and Storage:	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
Preparation:	<ul style="list-style-type: none">• DMSO ≤ 20 mM• Absolute ethanol ≤ 20 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 832 µL of DMSO. Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application. Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.

Published Applications

DIFFERENTIATION

- Stimulates colony growth of Thy-1.1^{low}Sca-1⁺Lin⁻ mouse hematopoietic stem and progenitor cells (HSPCs) in low concentrations in combination with growth factors but inhibits colony growth at high concentrations (Perry et al.).
- Promotes cell survival of neural stem and progenitor cells in mice through a calcineurin-independent pathway (Hunt et al.; Sachewsky et al.).

IMMUNOLOGY

- Suppresses T cell activation and proliferation through the inhibition of calcineurin/nuclear factor of activated T cells (NFAT) signaling (Ho et al.).
- Interferes with T cell maturation and selection in the thymus and with deletion of self-reactive cells (Damoiseaux et al; Jenkins et al.).
- Prevents rejection of allografts in experimental animal models (Lim & White).

DISEASE MODELING

- Inhibits formation of lesions and development of disease in mouse and rat models of rheumatoid arthritis (Pozo et al.; Takagishi et al.).
- Decreases levels of pro-inflammatory cytokines, reduces lymphocyte infiltration, and improves organ pathology in a mouse model of acute graft versus host disease (Corbett et al.).

References

- Amundsen R et al. (2012) Cyclosporine A- and tacrolimus-mediated inhibition of CYP3A4 and CYP3A5 in vitro. *Drug Metab Dispos* 40(4): 655–61.
- Corbett JM et al. (2021) Drug delivery formulation impacts cyclosporine efficacy in a humanized mouse model of acute graft versus host disease. *Transpl Immunol* 65: 101373.
- Damoiseaux JGMC et al. (1997) Multiple effects of cyclosporin A on the thymus in relation to T cell development and autoimmunity. *Clin Immunol Immunopathol* 82(3): 197–202.
- Handschumacher RE et al. (1984) Cyclophilin: a specific cytosolic binding protein for cyclosporin A. *Science* 226(4674): 544–7.
- Hermann-Kleiter N & Baier G. (2010) NFAT pulls the strings during CD4+ T helper cell effector functions. *Blood* 115(15): 2989–97.
- Ho S et al. (1996) The mechanism of action of cyclosporin A and FK506. *Clin Immunol Immunopathol* 80(3Pt 2): S40–5.
- Hunt J et al. (2010) Cyclosporin A has direct effects on adult neural precursor cells. *J Neurosci* 30(8): 2888.
- Lim SML & White DJG. (1989) Cyclosporin A and tolerance induction in experimental animals. *Cyclosporin* 82–111.
- Jenkins MK et al. (1988) Effects of cyclosporine A on T cell development and clonal deletion. *Science* 241(4873): 1655–8.
- Perry SS et al. (1999) Direct effects of cyclosporin A on proliferation of hematopoietic stem and progenitor cells. *Cell Transplant* 8(4): 339–44.
- Pozo del E et al. (1992) Prevention of adjuvant arthritis by cyclosporine in rats. *Semin Arthritis Rheum* 21(6): 23–9.
- Sachewsky N et al. (2014) Cyclosporin A enhances neural precursor cell survival in mice through a calcineurin-independent pathway. *Dis Model Mech* 7(8): 953–61.
- Takagishi K et al. (1992) Comparative study of effects of cyclosporins A and G on collagen arthritis in mice. *Agents Actions* 37(3-4): 284–9.

Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

Warning

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.