

# EPZ004777 (Formate)

Histone modifier; Inhibits disruptor of telomeric silencing 1-like (DOT1L)

Catalog #100-1655

10 mg

## Product Description

EPZ004777 (Formate) is a potent inhibitor of disruptor of telomeric silencing 1-like (DOT1L) histone lysine methyltransferase ( $IC_{50} = 400$  pM). By inhibiting DOT1L, EPZ004777 (Formate) prevents DOT1L-mediated methylation of lysine 79 of histone 3 (H3K79; Daigle et al.). H3K79 methylation impacts telomeric silencing, transcriptional regulation, DNA-damage repair, and cell cycle progression (Nguyen & Zhang). Dysregulation of DOT1L-mediated H3K79 methylation is involved in the development of mixed-lineage leukemia (MLL) and the pathophysiology of other cancers (McLean et al.).

**Alternative Names:** Not applicable

**CAS Number:** Not applicable


**Chemical Formula:**  $C_{28}H_{41}N_7O_4 \bullet CH_2O_2$

**Molecular Weight:** 585.7 g/mol

**Purity:**  $\geq 98\%$

**Chemical Name:** 7-[5-deoxy-5-[[[3-[[[4-(1,1-dimethylethyl)phenyl]amino]carbonyl]amino]propyl](1-methylethyl)amino]- $\beta$ -D-ribofuranosyl]-7H-pyrrolo[2,3-d]pyrimidin-4-amine, monoformic acid

**Structure:**



## Properties

<b>Product Format:</b>	A crystalline solid
<b>Stability and Storage:</b>	Product stable at -20°C as supplied. As a precaution, STEMCELL recommends storing all small molecules away from direct light. For long-term storage, store with a desiccant. Stable as supplied for 12 months from date of receipt.
<b>Preparation:</b>	<ul style="list-style-type: none"> <li>• DMSO <math>\leq</math> 17 mM</li> <li>• Absolute ethanol <math>\leq</math> 17 mM</li> </ul> <p>For example, to prepare a 10 mM stock solution in DMSO, resuspend 10 mg in 1.71 mL of DMSO.</p> <p>Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported; however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.</p> <p>Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO or absolute ethanol concentration above 0.1% due to potential cell toxicity.</p>

## Published Applications

### CANCER RESEARCH

- Reduces cell viability and colony formation in human prostate cancer cell lines and patient-derived xenograft model (Vatapalli et al.).
- Inhibits cell growth and upregulates pro-apoptotic genes in a mouse xenograft model of ovarian cancer and human ovarian cancer cell lines (Chava et al.).
- Reduces mouse pancreatic and colon cancer cell viability and colony formation and reduces tumor growth in a mouse xenograft model (Zhou et al.).
- Induces apoptosis, cell-cycle arrest, and terminal differentiation in human acute myeloid leukemia cells and suppresses leukemic tumor growth in a mouse xenograft model (Rau et al.).

## References

- Chava S et al. (2021) Disruptor of telomeric silencing 1-like promotes ovarian cancer tumor growth by stimulating pro-tumorigenic metabolic pathways and blocking apoptosis. *Oncogenesis* 10(7): 48.
- Daigle SR et al. (2011) Selective killing of mixed lineage leukemia cells by a potent small-molecule DOT1L inhibitor. *Cancer Cell* 20(1): 53–65.
- McLean CM et al. (2014) The emerging roles of DOT1L in leukemia and normal development. *Leukemia* 28(11): 2131–8.
- Nguyen AT & Zhang Y. (2011) The diverse functions of Dot1 and H3K79 methylation. *Genes Dev* 25(13): 1345–58.
- Rau RE et al. (2016) DOT1L as a therapeutic target for the treatment of DNMT3A-mutant acute myeloid leukemia. *Blood* 128(7): 971–81.
- Vatapalli R et al. (2020) Histone methyltransferase DOT1L coordinates AR and MYC stability in prostate cancer. *Nat Commun* 11(1): 4153.
- Zhou Z et al. (2019) Epigenetically modulated FOXM1 suppresses dendritic cell maturation in pancreatic cancer and colon cancer. *Mol Oncol* 13(4): 873–93.

## Related Products

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2024 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.