

## Anti-Tyrosine Hydroxylase Antibody, Clone TH-2



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Antibodies

Mouse monoclonal IgG1 antibody  
against human, rat, cow tyrosine  
hydroxylase, unconjugated

Catalog #60058

200 µL

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

## Product Description

The TH-2 antibody clone recognizes an epitope present in the N-terminal region of both rodent (~60 kDa) and human (62 - 68 kDa) tyrosine hydroxylase (TH) which catalyzes the hydroxylation of L-tyrosine to L-3,4 dihydroxyphenylalanine (L-dopa) in brain and adrenal medulla and can therefore be used to detect dopaminergic neurons. L-dopa is required for the biosynthesis of catecholamines (dopamine, norepinephrine and epinephrine) which function as neurotransmitters and hormones.

Target Antigen Name:	Tyrosine Hydroxylase
Alternative Names:	Tyrosine 3-monooxygenase, Tyrosine 3-hydroxylase
Gene ID:	7054 (human), 25085 (rat)
Species Reactivity:	Human, Rat, Cow, Guinea pig, Monkey, Rabbit, Sheep
Host Species:	Mouse
Clonality:	Monoclonal
Clone:	TH-2
Isotype:	IgG1
Immunogen:	Rat tyrosine hydroxylase
Conjugate:	Unconjugated

## Applications

Verified:	ICC
Reported:	ELISA, ICC, IF, IHC, IP, WB

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

## Properties

Formulation:	Mouse ascites fluid containing < 0.1% sodium azide
Purification:	The antibody was purified by column chromatography.
Stability and Storage:	Product stable at -20°C when stored undiluted. For product expiry date, please contact <a href="mailto:techsupport@stemcell.com">techsupport@stemcell.com</a> .
Directions for Use:	<p>For immunocytochemistry the suggested concentration of this antibody is 1:400 dilution. It is recommended that the antibody be titrated for optimal performance for each application.</p> <p>For instructions on how to use this antibody, refer to the Technical Manual: In Vitro Proliferation and Differentiation of Human Neural Stem and Progenitor Cells Using NeuroCult™ or NeuroCult™-XF (Document #28724) available on our website at <a href="http://www.stemcell.com">www.stemcell.com</a>.</p>

## Related Products

For a complete list of antibodies, including other conjugates, sizes and clones, as well as related products available from STEMCELL Technologies, please visit our website at [www.stemcell.com/antibodies](http://www.stemcell.com/antibodies) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

1. Beveridge TS et al. (2015) Anatomy of the nerves and ganglia of the aortic plexus in males. *J Anat* 226(1): 93–103. (IHC)
2. Broms J et al. (2015) Conserved expression of the GPR151 receptor in habenular axonal projections of vertebrates. *J Comp Neurol* 523(3): 359–80. (IF, IHC)
3. Nadella R et al. (2014) Transient transfection of human CDNF gene reduces the 6-hydroxydopamine-induced neuroinflammation in the rat substantia nigra. *J Neuroinflammation* 11: 209. (IF, IHC)
4. Arita DY et al. (2002) Purification and characterization of the active form of tyrosine hydroxylase from mesangial cells in culture. *J Cell Biochem* 87(1): 58–64.
5. Nagatsu T et al. (1964) Tyrosine hydroxylase. The initial step in norepinephrine biosynthesis. *J Biol Chem* 239: 2910–7.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and NeuroCult are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.