

## Anti-Human CD2 Antibody, Clone RPA-2.10



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Antibodies

Mouse monoclonal IgG1 antibody  
against human, rhesus, cynomolgus  
CD2, unconjugated

Catalog #60007

100 µg 0.5 mg/mL

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

## Product Description

The RPA-2.10 antibody reacts with CD2, an ~50 kDa type I transmembrane glycoprotein and a member of the immunoglobulin (Ig) superfamily; two Ig-like domains are located in its extracellular portion. CD2 is broadly expressed on peripheral T and NK cells, dendritic cells, erythrocytes, most thymic T cells, subsets of thymic B cells, and on the endothelium. CD2 expression appears early during T cell differentiation. Aberrant expression has been observed in some lymphomas and myeloid leukemias. CD2 is critically important for T cell activation and signaling, and lymphocyte adhesion. The primary ligand for CD2 is CD58 (LFA-3) located on antigen-presenting cells, with additional ligands comprising CD15 (SSEA-1), CD48 and CD59. Notably, the RPA-2.10 antibody blocks the mixed lymphocyte reaction.

|                      |  |
|----------------------|--|
| Target Antigen Name: | CD2  |
| Alternative Names:   | LFA-2, SRBC-R, T11   |
| Gene ID:             | 914  |
| Species Reactivity:  | Human, Rhesus, Cynomolgus, Baboon, Chimpanzee, Capuchin Monkey, Pigtailed Macaque, Pig |
| Host Species:        | Mouse  |
| Clonality:           | Monoclonal   |
| Clone:               | RPA-2.10   |
| Isotype:             | IgG1, kappa  |
| Immunogen:           | Human CD2 recombinant protein  |
| Conjugate:           | Unconjugated   |

## Applications


|                       |   |
|-----------------------|---|
| Verified:             | FC  |
| Reported:             | FA, FC, ICC, IF, IHC, WB  |
| Special Applications: | This antibody clone has been verified for purity assessments of cells isolated with EasySep™ kits, including EasySep™ HLA CD3 Positive Selection Kit (Catalog #18051HLA), EasySep™ Human CD3 Positive Selection Kit (Catalog #18051), EasySep™ HLA Whole Blood CD3 Positive Selection Kit (Catalog #18081HLA) and EasySep™ Human Whole Blood CD3 Positive Selection Kit (Catalog #18081). |

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; WB: Western blotting

## Properties

|                        |  |
|------------------------|--|
| Formulation:           | Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide   |
| Purification:          | The antibody was purified by affinity chromatography.  |
| Stability and Storage: | Product stable at 2 - 8°C when stored undiluted. Do not freeze. For product expiry date, please request a lot-specific Certificate of Analysis from <a href="mailto:techsupport@stemcell.com">techsupport@stemcell.com</a> . |
| Directions for Use:    | For flow cytometry the suggested use of this antibody is ≤ 0.25 µg per 1 × 10 <sup>6</sup> cells in 100 µL volume. It is recommended that the antibody be titrated for optimal performance for each application.             |

## Data



Flow cytometry analysis of human peripheral blood lymphocytes labeled with Anti-Human CD2 Antibody, Clone RPA-2.10, followed by anti-mouse IgG, FITC (solid line histogram) or a mouse IgG1, kappa isotype control antibody followed by anti-mouse IgG, FITC (dashed line histogram).

## Related Products

| PRODUCT NAME                                  | CATALOG # | SIZE      |
|---|-----------|-----------|
| Anti-Human CD2 Antibody, Clone RPA-2.10       | 60007     | 100 µg    |
| Anti-Human CD2 Antibody, Clone RPA-2.10, PE   | 60007PE   | 100 tests |
| Anti-Human CD2 Antibody, Clone RPA-2.10, PE   | 60007PE.1 | 25 tests  |
| Anti-Human CD2 Antibody, Clone RPA-2.10, FITC | 60007FI   | 100 tests |

## References

1. Aversa GG, et al. RPA-2.10: an anti-CD2 monoclonal antibody that inhibits alloimmune responses and monitors T cell activation. *Transplant Proc* 19(1 Pt 1): 277-78, 1987
2. Knapp W, et al. Eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford University Press. New York, pp. 628-34, 1989
3. Shaw AS, et al. Making the T cell receptor go the distance: a topological view of T cell activation. *Immunity* 6(4): 361-69, 1997
4. Scherthaner GH, et al. Expression, epitope analysis, and functional role of the LFA2 antigen detectable on neoplastic mast cells. *Blood* 98(13): 3784-92, 2001
5. Piriou-Guzylack L, Salmon H. Membrane markers of the immune cells in swine: an update. *Vet Res* 39(6): 54, 2008
6. Kap YS, et al. A monoclonal antibody selection for immunohistochemical examination of lymphoid tissues from non-human primates. *J Histochem Cytochem* 57(12): 1159-67, 2009 (IHC)
7. Zaretsky AG, et al. T follicular helper cells differentiate from Th2 cells in response to helminth antigens. *J Exp Med* 206(5): 991-99, 2009 (IHC)
8. Perona-Wright G, et al. Sustained signaling by canonical helper T cell cytokines throughout the reactive lymph node. *Nat Immunol* 11(6): 520-26, 2010 (FC)
9. Thummler K, et al. Immune regulation by peripheral suppressor T cells induced upon homotypic T cell/T cell interactions. *J Leukoc Biol* 88(5): 1041-50, 2010
10. Yoshino N, et al. Upgrading of flow cytometric analysis for absolute counts, cytokines and other antigenic molecules of cynomolgus monkeys (*Macaca fascicularis*) by using anti-human cross-reactive antibodies. *Exp Anim* 49(2): 97-110, 2000 (FC)

Copyright © 2013 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists and EasySep are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.