



EasySep™ Human CXCR4 Positive Selection Kit

Positive Selection

Catalog #18163

For processing 1 x 10⁹ cells



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Document #28068 | Version 1_0_2

Description

Isolate highly purified CXCR4+ cells derived from human pluripotent stem cells by immunomagnetic positive selection.

- Fast and easy-to-use
- Up to 99% purity
- No columns required

This kit targets CXCR4+ cells for positive selection with an antibody recognizing the CXCR4 (CD184) surface marker. The cocktail contains the B-R24 clone of anti-CXCR4 antibody which is known to inhibit HIV infection. Desired cells are labeled with antibodies and magnetic particles, and separated without columns using an EasySep™ magnet. Unwanted cells are simply poured off, while desired cells remain in the tube. Isolated cells are immediately available for downstream applications such as flow cytometry, culture, or DNA/RNA extraction.

Component Descriptions

COMPONENT NAME	COMPONENT #	QUANTITY	STORAGE	SHELF LIFE	FORMAT
EasySep™ Human CXCR4 Positive Selection Cocktail	18163C	1 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A combination of monoclonal antibodies in PBS.
EasySep™ Dextran RapidSpheres™ 50100	50100	1 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A suspension of magnetic particles in water.

PBS - phosphate-buffered saline

Components may be shipped at room temperature (15 - 25°C) and should be stored as indicated above.

Sample Preparation

Trypsin-EDTA is recommended for preparation of a single-cell suspension from human pluripotent stem cell-derived definitive endoderm.

1. Discard culture medium and rinse each well of a 6-well plate with 2 mL of DMEM/F-12 with 15 mM HEPES (Catalog #36254).
2. Dilute Trypsin-EDTA (0.25%; Catalog #07901) 10-fold in HBSS, Modified (Without Ca++ and Mg++; Catalog #37250).
3. Add 1 mL of diluted Trypsin-EDTA to each well.
4. Incubate at 37°C for 5 minutes until the cells are easily dislodged from the bottom of the plate with gentle tapping.
5. Pipette up and down 3 - 4 times with a 1 mL pipette tip to dissociate to single cells.
6. Combine the single-cell suspension from each well into a 15 mL conical tube (e.g. Catalog #38009).
7. Wash each well with 1 mL of Soybean Trypsin Inhibitor, ACF (Catalog #07457) and add to the 15 mL conical tube.
8. Centrifuge at 300 x g for 5 minutes.
9. Remove supernatant and resuspend the cells at 5 x 10⁷ cells/mL in recommended medium.


Recommended Medium

RoboSep™ Buffer 2 (Catalog #20164), or PBS containing 0.5% bovine serum albumin (BSA) and 2 mM EDTA. Medium should be free of Ca++ and Mg++.

Directions for Use

See page 1 for Sample Preparation and Recommended Medium. Refer to Table 1 for detailed instructions regarding the EasySep™ procedure.

Table 1. EasySep™ Human CXCR4 Positive Selection Kit Protocol

		EASYSEP™ MAGNETS	
STEP	INSTRUCTIONS	EasySep™ (Catalog #18000)	
1	Prepare sample at the indicated cell concentration within the volume range.	5 x 10 ⁷ cells/mL 0.1- 0.3 mL NOTE: If starting with fewer than 5 x 10 ⁶ cells, resuspend cells in 0.1 mL.	
	Add sample to required tube.	5 mL (12 x 75 mm) polystyrene round-bottom tube (e.g. Catalog #38007)	
2	Add Selection Cocktail to sample.	50 µL/mL of sample	
	Mix and incubate.	RT for 5 minutes	
3	Vortex RapidSpheres™. NOTE: Particles should appear evenly dispersed.	30 seconds	
4	Add RapidSpheres™ to sample.	25 µL/mL of sample	
	Mix and incubate.	RT for 5 minutes	
5	Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL	
	Place the tube (without lid) into the magnet and incubate.	RT for 5 minutes	
6	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring off the supernatant. Remove the tube from the magnet; the tube contains the isolated cells.	Discard supernatant	
7	Repeat steps as indicated.	Steps 5 and 6, two more times (total of 3 x 5-minute separations)	
8	Resuspend cells in desired medium. Be sure to collect cells from the sides of the tube.	Isolated cells are ready for use	

RT - room temperature (15 - 25°C)

* Leave the magnet and tube inverted for 2 - 3 seconds, then return upright. Do not shake or blot off any drops that may remain hanging from the mouth of the tube.

Notes and Tips

ASSESSING PURITY


EasySep™ Human CXCR4 Positive Selection Cocktail uses the anti-CXCR4 antibody clone B-R24. For purity assessment by flow cytometry use the following fluorochrome-conjugated antibody clones:

- Anti-Human CD184 (CXCR4) Antibody, Clone 12G5 (partial blocking; Catalog #60089)
- Anti-Human CD117 (c-Kit) Antibody, Clone 104D2 (Catalog #60087)
- Anti-human SOX17 antibody

The following method can also be used:

- Use an alternative fluorochrome-conjugated antibody marker.

Data



Starting with human pluripotent stem cell-derived definitive endoderm, the CXCR4+SOX17+ cell content of the isolated fraction is typically $92.9 \pm 3.1\%$ (mean \pm SD). In the above example, the purities of the start and final isolated fractions are 24.9% and 95.5%, respectively.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, RapidSpheres, and RoboSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.