

## Small Molecules

LY294002

PI3K/AKT pathway inhibitor; Inhibits PI3K

Catalog # 72152  
72154

5 mg  
25 mg



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

LY294002 is a phosphoinositide 3-kinase (PI3K) inhibitor that has greater potency and selectivity than quercetin, the structure on which it is based. LY294002 inhibits PI3K ( $IC_{50} = 1.4 \mu M$ ) and also shows activity against casein kinase 2 (CK2), but not PI4K, EGFR, PDGFR, MAPK, PKA, or PKC (Davies et al.; Vlahos et al.).

Molecular Name:	LY294002
Alternative Names:	Not applicable
CAS Number:	154447-36-6
Chemical Formula:	$C_{19}H_{17}NO_3$
Molecular Weight:	307.3 g/mol
Purity:	$\geq 98\%$
Chemical Name:	2-(4-morpholinyl)-8-phenyl-4H-1-benzopyran-4-one
Structure:	



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at $-20^{\circ}C$ as supplied. Protect from prolonged exposure to light. For product expiry date, please contact <a href="mailto:techsupport@stemcell.com">techsupport@stemcell.com</a> .
Solubility:	<ul style="list-style-type: none"><li>· DMSO <math>\leq 50</math> mM</li><li>· Absolute ethanol <math>\leq 50</math> mM</li><li>· PBS (pH 7.2) <math>\leq 160 \mu M</math></li></ul> For example, to prepare a 10 mM stock solution in DMSO, resuspend 5 mg in 1.63 mL of fresh DMSO.

Prepare stock solution fresh before use. Stock solutions in DMSO or ethanol are stable for up to 6 months if stored at  $-20^{\circ}C$ . Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### MAINTENANCE AND SELF-RENEWAL

- Suppresses proliferation and self-renewal of mouse embryonic stem (ES) cells (Lianguzova et al.; Paling et al.).

### DIFFERENTIATION

- Promotes differentiation to insulin-producing cells from mouse ES cells (Hori et al.).
- Inhibits myotube formation from myoblasts (Coolican et al.; Jiang et al.).

## References

- Coolican SA et al. (1997) The mitogenic and myogenic actions of insulin-like growth factors utilize distinct signaling pathways. *J Biol Chem* 272(10): 6653–62.
- Davies SP et al. (2000) Specificity and mechanism of action of some commonly used protein kinase inhibitors. *Biochem J* 351(Pt 1): 95–105.
- Hori Y et al. (2002) Growth inhibitors promote differentiation of insulin-producing tissue from embryonic stem cells. *Proc Natl Acad Sci USA* 99(25): 16105–10.
- Jiang BH et al. (1998) An essential role of phosphatidylinositol 3-kinase in myogenic differentiation. *Proc Natl Acad Sci USA* 95(24): 14179–83.
- Lianguzova MS et al. (2007) Phosphoinositide 3-kinase inhibitor LY294002 but not serum withdrawal suppresses proliferation of murine embryonic stem cells. *Cell Biol Int* 31(4): 330–7.
- Paling NRD et al. (2004) Regulation of embryonic stem cell self-renewal by phosphoinositide 3-kinase-dependent signaling. *J Biol Chem* 279(46): 48063–70.
- Vlahos CJ et al. (1994) A specific inhibitor of phosphatidylinositol 3-kinase, 2-(4-morpholinyl)-8-phenyl-4H-1-benzopyran-4-one (LY294002). *J Biol Chem* 269(7): 5241–8.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, please visit our website at [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.