

PRODUCT DESCRIPTION

StemAdhere™ Defined Matrix for hESCs and hiPSCs is a matrix that supports the culture of human embryonic stem (ES) cells and human induced pluripotent stem (iPS) cells. It contains a single recombinant protein composed of entirely human sequences. It can be used in conjunction with TeSR™2 feeder-free maintenance medium (Catalog #05860) for long term passaging while maintaining a pluripotent state.

StemAdhere™ Defined Matrix for hESCs and hiPSCs is developed and manufactured by Primorigen Biosciences®, Inc.


COMPONENTS

07160 StemAdhere™ 250 µg/mL, 2 mL

StemAdhere™ contains a single recombinant protein composed of entirely human sequences.

STABILITY AND STORAGE

Product stable at -80°C or -20°C until expiry date as indicated on label. Thaw entire vial at room temperature (15 - 25°C) or 2 - 8°C, prior to use. Store thawed StemAdhere™ at 2 - 8°C and use within one week. Product can be aliquoted immediately after thaw and refrozen at -80°C or -20°C. Avoid repeated freezing and thawing.



DIRECTIONS FOR USE

For complete instructions on maintaining high quality human ES cells and iPS cells using StemAdhere™, please see the Technical Manual Maintenance of human pluripotent stem cells on StemAdhere™ (Document #28795) available on our website at www.stemcell.com.

RELATED PRODUCTS

PRODUCT	CATALOG #
StemAdhere™ Kit	07170
StemAdhere™ Dilution Buffer	07163
Gentle Cell Dissociation Reagent	07174
Non Tissue Culture-Treated 6-well Plates	27147

Copyright © 2013 by STEMCELL Technologies, Inc. All rights reserved, including graphics and images. STEMCELL Technologies & Design, STEMCELL shield, Scientists Helping Scientists are trademarks of STEMCELL Technologies, Inc. StemAdhere™ is developed and manufactured by Primorigen Biosciences®, Inc., and StemAdhere™ is a trademark of Primorigen Biosciences®, Inc. TeSR is a trademark of WARF. All other trademarks and registered trademarks are the property of their respective holders.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.
FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.



TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM • FOR GLOBAL CONTACT DETAILS VISIT WWW.STEMCELL.COM

VERSION 1.2.0
DOCUMENT #28796