

SpinSep™ Density Medium

Density gradient medium

Catalog # 17531

100 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

SpinSep™ Density Medium is a density gradient medium that is designed for the enrichment of specific human cell subsets.

Properties

- Storage:** Store at 15 - 25°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- Iodixanol
 - Hetastarch
 - Sodium chloride
 - Sodium lactate
 - Dextrose
 - Calcium chloride
 - Potassium chloride
 - Magnesium chloride
 - Water

Handling / Directions For Use

Storage at 2 - 8°C is acceptable; ensure that the medium warms to room temperature (15 - 25°C) and invert bottle to mix contents before use.

Use with EasySep™ Human Myeloid Positive Selection Kit (Catalog #18653). Detailed information is included in the Product Information Sheet (Document #29009), available at www.stemcell.com or contact us to request a copy.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, and SpinSep are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.