

## Small Molecules

### Garcinol

Epigenetic modifier; Inhibits histone acetyltransferases (HATs) p300 and pCAF

Catalog # 72452

5 mg



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Garcinol, a polyisoprenylated benzophenone isolated from *Garcinia indica*, is an inhibitor of the histone acetyltransferases (HATs) p300 and pCAF ( $IC_{50} = 7$  and  $5 \mu M$ , respectively; Balasubramanyam et al.). It also inhibits the HAT Gcn5 in *Cryptococcus neoformans*, inducing temperature sensitivity and impairing growth (O'Meara et al.).

Molecular Name:	Garcinol
Alternative Names:	Camboginol
CAS Number:	78824-30-3
Chemical Formula:	$C_{38}H_{50}O_6$
Molecular Weight:	602.8 g/mol
Purity:	$\geq 95\%$
Chemical Name:	3-(3,4-dihydroxybenzoyl)-4-hydroxy-8,8-dimethyl-1,7-bis(3-methyl-2-buten-1-yl)-5-[(2S)-5-methyl-2-(1-methylethenyl)-4-hexen-1-yl]-bicyclo[3.3.1]non-3-ene-2,9-dione

Structure:



## Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at $-20^{\circ}C$ as supplied. Protect from prolonged exposure to light. Stable as supplied for 12 months from date of receipt.
Solubility:	· Absolute ethanol $\leq 30$ mM · DMSO $\leq 30$ mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 166 $\mu L$ of fresh DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at  $-20^{\circ}C$ . Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### MAINTENANCE AND SELF-RENEWAL

- Promotes ex vivo expansion of human hematopoietic stem cells (Nishino et al.).

### DIFFERENTIATION

- Promotes neurogenesis in rat cortical progenitor cells (Weng et al.).

### CANCER RESEARCH

- Induces apoptosis in several types of cancer cells and has anti-inflammatory actions (Koeberle et al.; Prasad et al.).

## References

- Balasubramanyam K et al. (2004) Polyisoprenylated benzophenone, garcinol, a natural histone acetyltransferase inhibitor, represses chromatin transcription and alters global gene expression. *J Biol Chem* 279(32): 33716–26.
- Koeberle A et al. (2009) Identification of 5-lipoxygenase and microsomal prostaglandin E2 synthase-1 as functional targets of the anti-inflammatory and anti-carcinogenic garcinol. *Biochem Pharmacol* 77(9): 1513–21.
- Nishino T et al. (2011) Ex vivo expansion of human hematopoietic stem cells by garcinol, a potent inhibitor of histone acetyltransferase. P. Rameshwar (Ed.). *PLoS One* 6(9): e24298.
- O'Meara TR et al. (2010) *Cryptococcus neoformans* histone acetyltransferase Gcn5 regulates fungal adaptation to the host. *Eukaryot Cell* 9(8): 1193–202.
- Prasad S et al. (2010) Garcinol potentiates TRAIL-induced apoptosis through modulation of death receptors and antiapoptotic proteins. *Mol Cancer Ther* 9(4): 856–68.
- Weng M-S et al. (2011) Garcinol promotes neurogenesis in rat cortical progenitor cells through the duration of extracellular signal-regulated kinase signaling. *J Agric Food Chem* 59(3): 1031–40.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.