

Small Molecules

Kartogenin

RUNX1 transcriptional activator; Binds filamin A

Catalog # 72572

5 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Kartogenin induces chondrogenesis by binding the actin-binding protein filamin A, which disrupts its interaction with the transcription factor core-binding factor β subunit (CBF β). When dissociated from filamin A, CBF β translocates to the nucleus and forms a transcriptional complex with the runt-related transcription factor RUNX1, which enables chondrocyte differentiation (Johnson et al.).

| | |
|--------------------|---|
| Molecular Name: | Kartogenin |
| Alternative Names: | Not applicable |
| CAS Number: | 4727-31-5 |
| Chemical Formula: | C ₂₀ H ₁₅ NO ₃ |
| Molecular Weight: | 317.3 g/mol |
| Purity: | ≥ 98% |
| Chemical Name: | 2-[(4-phenylphenyl)carbamoyl]benzoic acid |
| Structure: | |



Properties

| | |
|----------------------|--|
| Physical Appearance: | A crystalline solid |
| Storage: | Product stable at -20°C as supplied. Protect from prolonged exposure to light. Stable as supplied for 12 months from date of receipt. |
| Solubility: | · DMSO ≤ 90 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 315 μ L of DMSO. |

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

DIFFERENTIATION

- Promotes differentiation of human bone marrow mesenchymal stem cells into chondrocytes (Johnson et al.; Zhang et al.).
- Promotes cartilage formation/repair in mouse models of osteoarthritis or when injected into mouse tendon-bone junctions (Johnson et al.; Zhang et al.).
- Promotes type-I collagen synthesis in human dermal fibroblasts in vitro and in the dermis of mice through activation of the SMAD4/SMAD5 pathway (Wang et al.).

References

Johnson K et al. (2012) A stem cell-based approach to cartilage repair. *Science* 336(6082): 717–21.

Wang J et al. (2014) A heterocyclic molecule kartogenin induces collagen synthesis of human dermal fibroblasts by activating the smad4/smad5 pathway. *Biochem Biophys Res Commun* 450(1): 568–74.

Zhang J & Wang JH-C. (2014) Kartogenin induces cartilage-like tissue formation in tendon–bone junction. *Bone Res* 2: 14008.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.