

## Small Molecules

AM580

Retinoid pathway activator; Activates retinoic acid receptor (RAR) alpha

Catalog # 72962  
72964

1 mg  
10 mg



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

AM580 is a retinoic acid receptor (RAR) agonist that is selective for RAR $\alpha$  ( $EC_{50}$  = 0.36 nM) compared to RAR $\beta$  ( $EC_{50}$  = 24.6 nM) and RAR $\gamma$  ( $EC_{50}$  = 27.9 nM; Bernard et al.) It is a derivative of retinoic acid (RA), however it demonstrates greater specific binding to RAR $\alpha$  compared to RA, which exhibits little selectivity across RAR $\alpha$ ,  $\beta$ , or  $\gamma$  (Gianní et al.; Bernard et al.; Kim et al; Rochette-Egly & Germain).

**Molecular Name:** AM580  
**Alternative Names:** CD336; NSC 608001; Ro 40-6055  
**CAS Number:** 102121-60-8  
**Chemical Formula:** C<sub>22</sub>H<sub>25</sub>NO<sub>3</sub>  
**Molecular Weight:** 351.4 g/mol  
**Purity:**  $\geq$  98%  
**Chemical Name:** 4-[[[(5,6,7,8-tetrahydro-5,5,8,8-tetramethyl-2-naphthalenyl)carbonyl]amino]-benzoic acid  
**Structure:**



## Properties

**Physical Appearance:** A crystalline solid  
**Storage:** Product stable at -20°C as supplied. Protect from prolonged exposure to light. For product expiry date, please contact [techsupport@stemcell.com](mailto:techsupport@stemcell.com).  
**Solubility:** · DMSO  $\leq$  55 mM  
· Absolute ethanol  $\leq$  25 mM  
For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 285  $\mu$ L of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

## Published Applications

### REPROGRAMMING

- Promotes reprogramming of somatic cells to induced pluripotent stem cells (Wang et al.).

### DIFFERENTIATION

- Induces differentiation of human induced pluripotent stem cells into intermediate mesoderm, in combination with the GSK3 $\beta$  inhibitor CHIR99021 (Catalog #72052; Araoka et al.).

### CANCER RESEARCH

- Inhibits tumor cell proliferation and survival signaling pathways, and induces apoptosis, leading to inhibition of mouse mammary tumor virus (MMTV)-neu- and MMTV-wnt1-induced mammary gland hyperplasia (Lu et al.).
- Inhibits tumor growth in MMTV-Myc mice (Bosch et al.).
- Inhibits endometrial cancer cell proliferation (Cheng et al.).
- Induces differentiation in acute promyelocytic leukemia cells (Gianní et al.).

## References

- Araoka T et al. (2014) Efficient and Rapid Induction of Human iPSCs/ESCs into Nephrogenic Intermediate Mesoderm Using Small Molecule-Based Differentiation Methods. PLoS One 9(1): e84881.
- Bernard BA et al. (1992) Identification of synthetic retinoids with selectivity for human nuclear retinoic acid receptor  $\gamma$ . Biochem Biophys Res Commun 186(2): 977–83.
- Bosch A et al. (2012) Reversal by RAR $\alpha$  agonist Am580 of c-Myc-induced imbalance in RAR $\alpha$ /RAR $\gamma$  expression during MMTV-Myc tumorigenesis. Breast Cancer Res 14(4): R121.
- Cheng Y-H et al. (2011) Retinoic acid inhibits endometrial cancer cell growth via multiple genomic mechanisms. J Mol Endocrinol 46(2): 139–53.
- Gianní M et al. (1996) AM580, a stable benzoic derivative of retinoic acid, has powerful and selective cyto-differentiating effects on acute promyelocytic leukemia cells. Blood 87(4): 1520–31.
- Kim M-J et al. (2000) The Role of Specific Retinoid Receptors in Sebocyte Growth and Differentiation in Culture1. J Invest Dermatol 114(2): 349–53.
- Lu Y et al. (2010) Mechanism of inhibition of MMTV-neu and MMTV-wnt1 induced mammary oncogenesis by RARalpha agonist AM580. Oncogene 29(25): 3665–76.
- Rochette-Egly C & Germain P. (2009) Dynamic and combinatorial control of gene expression by nuclear retinoic acid receptors (RARs). Nucl Recept Signal 7: e005.
- Wang W et al. (2011) Rapid and efficient reprogramming of somatic cells to induced pluripotent stem cells by retinoic acid receptor gamma and liver receptor homolog 1. Proc Natl Acad Sci USA 108(45): 18283–8.

## Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, please visit our website at [www.stemcell.com/smallmolecules](http://www.stemcell.com/smallmolecules) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.