

Antibodies

Goat Anti-Rat IgM (Heavy Chain) Antibody, Polyclonal, APC

Goat polyclonal antibody against rat IgM (heavy chain), APC-conjugated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #60140AZ

250 µg

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

APC-conjugated goat anti-rat IgM, μ chain specific, F(ab')₂ fragment antibody reacts with the heavy chain of rat IgM but not with rat IgG. This information on reactivity was based on immunoelectrophoresis and/or ELISA data. This antibody has minimal cross-reaction with human, cow, and horse serum proteins, but it may cross-react with IgM from other species. The antibody was pepsin digested to obtain the F(ab')₂ fragment. Fc fragments and whole IgG molecules have been removed.

Target Antigen Name:	IgM (Heavy Chain), Rat
Alternative Names:	Not applicable
Gene ID:	Not applicable
Species Reactivity:	Rat IgM, μ chain specific. Minimal cross reactivity to human, cow, and horse serum proteins
Host Species:	Goat
Clonality:	Polyclonal
Clone:	Not applicable
Isotype:	Not applicable
Immunogen:	Not applicable
Conjugate:	APC

Applications

Verified:	FC
Reported:	FC, ICC, IF, IHC

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Formulation:	Lyophilized from a solution containing sodium phosphate, sodium chloride, bovine serum albumin, and sodium azide
Purification:	The antibody was prepared by a combination of pepsin digestion and antigen affinity chromatography. The preparation is free of Fc fragments and whole IgG molecules.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. For product expiry date, please contact techsupport@stemcell.com .
Directions for Use:	<p>Centrifuge vial before opening. Resuspend the product in 0.5 mL deionized water, this is the stock dilution. Centrifuge if solution is not clear. Prepare working dilution fresh each day.</p> <p>NOTE: Once resuspended, store stock dilution at 2 - 8°C and use within 6 months. Do not freeze.</p> <p>For flow cytometry the suggested concentration of this antibody is 1:100 dilution. It is recommended that the antibody be titrated for optimal performance for each application.</p>

Related Products

For a complete list of antibodies, including other conjugates, sizes and clones, as well as related products available from STEMCELL Technologies, please visit our website at www.stemcell.com/antibodies or contact us at techsupport@stemcell.com.

References

1. Kummitha CM et al. (2012) HECA-452 is a non-function blocking antibody for isolated sialyl Lewis x adhesion to endothelial expressed E-selectin under flow conditions. *J Immunol Methods* 384(1-2): 43–50. (FC)
2. Barbaric I et al. (2011) High-content screening for chemical modulators of embryonal carcinoma cell differentiation and survival. *J Biomol Screen* 16(6): 603–17. (ICC, IF)
3. Lawson BR et al. (2007) Immunomodulation of murine collagen-induced arthritis by N, N-dimethylglycine and a preparation of Perna canaliculus. *BMC Complement Altern Med* 7(1): 20. (FA)
4. Haspot F et al. (2005) Anti-CD28 antibody-induced kidney allograft tolerance related to tryptophan degradation and TCR- class II- B7+ regulatory cells. *Am J Transplant* 5(10): 2339–48. (IF, IHC)
5. Johansson-Lindbom B et al. (2003) Selective generation of gut tropic T cells in gut-associated lymphoid tissue (GALT): Requirement for GALT dendritic cells and adjuvant. *J Exp Med* 198(6): 963–9. (FC)

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.