

## Anti-Biotin Antibody, Clone C6D5.1.1

## Antibodies

Mouse monoclonal IgG1 antibody  
against biotin, unconjugated



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #60145

100 µg 1 mg/mL

## Product Description

The C6D5.1.1 antibody reacts with biotin (vitamin H, vitamin B7, coenzyme R), a 244.3 Da protein that is present in all living cells and is essential for cell growth and the citric acid cycle. Biotin is commonly coupled to proteins and other molecules without altering their biological activity for use in a number of applications. Both avidin and streptavidin form strong non-covalent bonds with biotin in a 1:4 ratio and thus efficiently cross-link biotinylated molecules. Both avidin and streptavidin can, however, exhibit nonspecific binding that can be problematic in certain applications. Cross-linking of biotinylated molecules, such as immunoglobulins, can also be achieved using anti-biotin antibodies such as clone C6D5.1.1. This approach may be preferable to eliminate the non-specific binding problems associated with streptavidin, or when a lesser degree of cross-linking is required.

Target Antigen Name:	Biotin
Alternative Names:	Coenzyme R, Vitamin B7, Vitamin H
Gene ID:	Not applicable
Species Reactivity:	Not applicable
Host Species:	Mouse (BALB/c)
Clonality:	Monoclonal
Clone:	C6D5.1.1
Isotype:	IgG1
Immunogen:	Biotin conjugated to keyhole limpet hemocyanin
Conjugate:	Unconjugated

## Applications

Verified:	CellSep
Reported:	CellSep, Cross-linking, ELISA

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

## Properties

Formulation:	Phosphate-buffered saline
Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. Addition of 0.1% sodium azide (final) is recommended once the vial has been opened. For product expiry date, please contact <a href="mailto:techsupport@stemcell.com">techsupport@stemcell.com</a> .
Directions for Use:	It is recommended that the antibody be titrated for optimal performance for each application.

## Related Products

For a complete list of antibodies, including other conjugates, sizes and clones, as well as related products available from STEMCELL Technologies, please visit our website at [www.stemcell.com/antibodies](http://www.stemcell.com/antibodies) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.