

Anti-Human CD25 Antibody, Clone 2A3, PE

Antibodies

Mouse monoclonal IgG1 antibody
against human, cynomolgus CD25, PE-
conjugated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #60153PE

50 tests 20 µL/test

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

The 2A3 antibody reacts with human CD25, a 55 kDa type I integral transmembrane glycoprotein also known as low-affinity interleukin-2 receptor alpha chain (IL-2R α) or Tac antigen. CD25 is present on activated T and B cells, monocytes, macrophages, and regulatory T cells. CD25 has a role in T-lymphocyte activation and can be used in the identification of T-cell subsets in viral infections like HIV and HTLV. CD25 associates with IL-2 receptor β (CD122) and common γ (CD132) chains to form a high affinity IL-2R complex. The 2A3 antibody binds to an epitope on p55 at or near the IL-2 binding site and prevents high-affinity IL-2 binding to T lymphocytes. Daclizumab, a therapeutic humanised monoclonal antibody, blocks the binding of 2A3 antibody.

Target Antigen Name:	CD25
Alternative Names:	Interleukin 2 receptor alpha (IL-R2 α), Low affinity IL-2R, p55, Tac
Gene ID:	100
Species Reactivity:	Human, Cynomolgus, Common marmoset, Pigtailed macaque, Squirrel monkey
Host Species:	Mouse (BALB/c)
Clonality:	Monoclonal
Clone:	2A3
Isotype:	IgG1, kappa
Immunogen:	Human phytohemagglutinin-activated T cells
Conjugate:	PE

Applications

Verified:	FC
Reported:	FACS, FC

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; RIA: Radioimmunoassay; WB: Western blotting

Properties

Formulation:	Phosphate-buffered saline containing 0.1% sodium azide and gelatin
Purification:	The antibody was purified by column chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. Protect product from prolonged exposure to light. For product expiry date, please contact techsupport@stemcell.com .
Directions for Use:	For flow cytometry the suggested use of this antibody is 20 µL per 1 x 10 ⁶ cells in 100 µL volume. It is recommended that the antibody be titrated for optimal performance for each application.

Related Products

For a complete list of antibodies, including other conjugates, sizes and clones, as well as related products available from STEMCELL Technologies, please visit our website at www.stemcell.com/antibodies or contact us at techsupport@stemcell.com.

References

1. Rech AJ et al. (2012) CD25 blockade depletes and selectively reprograms regulatory T cells in concert with immunotherapy in cancer patients. *Sci Transl Med* 4(134): 134ra62. (FC)
2. Ahmadzadeh M & Rosenberg SA. (2006) IL-2 administration increases CD4+ CD25(hi) Foxp3+ regulatory T cells in cancer patients. *Blood* 107(6): 2409–14. (FACS)
3. Krueger JG et al. (2000) Successful in vivo blockade of CD25 (high-affinity interleukin 2 receptor) on T cells by administration of humanized anti-Tac antibody to patients with psoriasis. *J Am Acad Dermatol* 43(3): 448–58. (FC)
4. Anasetti C et al. (1990) A phase I-II study evaluating the murine anti-IL-2 receptor antibody 2A3 for treatment of acute graft-versus-host disease. *Transplantation* 50(1): 49–54.
5. Urdal DL et al. (1984) Purification and chemical characterization of the receptor for interleukin 2 from activated human T lymphocytes and from a human T-cell lymphoma cell line. *Proc Natl Acad Sci USA* 81(20): 6481–5.

Please refer to the Safety Data Sheet (SDS) for hazard information.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.