

Small Molecules

Chaetocin

Epigenetic modifier; Inhibits histone methyltransferase SU(VAR)3-9

Catalog # 73592

500 µg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Chaetocin is a fungal mycotoxin, originally produced by Chaetomium species, that inhibits the lysine-specific histone methyltransferases (HMTs) SU(VAR)3-9 ($IC_{50} = 0.8 \mu M$). HMTs are responsible for methylation of histones, thereby affecting heterochromatinization. Chaetocin is also known to inhibit G9a and DIM5 ($IC_{50} = 2.5 \mu M$ and $3 \mu M$ respectively; Cherblanc et al.; Greiner et al.). Chaetocin also acts as a competitive substrate and inhibitor of the central oxidative stress remediation enzyme thioredoxin reductase-1 (TrxR1; $K_m = 4.6 \mu M$) more potently than glutathione reductase or thioredoxin (Tibodeau et al.).

| | |
|--------------------|---|
| Molecular Name: | Chaetocin |
| Alternative Names: | Chaetocin A |
| CAS Number: | 28097-03-2 |
| Chemical Formula: | $C_{30}H_{28}N_6O_6S_4$ |
| Molecular Weight: | 696.8 g/mol |
| Purity: | ≥ 95% |
| Chemical Name: | (3S,3'S,5aR,5aR',10bR,10'bR,11aS,11'aS)-2,2',3,3',5a,5'a,6,6'-octahydro-3,3'-bis(hydroxymethyl)-2,2'-dimethyl-[10b,10'b(11H,11'H)-bi3,11a-epidithio-11aH-pyrazino[1',2':1,5]pyrrolo[2,3-b]indole]-1,1',4,4'-tetrone |

Structure:



Properties

| | |
|----------------------|--|
| Physical Appearance: | A crystalline solid |
| Storage: | Product stable at $-20^{\circ}C$ as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. For product expiry date, please contact techsupport@stemcell.com . |
| Solubility: | · DMSO ≤ 20 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 0.5 mg in 72 µL of DMSO. |

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at $-20^{\circ}C$. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

CANCER RESEARCH

- Induces cellular oxidative stress, selectively killing cancer cells and rapidly-proliferating primary cells (Isham et al.; Spannhoff et al.).
- Attenuates the growth of glioma xenografts accompanied by an increase in reactive oxygen species (ROS) production (Dixit et al.).

References

- Cherblanc FL et al. (2013) Chaetocin is a nonspecific inhibitor of histone lysine methyltransferases. *Nat Chem Biol* 9(3): 136–7.
- Dixit D et al. (2014) Chaetocin-induced ROS-mediated apoptosis involves ATM-YAP1 axis and JNK-dependent inhibition of glucose metabolism. *Cell Death Dis* 5: e1212.
- Greiner D et al. (2005) Identification of a specific inhibitor of the histone methyltransferase SU(VAR)3-9. *Nat Chem Biol* 1(3): 143–5.
- Isham CR et al. (2007) Chaetocin: a promising new antimaloma agent with in vitro and in vivo activity mediated via imposition of oxidative stress. *Blood* 109(6): 2579–88.
- Spannhoff A et al. (2009) Cancer treatment of the future: inhibitors of histone methyltransferases. *Int J Biochem Cell Biol* 41(1): 4–11.
- Tibodeau JD et al. (2009) The anticancer agent chaetocin is a competitive substrate and inhibitor of thioredoxin reductase. *Antioxid Redox Signal* 11(5): 1097–106.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, please visit our website at www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.