

Small Molecules

EX527

Epigenetic modifier; Inhibits SIRT1 histone deacetylase

Catalog # 73652
73654

1 mg
10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

EX527 is a cell-permeable, selective inhibitor of mammalian sirtuin 1 (SIRT1; $IC_{50} = 98 \text{ nM}$) over SIRT2 and SIRT3 and has no effect on other histone deacetylases (HDACs; Nayagam et al.). SIRT1 is a nicotinamide adenine dinucleotide (NAD)-dependent deacetylase with roles in energy metabolism and inflammation. Studies have shown that EX527 inhibits sirtuins by forming a trimeric sirtuin complex with an NAD⁺-derived coproduct (Gertz M et al.).

| | |
|--------------------|--|
| Molecular Name: | EX527 |
| Alternative Names: | Selisistat |
| CAS Number: | 49843-98-3 |
| Chemical Formula: | $C_{13}H_{13}ClN_2O$ |
| Molecular Weight: | 248.7 g/mol |
| Purity: | ≥ 98% |
| Chemical Name: | 6-chloro-2,3,4,9-tetrahydro-1H-carbazole-1-carboxamide |
| Structure: | |



Properties

| | |
|----------------------|---|
| Physical Appearance: | A crystalline solid |
| Storage: | Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. For product expiry date, please contact techsupport@stemcell.com . |
| Solubility: | <ul style="list-style-type: none">· DMSO ≤ 80 mM· Ethanol ≤ 20 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 402 µL of DMSO. |

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

DIFFERENTIATION

- Increases the production of oligodendrocytes from differentiating neural stem cells and neural progenitor cells in vitro (Rafalski et al.).

IMMUNOLOGY

- Restores the microvascular response during the hypoinflammatory phase in a mouse model of sepsis, and enhances the systemic innate immune response (Vachharajani et al.).

DISEASE MODELING

- Delays cyst growth in kidneys of PKD1 knockout mouse models (Zhou et al.).

References

Gertz M et al. (2013) Ex-527 inhibits Sirtuins by exploiting their unique NAD⁺-dependent deacetylation mechanism. *Proc Natl Acad Sci U S A* 110(30): E2772–81.

Nayagam VM et al. (2006) SIRT1 modulating compounds from high-throughput screening as anti-inflammatory and insulin-sensitizing agents. *J Biomol Screen* 11(8): 959–67.

Rafalski VA et al. (2013) Expansion of oligodendrocyte progenitor cells following SIRT1 inactivation in the adult brain. *Nat Cell Biol* 15(6): 614–24.

Vachharajani VT et al. (2014) SIRT1 inhibition during the hypoinflammatory phenotype of sepsis enhances immunity and improves outcome. *J Leukoc Biol* 96(5): 785–96.

Zhou X et al. (2013) Sirtuin 1 inhibition delays cyst formation in autosomal-dominant polycystic kidney disease. *J Clin Invest* 123(7): 3084–98.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, please visit our website at www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. All other trademarks are While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.