

Small Molecules

(Z)-Guggulsterone

Retinoic acid receptor (RAR) pathway inhibitor; Inhibits farnesoid X receptor (FXR)

1 mg

Catalog # 73702



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

(Z)-Guggulsterone is a plant steroid found in the resin of the guggul plant *Commiphora mukul* that acts as a selective antagonist of farnesoid X receptor (FXR; Cui et al.). It decreases chenodeoxycholic acid (CDCA)-induced FXR activation ($IC_{50} = 10 \mu M$) in the presence of $100 \mu M$ CDCA (Urizar et al.; Cui et al.).

Molecular Name:	(Z)-Guggulsterone
Alternative Names:	Not applicable
CAS Number:	39025-23-5
Chemical Formula:	$C_{21}H_{28}O_2$
Molecular Weight:	312.5 g/mol
Purity:	$\geq 95\%$
Chemical Name:	(8R,9S,10R,13S,14S,17Z)-17-ethylidene-10,13-dimethyl-1,2,6,7,8,9,11,12,14,15-decahydrocyclopenta[a]phenanthrene-3,16-dione

Structure:



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at $-20^{\circ}C$ as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. For product expiry date, please contact techsupport@stemcell.com .
Solubility:	<ul style="list-style-type: none">· DMSO $\leq 800 \mu M$· Absolute ethanol $\leq 3.2 \mu M$· DMF $\leq 30 mM$ For example, to prepare a 10 mM stock solution in DMF, resuspend 1 mg in 320 μL of DMF.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at $-20^{\circ}C$. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

DIFFERENTIATION

- Induces differentiation of human pluripotent stem cell-derived neural stem cells into dopaminergic neurons (Gonzalez et al.).

METABOLISM

- Decreases hepatic low-density lipoprotein cholesterol and triglyceride levels in rodents fed with a high cholesterol diet (Urizar et al.).

CANCER RESEARCH

- Regulates growth and metastasis of tumor cells in human multiple myeloma cells by the expression of STAT3-regulated antiapoptotic (Bcl-2, Bcl-xL, and Mcl-1), proliferative (cyclin D1), and angiogenic (VEGF) gene products (Ahn et al.).
- Inhibits tube formation in HUVEC and migration in HUVEC or a human prostate cancer cell line (Xiao et al.).

References

- Ahn KS et al. (2008) Guggulsterone, a farnesoid X receptor antagonist, inhibits constitutive and inducible STAT3 activation through induction of a protein tyrosine phosphatase SHP-1. *Cancer Res* 68(11): 4406–15.
- Cui J et al. (2003) Guggulsterone is a farnesoid X receptor antagonist in coactivator association assays but acts to enhance transcription of bile salt export pump. *J Biol Chem* 278(12): 10214–20.
- Gonzalez R et al. (2013) Deriving dopaminergic neurons for clinical use. A practical approach. *Sci Rep* 3: 1463.
- Urizar NL et al. (2002) A natural product that lowers cholesterol as an antagonist ligand for FXR. *Science* 296(5573): 1703–6.
- Xiao D & Singh S V. (2008) z-Guggulsterone, a constituent of Ayurvedic medicinal plant *Commiphora mukul*, inhibits angiogenesis in vitro and in vivo. *Mol Cancer Ther* 7(1): 171–80.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, please visit our website at www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.