

Small Molecules

Purvalanol A

Cyclin/CDK pathway inhibitor; Inhibits CDKs

Catalog # 73772
73774

1 mg
10 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Purvalanol A is a cell-permeable, potent, and selective inhibitor of cyclin-dependent kinases (CDKs). CDKs and cyclins form a stoichiometric complex, which is necessary for the CDK subunit to gain its protein kinase activity. It has been shown that these CDK/cyclin complexes play a key role in initiating G2/M transitions of the cell cycle (Jackman & Pines). Purvalanol A acts through competitive inhibition of ATP binding, to inhibit CDK1/cyclin B (IC_{50} = 4 nM), CDK2/cyclin A (IC_{50} = 70 nM), CDK2/cyclin E (IC_{50} = 35 nM), CDK4/cyclin D1 (IC_{50} = 850 nM), and CDK5-p35 (IC_{50} = 75 nM; Bain et al.; Gray et al.), thereby arresting cells in G1 and G2.

Molecular Name:	Purvalanol A
Alternative Names:	NG 60
CAS Number:	212844-53-6
Chemical Formula:	$C_{19}H_{25}ClN_6O$
Molecular Weight:	388.9 g/mol
Purity:	≥ 98%
Chemical Name:	(2R)-2-[[6-[(3-chlorophenyl)amino]-9-(1-methylethyl)-9H-purin-2-yl]amino]-3-methyl-1-butanol
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 75 mM · Absolute ethanol ≤ 25 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 257 µL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

CANCER RESEARCH

- Inhibits proliferation in exponentially growing cancer cell lines and reversibly arrests synchronised cells in G1 and G2 phase of cell cycle (Villerbu et al.).
- Induces apoptosis in MCF-7 estrogen receptor positive breast cancer cells (Obakan et al.).
- Suppresses cancer progression associated with Src up-regulation by the coordinated inhibition of cell cycle progression and tyrosine kinase signaling (Hikita et al.).

References

- Bain J et al. (2003) The specificities of protein kinase inhibitors: an update. *Biochem J* 371(Pt 1): 199–204.
- Gray NS et al. (1998) Exploiting chemical libraries, structure, and genomics in the search for kinase inhibitors. *Science* 281(5376): 533–8.
- Hikita T et al. (2010) Purvalanol A, a CDK inhibitor, effectively suppresses Src-mediated transformation by inhibiting both CDKs and c-Src. *Genes Cells* 15(10): 1051–62.
- Jackman MR & Pines JN. (1997) Cyclins and the G2/M transition. *Cancer Surv* 29: 47–73.
- Obakan P et al. (2014) Purvalanol A is a strong apoptotic inducer via activating polyamine catabolic pathway in MCF-7 estrogen receptor positive breast cancer cells. *Mol Biol Rep* 41(1): 145–54.
- Villerbu N et al. (2002) Cellular effects of purvalanol A: a specific inhibitor of cyclin-dependent kinase activities. *Int J Cancer* 97(6): 761–9.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.