

Small Molecules

Splitomicin

Epigenetic modifier; Inhibitor of Sir2p

Catalog # 73842

5 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Splitomicin is a cell-permeable, selective inhibitor of Sir2p NAD⁺-dependent histone deacetylase (HDAC) activity. Sir2p negatively regulates gene expression and initiation of DNA replication. Splitomicin induces dose-dependent inhibition of the HDAC activity of Sir2p in vitro (IC₅₀ = 60 μM) and in vivo (minimal inhibitory concentration = 0.49 μM) in yeast extract (Bedalov et al.; Hirao et al.). It was shown that the hydrolytically unstable lactone ring of Splitomicin is critically important for its activity (Posakony et al.).

Molecular Name:	Splitomicin
Alternative Names:	1-Naphthalenepropanoic Acid
CAS Number:	5690-03-9
Chemical Formula:	C ₁₃ H ₁₀ O ₂
Molecular Weight:	198.2 g/mol
Purity:	≥ 98%
Chemical Name:	1,2-dihydro-3H-naphtho[2,1-b]pyran-3-one
Structure:	



Properties

Physical Appearance:	A crystalline solid
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· DMSO ≤ 75 mM · Absolute ethanol ≤ 75 mM For example, to prepare a 10 mM stock solution in DMSO, resuspend 1 mg in 505 μL of DMSO.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in DMSO at -20°C. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

Compound has low solubility in aqueous media. For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use. Avoid final DMSO concentration above 0.1% due to potential cell toxicity.

Published Applications

IMMUNOLOGY

- Inhibits platelet aggregation induced by thrombin, collagen, arachidonic acid, and U46619 (Liu et al.).

DISEASE MODELING

- Reactivates FMR1 expression in neurons of Fragile X patients in vitro (Biacsi et al.).

CANCER RESEARCH

- Sensitizes cancer cells to a variety of DNA-damaging agents by abrogating Sir2p downregulation of p53 expression (Botta et al.).

References

- Bedalov A et al. (2001) Identification of a small molecule inhibitor of Sir2p. *Proc Natl Acad Sci USA* 98(26): 15113–8.
- Biacsi R et al. (2008) SIRT1 inhibition alleviates gene silencing in Fragile X mental retardation syndrome. *PLoS Genet* 4(3): e1000017.
- Botta G et al. (2012) Current advances in the synthesis and antitumoral activity of SIRT1-2 inhibitors by modulation of p53 and pro-apoptotic proteins. *Curr Med Chem* 19(34): 5871–84.
- Hirao M et al. (2003) Identification of selective inhibitors of NAD⁺-dependent deacetylases using phenotypic screens in yeast. *J Biol Chem* 278(52): 52773–82.
- Liu F-C et al. (2009) Splitomicin suppresses human platelet aggregation via inhibition of cyclic AMP phosphodiesterase and intracellular Ca⁺⁺ release. *Thromb Res* 124(2):199–207.
- Posakony J et al. (2004) Inhibitors of Sir2: evaluation of splitomicin analogues. *J Med Chem* 47(10): 2635–44.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2017 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.