

Cytokines

Human Recombinant IL-8 (CXCL8)

Interleukin 8



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 78084
78084.1

5 µg
25 µg

Product Description

Interleukin-8 (IL-8) is a member of the CXC subfamily of chemokines and is produced by leukocytic cells (monocytes, T cells, neutrophils, and natural killer cells) and non-leukocytic somatic cells (endothelial cells, fibroblasts, and epithelial cells), with the most prominent source being monocytes and macrophages. Its production is induced by inflammatory stimuli, such as IL-1. IL-8, also known as CXCL8, activates neutrophils inducing chemotaxis, exocytosis, and the respiratory burst (Baggiolini & Clark-Lewis; Mukaida). IL-8 is considered one of the most potent neutrophil chemoattractants in inflammation and binds to two different chemokine receptors on leukocytes: the G protein-coupled receptors CXCR1 and CXCR2 (Hoffmann et al.; de Oliveira et al.). IL-8 has angiogenic effects on human intestinal microvascular endothelial cells in vitro that are mediated by CXCR2 (Heidemann et al.). IL-8 is reported to promote breast cancer progression by increasing cell invasion, angiogenesis, and metastasis and has been reported to be involved in regulating breast cancer stem-like cells (Singh et al.). IL-8 also has proangiogenic properties in inflammatory diseases of the conjunctiva, cornea, iris, retina, and orbit (Ghasemi et al.). It was also shown that a major T cell effector function in human newborns is IL-8 production, which has the potential to activate antimicrobial neutrophils and gamma/delta T cells (Gibbons et al.). A variety of human pathogens, such as HIV and Mycobacterium tuberculosis, have been shown to induce IL-8 production by monocytes and macrophages (Friedland et al.; Meddows-Taylor et al.).

Product Information

Alternative Names: CXC motif ligand 8, GCP-1, Granulocyte chemotactic protein 1, Interleukin-8, MDNCF, Monocyte-derived neutrophil chemotactic factor, NAF, NAP-1, Neutrophil activating factor, SCYB8, Small inducible cytokine subfamily B member 8

Accession Number: P10145

Amino Acid Sequence: AVLPRSAKEL RCQCIKTYSK PFHPKFIKEL RVIESGPHCA NTEIIVKLSD GRELCCLDPKE NWWQRVVEKF LKRAENS

Predicted Molecular Mass: 8.9 kDa

Species: Human

Cross Reactivity: Mouse

Formulation: Lyophilized after dialysis against phosphate-buffered saline.

Source: E. coli

Specifications

Activity: The specific activity is $\geq 6.7 \times 10^3$ units/mg ($EC_{50} \leq 0.15$ µg/mL) as determined by Ca^{2+} mobilization assay in CHO-K1/Gα15/hCXCR1 cells (human Gα15 and human CXCR1 stably expressed in CHO-K1 cells).

Purity: ≥ 95 %

Endotoxin Level: Measured by kinetic limulus amebocyte lysate (LAL) analysis and is ≤ 0.2 EU/µg protein.


Preparation and Storage

Storage: Store at -80°C.

Stability: Stable as supplied for 12 months from date of receipt.

Preparation: Centrifuge vial before opening. Resuspend the product in sterile water containing 0.1% bovine serum albumin (BSA) to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex. Store at 2 - 8°C for up to 2 weeks or at -20°C to -80°C for up to 3 months. Avoid repeated freeze-thaw cycles. NOTE: If reconstituted product will be used immediately BSA is not required.

Data



(A) The biological activity of Human Recombinant IL-8 (CXCL8) was tested by its ability to mobilize Ca²⁺ in CHO-K1/Gα15/hCXCR1 cells (human Gα15 and human CXCR1 stably expressed in CHO-K1 cells). Ca²⁺ mobilization was measured using a fluorometric assay method. The EC₅₀ is defined as the effective concentration of the growth factor at which Ca²⁺ mobilization is at 50% of maximum. The EC₅₀ in the example above is less than 0.15 µg/mL.

(B) 2 µg of Human Recombinant IL-8 (CXCL8) was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Human Recombinant IL-8 (CXCL8) has a predicted molecular mass of 8.9 kDa.

Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, please visit our website at www.stemcell.com/cytokines or contact us at techsupport@stemcell.com.

References

- Baggiolini M & Clark-Lewis I. (1992) Interleukin-8, a chemotactic and inflammatory cytokine. *FEBS Lett* 307(1): 97–101.
- de Oliveira S et al. (2013) Cxcl8 (IL-8) mediates neutrophil recruitment and behavior in the zebrafish inflammatory response. *J Immunol* 190(8): 4349–59.
- Friedland JS et al. (1992) Secretion of interleukin-8 following phagocytosis of *Mycobacterium tuberculosis* by human monocyte cell lines. *Eur J Immunol* 22(6): 1373–8.
- Ghasemi H et al. (2011) Roles of IL-8 in ocular inflammations: a review. *Ocul Immunol Inflamm* 19(6): 401–12.
- Gibbons D et al. (2014) Interleukin-8 (CXCL8) production is a signatory T cell effector function of human newborn infants. *Nat Med* 20(10): 1206–10.
- Heidemann J et al. (2003) Angiogenic effects of interleukin 8 (CXCL8) in human intestinal microvascular endothelial cells are mediated by CXCR2. *J Biol Chem* 278(10): 8508–15.
- Hoffmann E et al. (2002) Multiple control of interleukin-8 gene expression. *J Leukoc Biol* 72(5): 847–55.
- Meddows-Taylor S et al. (1999) Dysregulated production of interleukin-8 in individuals infected with human immunodeficiency virus type 1 and *Mycobacterium tuberculosis*. *Infect Immun* 67(3): 1251–60.
- Mukaida N. (2003) Pathophysiological roles of interleukin-8/CXCL8 in pulmonary diseases. *Am J Physiol Lung Cell Mol Physiol* 284(4): L566–77.
- Singh JK et al. (2013) Recent advances reveal IL-8 signaling as a potential key to targeting breast cancer stem cells. *Breast Cancer Res* 15(4): 210.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.